

Bion, Commode, 306

BION Code: 02405



Overall Width	22¾" / 58cm
Overall Height	38½" / 97.5cm
Overall Length	32¼" / 82cm
Distance between Armrest to Armrest	18½" / 47cm
Seat Width	15¾" / 40cm
Seat Depth	16½" / 42cm
Seat Height	21¼" / 54cm
Backrest Height	14¼" / 36cm
Clearance for Toilet Bowl	17" / 43cm
Net Weight	9.5kg
Weight Capacity	100kg
Wheels Diameter	4"
Number of Lockable Wheels	4
Armrests Type	Height-Adjustable, Detachable
Footrests	Flip-up
Frame Material	Aluminium
Foldable Frame	No
Seat Cover	Yes
Toilet Bowl Friendly	Yes
Rear Anti-Tipper	No
Front Anti-tipper	Yes
Back Lift	No
Safety Belt	Optional