

## Bion, Blood Pressure Monitor MA350f

Item code: 02406



movement detection



hypertension risk indication



irregular heartbeat detection (IHB)



latex-free patented universal cone cuff

**90M**

90 memories



average of last 3 readings



### Clinically Validated By Multiple Medical Associations

The Bion Blood Pressure Monitor MA350f is a basic blood pressure monitor developed by Bion. It has been clinically validated by multiple medical associations, including BHIS (British and Irish Hypertension Society) and ESH (European Society of Hypertension).

### Specifications:

<b>Device Weight</b>	220 g (without batteries)
<b>Power Source</b>	DC 6V, ≥600mA (Plug size: outer(-) is Ø4.0, inner(+) is Ø1.7) or DC 6V four AAA Batteries
<b>Arm Cuff Circumference</b>	24-40 cm (9.4"~15.7")
<b>Ingress Protection Rating</b>	IP 21
<b>Pressure Measurement Range</b>	30~260mmHg
<b>Pressure Accuracy</b>	±3 mmHg
<b>Pulse Measurement Range</b>	40~199 bpm
<b>Pulse Measurement Accuracy</b>	±5%
<b>Dimensions</b>	129.8mm(L) X 85mm(W) X 68.2mm(H)

### Product Diagram:



Scan QR code for more information on product and warranty

