

BION Non-Contact Temple Thermometer

Model: HB500

- 1 Second Measurement
- Temple Temperature
- Object Temperature
- 9 Memories
- LCD Display with Backlight



- 1 second measurement
- Temple temperature
- Object temperature
- Fever alarm
- 9 memories
- Backlight
- Auto shut-off
- °C / °F switchable
- Low battery indication
- Self-diagnosis for malfunction



1 Second Measurement



Temple Temperature



Object Temperature



Fever Alarm



9 Memories



LCD Display with Backlight



Auto Shut-off



°C / °F Switchable



Low battery indication

Clinically Accurate

ASTM E1965-98:2009

This Temple thermometer has been clinically tested in several big teaching hospitals based on the protocol of ASTM E1965-98:2009 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading.

Model	Qty per carton	Carton volume
HB500	50 pcs	0.079 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
17700 pcs	35400pcs	41100 pcs

Storage and Transportation Condition
 -25°C~-55°C(-13°F~131°F)
 ≤ 95% RH, 700~1060 hPa