

Bion Standard Wheelchair

BION Code: 02424



Seat Width	17³/₄" (45cm)
Seat Depth	16¹/₂" (42cm)
Seat Height	18¹/₂" (47cm)
Backrest Height (From Seat)	15³/₄" (40cm)
Distance between Armrest	17" (43cm)
Overall Width	26" (66cm)
Overall Height	34¹/₄" (87cm)
Overall Length	41¹/₂" (106cm)
Net Weight	16kg
Weight Capacity	100kg
Front Wheel Size	7"
Rear Wheel Size	22" x 1³/₈" Solid Tyres, Spokes Rim
Frame Material	Epoxy Coated Steel
Armrest	Full Length, Padded, Fixed
Footrest	Composite Footplate, Fixed
Foldable Backrest	x
Removable Backrest Cover	x
Removable Seat Cover	x
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	No

1 year warranty on frame only. Wear and tear parts are not covered under warranty.