

bion

Bion Comfy Wheelchair 4G

Item code: 02425



Comfortable & Compact

Bion Comfy Wheelchair 4G is part of our Compact Series. With a foldable footrest and backrest, the unit can be further compacted into a more portable form.


In addition to being compact, this unit has a tension back with seat cushions to improve comfort, making it a luxury portable wheelchair.



Specifications:

Seat Width	15.7" / 40cm
Seat Depth	15.7" / 40cm
Seat Height	17.7" / 45cm
Backrest Height (From Seat)	16.9" / 43cm
Distance between Armrest	17.9" / 45.5cm
Overall Width (Open)	24" / 61cm
Overall Height	34.4" / 87.5cm
Overall Length (Open)	40.2" / 102cm
Net Weight	11.9kg
Weight Capacity	100kg
Front Wheel Size	6"
Rear Wheel Size	20" X 1 1/8" Solid Tyres, Spokes Rim
Frame Material	Aluminium
Armrests	Half-Length, PU, Fixed
Footrests	Composite Footplate, Flip-Up
Foldable Backrest	✓
Removable Backrest Cover	✓
Removable Seat Cover	✓
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	✓

Dimensions:



Warranty Coverage

- 1 Year off-site warranty on manufacturing defects
- Wear & tear parts are not covered under warranty