



Bion Homecare Chair J100

Item Code: 02433


Comfortable and Lighter to Handle

The Bion Homecare Chair J100 is a lighter, more comfortable alternative to traditional geriatric chairs. It features plush cushioning that offers excellent support for the spine and neck, making it ideal for home use as well as clinical settings like dialysis and chemotherapy.

General Specs	
Headrest	Flat
Meal Tray	Blow Molded
Weight Capacity	120kg
Armrests	Fixed, Waterproof, Padded
Frame Material	Non-Toxic Epoxy Coating Steel
Upholstery Material	Fire Retardant
Safety Belt	Optional
Net Weight	20kg
Seat	
Distance Between Armrests	18.9" (48cm)
Seat Width	18.9" (48cm)
Seat Depth	19.3"-21.7" (49cm-55cm)
Seat Height	14.6"-19.3" (37cm-49cm)
Backrest Height (From Seat)	27.2" (69cm)
Dimensions	
Overall Width	25.6" (65cm)
Overall Height	40.6"-45.7" (103cm-116cm)
Overall Length	25.6" (65cm)



Technical Drawings



Scan QR Code for more information on product and warranty