

bion

Bion Comfy Pushchair 4G

Item code: 02435



Comfortable & Compact

Bion Comfy Pushchair 4G is part of our Compact Series. With a foldable footrest and backrest, the unit can be further compacted into a more portable form.


In addition, this unit has a tension back with seat cushions to enhance comfort, making it a luxury portable pushchair.



Specifications:

Seat Width	15.7" / 40cm
Seat Depth	16.1" / 41cm
Seat Height	16.9" / 43cm
Backrest Height (From Seat)	16.7" / 42.5cm
Distance between Armrest	17.5" / 44.5cm
Overall Width (Open)	21.7" / 55cm
Overall Height	34.3" / 87cm
Overall Length (Open)	39.8" / 101cm
Net Weight	11.2kg
Weight Capacity	100kg
Front Wheel Size	6"
Rear Wheel Size	12" Solid Tyres, Spokes Rim
Frame Material	Aluminium
Armrests	Half-Length, PU, Fixed
Footrests	Composite Footplate, Flip-Up
Foldable Backrest	✓
Removable Backrest Cover	✓
Removable Seat Cover	✓
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	✓

Dimensions:



Warranty Coverage

- 1 Year off-site warranty on manufacturing defects
- Wear & tear parts are not covered under warranty