

bion

Bion Traveller Elite

Item code: 02442




Lightweight & Compact

Bion Traveller Elite is our premium travel companion. Being lightweight and compact with a royal yellow finish, you will be glad to bring this around anywhere and anytime as you travel.



Dimensions:



Specifications:

Seat Width	15.4" / 39cm
Seat Depth	14.8" / 37.5cm
Seat Height	18.5" / 47cm
Backrest Height (From Seat)	15" / 38cm
Distance between Armrest	15.7" / 40cm
Overall Width (Open)	18.9" / 48cm
Overall Height	34" / 86.5cm
Overall Length (Open)	37" / 94cm
Net Weight	9.6kg
Weight Capacity	80kg
Front Wheel Size	6"
Rear Wheel Size	8" Solid Tyres, Mag Rim
Frame Material	Aluminium
Armrests	Half-Length, PU, Flip-Up
Footrests	Flip-Up
Foldable Backrest	✓
Removable Backrest Cover	X
Removable Seat Cover	X
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	X

Warranty Coverage

- 1 Year off-site warranty on manufacturing defects
- Wear & tear parts are not covered under warranty