



Bion Finger Pulse Oximeter P100

Item code: 02506 / 02478



oxygen saturation



perfusion index



singapore brand



pulse rate




threshold alerts



1 year warranty

Product Diagram:



Simple Assuring

The Bion Finger Pulse Oximeter P100 is utilised for measuring arterial oxygen saturation (% SpO₂) of hemoglobin and pulse rate, which are important indicators of your respiratory function.

Specifications:

Display Range	Oxygen Saturation	35% to 99%
	Pulse Rate	25 bpm to 250 bpm
Resolution	Oxygen Saturation	1%
	Pulse Rate	1 bpm
Measurement Precision	Oxygen Saturation	±2% (70% to 99%) no requirement (≤69%)
	Pulse Rate	±2bpm
Alert Range	Oxygen Saturation	Upper limit: 50% to 100% Lower limit: 50% to 98%
	Pulse Rate	Upper limit: 25 bpm to 250 bpm Lower limit: 25 bpm to 250 bpm
Alert Error	Oxygen Saturation	±1% of the preset value
	Pulse Rate	The greater of ±10% of the preset value and ±5 bpm
PI	Weak PI	Min. 0.3%

Scan QR code for more information on product and warranty

