

BION Automatic Upper Arm Blood Pressure Monitor

Model: B100



clinically validated



movement detection



99 memories x 2 users



irregular heartbeat detection



hypertension risk indication



universal cuff*



one touch operation



5 years warranty


SG

singapore brand


Specifications

Measurement Method	Oscillometric
Measurement Range	Pressure: 0-295mmHg; Pulse: 40-199bpm
Accuracy	Pressure: ± 3mmHg; Pulse: ± 5% of reading
Memory Capacity	99 Memories x 2 Users
Auto Shut-off	30 seconds after last operation
Operation Environment	5°C-40°C (41°F~104°F); 15%-85% RH; 70-106kPa
Storage & Transportation Environment	-20°C-55°C (-4°F~131°F); 10%~93% RH; 70-10kPa
Power Source	4x AA Dry Batteries or DC 5V USB external power supply (not included)
Monitor Dimensions	138.5(L) X 118.5(W) X 56.0(H) mm
Device Weight (Without Batteries)	362g
Arm cuff circumference	22-36cm
Ingress Protection Rating	IP 21
Classification	Type BF Applied Part

Product Diagram



Display Interface



Packing List

Description	Quantity
Electronic Blood Pressure Monitor	1
Cuff	1
Dry Battery (AA)	4
Instruction Manual	1
Zippered Carrying Bag	1