

BION Automatic Upper Arm Blood Pressure Monitor

Model: B100



clinically validated



movement detection



99 memories x 2 users



irregular heartbeat detection



hypertension risk indication



universal cuff*



one touch operation



5 years warranty


SG

singapore brand

Specifications

| | |
|--------------------------------------|---|
| Measurement Method | Oscillometric |
| Measurement Range | Pressure: 0-295mmHg; Pulse: 40-199bpm |
| Accuracy | Pressure: ± 3mmHg; Pulse: ± 5% of reading |
| Memory Capacity | 99 Memories x 2 Users |
| Auto Shut-off | 30 seconds after last operation |
| Operation Environment | 5°C-40°C (41°F~104°F); 15%-85% RH; 70-106kPa |
| Storage & Transportation Environment | -20°C-55°C (-4°F~131°F); 10%~93% RH; 70-10kPa |
| Power Source | 4x AA Dry Batteries or DC 5V USB external power supply (not included) |
| Monitor Dimensions | 138.5(L) X 118.5(W) X 56.0(H) mm |
| Device Weight (Without Batteries) | 362g |
| Arm cuff circumference | 22-42 cm |
| Ingress Protection Rating | IP 21 |
| Classification | Type BF Applied Part |

Product Diagram



Packing List

| Description | Quantity |
|-----------------------------------|----------|
| Electronic Blood Pressure Monitor | 1 |
| Cuff | 1 |
| Dry Battery (AA) | 4 |
| Instruction Manual | 1 |
| Zippered Carrying Bag | 1 |

Display Interface

