



# Bion Overbed Table, 001

Item code: 02295



## Specifications:

Overall Width	15.7" / 40cm
Overall Height	28.5"-38.2" / 72.50-97cm
Overall Length	31.5" / 80cm
Table Top Width	15.7" / 40cm
Table Top Height	24.4"-33.9" / 62-86cm
Table Top Length	31.5" / 80cm
Table Top Thickness	0.8" / 2cm
Base Width	15.4" / 39cm
Base Height	3.9" / 10cm
Base Length	29.5" / 75cm
Height Adjustable	Yes
Table Material	Laminated Wood
Frame Material	Aluminium Alloy
Number of Lockable Wheels	4




Side View



Elevated View

## Technical Drawings



Ver20240102

**Warranty Coverage**  
• 1 Year off-site warranty on manufacturing defects