



Bion, Homecare Bed, H100

BION Code: 02525

Features

- Versatile homecare bed that comes with wooden bed frame
- Remote controlled adjustable positions:
 1. Backrest Angle
 2. Elevating Leg Rest
 3. Bed Height
- Wired and wireless remote controller included
- Quad Side Rails
- Four lockable castors

Maximum Weight Capacity	150 kg
Dimensions (L X W)	210cm x 102cm
Gross Weight	120kg
Mattress Weight	15 kg
Mattress Thickness	3 inch
Lowest Height (from ground to top of mattress)	43cm
Highest Height (from ground to top of mattress)	81 cm
Back Recline Angle	75°
Leg Elevating Angle	30°
Frame Finishing	Powder Coated Steel
Locking System	4 Lockable Castors 4"

1 year warranty on frame only. Wear and tear parts are not covered under warranty.