

## Bion, Commode, C100

**BION Code: 02526**



<b>Overall Width</b>	<b>27½" / 70cm</b>
<b>Overall Height</b>	<b>30¼" - 34¾" / 77-88cm</b>
<b>Overall Length</b>	<b>20¾" / 53cm</b>
<b>Distance between Armrest to Armrest</b>	<b>17" / 43cm</b>
<b>Seat Width</b>	<b>17" / 43cm</b>
<b>Seat Depth</b>	<b>15" / 38cm</b>
<b>Seat Height</b>	<b>16½" - 21" / 42-53cm</b>
<b>Backrest Height</b>	<b>13¾" / 35cm</b>
<b>Clearance for Toilet Bowl</b>	<b>5" - 9.5" / 13-24cm</b>
<b>Net Weight</b>	<b>5.6kg</b>
<b>Weight Capacity</b>	<b>120kg</b>
<b>Wheels Diameter</b>	<b>Nil</b>
<b>Number of Lockable Wheels</b>	<b>Nil</b>
<b>Armrests Type</b>	<b>ABS Armrest Fixed</b>
<b>Footrests</b>	<b>Nil</b>
<b>Frame Material</b>	<b>ABS/Aluminium (Legs)</b>
<b>Foldable Frame</b>	<b>No</b>
<b>Seat Cover</b>	<b>Yes</b>
<b>Toilet Bowl Friendly</b>	<b>Yes</b>
<b>Rear Anti-Tipper</b>	<b>No</b>
<b>Front Anti-tipper</b>	<b>No</b>
<b>Back Lift</b>	<b>No</b>
<b>Safety Belt</b>	<b>No</b>