



Bion Overbed Table D100

Item code: 02582




Specifications:

Overall Width	17.3" / 44cm
Overall Height	31.1"-44.9" / 79-114cm
Overall Length	29.9" / 76cm
Table Top Width	15" / 38cm
Table Top Height	27.2"-40.9" / 69-104cm
Table Top Length	29.9" / 76cm
Table Top Thickness	0.8" / 2cm
Base Width	17.3" / 44cm
Base Height	3.9" / 10cm
Base Length	28" / 71cm
Height Adjustable	Yes
Table Material	Waterproof Desktop
Frame Material	High Quality Powder Coated Steel Frame
Number of Lockable Wheels	2



Side View

Technical Drawings



Ver20240102

Warranty Coverage
• 1 Year off-site warranty on manufacturing defects