

bion

Bion, Postur Pushchair, 16", S310

Item code: 02589




Specifications

Seat Width	16" / 40.6cm
Seat Depth	16.5" / 42cm
Seat Height	19.7" / 50cm
Backrest Height (From Seat)	17.5" / 44.5cm
Distance between Armrest	18.7" / 47.5cm
Overall Width (Open)	23.4" / 59.5cm
Overall Height	34.8" / 88.5cm
Overall Length (Open)	39.8" / 101cm
Net Weight	11.8kg
Weight Capacity	100kg
Front Wheel Size	6" PU
Rear Wheel Size	16" X 1.75 PU Tyres, Spokes Rim
Frame Material	Aluminium Alloy
Armrests	Half-Length, Aluminium Alloy / PU Pad, Flip-Up
Footrests	Plastic Composite Footplate, Detachable
Foldable Backrest	✓
Removable Backrest Cover	✓
Removable Seat Cover	✓
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	✓

Postur Seating

Ergonomically designed with your comfort as our top priority with every aspect tailored to support your body's natural alignment.

Curve Seating System



Tilted Backrest & Seat Angle

Warranty Coverage

- 1 Year off-site warranty on manufacturing defects
- Wear & tear parts are not covered under warranty

Dimensions:

