

Bion iLight Pushchair
L210

Item Code: 02630





Mobility, Light as a Feather

Travel effortlessly with Bion iLight Pushchair L210. At a mere 8.6kg, this pushchair with 16" rear wheels is the lightest in its class worldwide. Featuring an ultra-lightweight magnesium alloy frame and detachable footrests for added convenience, this pushchair is easy to lift and transport—making it perfect for transfers and everyday use.

Specifications

General Specs	
Net Weight	8.6kg
Weight Capacity	100kg
Frame Material	Magnesium Alloy
Armrest Material	Magnesium Alloy/PU Pad
Colour	Sky Blue
Foldable Backrest	Yes
Seat	
Seat Width	18.1" (46cm)
Seat Depth	15" (38cm)
Seat Height	19.3" (49cm)
Backrest Height (From Seat)	16.9" (43cm)
Distance between Armrest	17.9" (45.5cm)
Wheel	
Front Wheel Size	6" PVC
Rear Wheel Size	16" X 1.75, Spoke PU
Anti-Tipper	NA
Dimensions	
Overall Width	25.6" (65cm)
Overall Height	35.4" (90cm)
Overall Length	38.8" (98.5cm)
Folded Length (With Legrest)	35.8" (91cm)
Folded Width	11.4" (29cm)
Folded Height	30.3" (77cm)

Technical Drawings



Scan QR Code for more information on product and warranty