

Bion Code	02318	02437	00646	02311	02435	02353	01224	02297	02442
Series	iLight	iLight	iLight	Comfy	Comfy	iLight Detach	iLight Detach	Traveller	Traveller
Intended Purpose	General Purpose	General Purpose	General Purpose	Urban Travel	Urban Travel	Patient-Use	Patient-Use	Overseas Travelling	Overseas Travelling
Unique Feature	Lightweight. Economy. -.	Lightweight. Economy. -.	Comfortable. Lightweight. -.	Lightweight. Compact. -.	Lightweight. Tension Back. Anti-Tippers.	Flip-Up Armrests. Removable Leg Rests. Anti-Tippers.	Flip-Up Armrests. Removable Leg Rests. Anti-Tippers.	Lightweight. Carrying Bag. Compact.	Lightweight. Carrying Bag. Flip-Up Armrests.
Image									
Description	Bion, iLight Pushchair, EZ	Bion, iLight Pushchair, 3G	Bion, iLight Pushchair	Bion, Comfy Pushchair, 3G	Bion, Comfy Pushchair, 4G	Bion, iLight Pushchair, Detachable, 16"	Bion, iLight Pushchair, Detachable, 18"	Bion Traveller	Bion Traveller Elite
Seat Width	17½" (45cm)	16½" (42cm)	15½" (40cm)	15¾" (40cm)	15¾" (40cm)	15½" (39cm)	18" (46cm)	14" (36cm)	15¾" (39cm)
Seat Depth	16¾" (43cm)	15½" (39cm)	16" (41cm)	15¾" (40cm)	16" (41cm)	17" (43cm)	17" (43cm)	14¾" (37.5cm)	14¾" (37.5cm)
Seat Height	18" (46cm)	18¾" (47.5cm)	19" (48cm)	17½" (44.5cm)	17" (43cm)	18.5" (47cm)	18" (46cm)	18½" (47cm)	18½" (47cm)
Backrest Height (from Seat)	15¼" (39cm)	17¼" (44cm)	16" (41cm)	16¼" (41.5cm)	16¾" (42.5cm)	15" (38cm)	16¼" (42cm)	16" (41cm)	15" (38cm)
Distance between Armrest	16¾" (43cm)	15¾" (40cm)	17¾" (45cm)	18" (45.5cm)	17½" (44.5cm)	16" (41cm)	18½" (47cm)	14¾" (37.5cm)	15¾" (40cm)
Overall Width	25" (65cm)	24¾" (61.5cm)	23½" (60cm)	21½" (54.5cm)	21¾" (55cm)	21" (53cm)	24" (61cm)	17" (43cm)	19" (48cm)
Overall Height	33½" (85cm)	35" (89cm)	35" (89cm)	34½" (87.5cm)	34¾" (87cm)	34" (86cm)	33½" (85cm)	34½" (87.5cm)	34" (86.5cm)
Overall Length	40½" (103cm)	39½" (100cm)	39½" (100cm)	39" (99cm)	39¾" (101cm)	40½" (103cm)	40½" (103cm)	33" (83.75cm)	37" (94cm)
Net Weight	10.65kg	10.6kg	11kg	10.4kg	11.2kg	13kg	13kg	7.3kg	9.6kg
Weight Capacity	100kg	100kg	100kg	100kg	100kg	100kg	100kg	90kg	100kg
Front Wheel Size	6"	7"	6"	6"	6"	6"	6"	5"	6"
Rear Wheel Size	12" Solid Tyres, Mag Rim	12" Solid Tyres, Mag Rim	12" Solid Tyres, Mag Rim	12" Solid Tyres, Mag Rim	12" Solid Tyres, Spokes Rim	12" Solid Tyres, Mag Rim	12" Solid Tyres, Mag Rim	8" Solid Tyres, Mag Rim	8" Solid Tyres, Mag Rim
Frame Material	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium
Armrest	Half-Length, Padded, Fixed	Half-Length, PU, Fixed	Half-Length, PU, Fixed	Half-Length, Padded, Fixed	Half-Length, PU, Fixed	Half-Length, Padded, Flip-Up	Half-Length, Padded, Flip-Up	Half-Length, Rubber, Flip-Up	Half-Length, PU, Flip-Up
Footrest	Composite Footplate, Fixed	Composite Footplate, Fixed	Composite Footplate, Fixed	Composite Footplate, Flip-Up	Composite Footplate, Flip-Up	Composite Footplate, Detachable	Composite Footplate, Detachable	Flip-Up	Flip-Up
Foldable Backrest	✓	✓	✓	✓	✓	✓	✓	✓	✓
Removable Backrest Cover	✗	✗	✓	✗	✓	✓	✓	✗	✗
Removable Seat Cover	✗	✗	✓	✗	✓	✓	✓	✗	✗
User Brakes	✓	✓	✓	✓	✓	✓	✓	✗	✓
Attendant Brakes	✓	✓	✓	✓	✓	✓	✓	✓	✓
Safety Belt	✓	✓	✓	✓	✓	✓	✓	✓	✓
Anti-Tipper	Optional	Optional	No	Optional	Yes	Yes	Yes	No	No
Tension Backrest	✗	✗	✗	✗	✓	✗	✗	✗	✗
Quick Release Wheels	✗	✗	✗	✗	✗	✗	✗	✗	✗
Updated as of	11/10/2018	30/08/2019	6/10/2020	15/02/2019	30/08/2019	19/08/2020	6/07/2020	18/03/2020	18/03/2020