



Bion, Homecare Bed, H110

Item code: 02581



Specifications:

Max. Weight Capacity	150 kg
Dimensions (L x W)	210 cm x 102 cm
Gross Weight	120 kg
Mattress (Single)	192cm (L) x 85cm (W) x 7.5cm (H)
Mattress Weight	15 kg
Mattress Thickness	3 inch
Lowest Height (from ground to top of mattress)	41 cm
Highest Height (from ground to top of mattress)	80 cm
Back Recline Angle	75°
Leg Elevating Angle	30°
Frame Finishing	Powder Coated Steel
Locking System	4 Lockable Castors 4"
Bed Base to Side Rail Height	11.6" / 29.5cm
Mattress to Side Rail Height	8.7" / 22cm



Lower Body Elevation



Upper Body Elevation



Adjustable Bed Height

Warranty Coverage

- Wired controller and electronics (1 year on-site warranty)
- Motor and actuator (2 years warranty)
- Frame (4 years warranty)