

bion

Bion, Portable Water Flosser W100


Item code: 02503



Specifications:

Power Consumption	5V, 500mA
Battery Capacity	2.4v 1300mA/h
Water Pressure	Average - approx. 390 kPa (57psi) Strong - approx. 590 kPa (86psi)
Water Jet Pulsation Frequency	1500 pulses/min
Built-In Container Capacity	160 ml
Waterproof Rating	IPX5
Weight (without packaging)	304 g
Noise Level	Less than 65 dB
Package Contents	Main Unit, Connector, Container, 3 Different Nozzles, Long Tube With Holder Short Tube, Power Cord (USB), Instruction Manual, Individual Carton Package.

Product Diagram:



Water Flosser Nozzles:

Standard Nozzle W101




- For general use
- Made of plastic
- For individual use

Subgingival Nozzle W102



- For cleaning difficult areas
- Made of plastic with flexible nozzle
- For individual use

Orthodontic Nozzle W103



- For clearing dental brackets
- Made of plastic with brush head
- For individual use



Scan QR code for more information on product and warranty