

Bion, Commode, with Wheels, 103

Item code: 02382



Do not use on uneven surfaces





Do not push over ramps



Specifications:

Overall Width	20.9" / 53cm
Overall Height	36.4"-40.2" / 92.5-102cm
Overall Length	22.8" / 58cm
Distance between Armrest to Armrest	17.3" / 44cm
Seat Width	15.4" / 39cm
Seat Depth	16.5" / 42cm
Seat Height	18.3"-22.2" / 46.5-56.5cm
Backrest Height	18.5" / 47cm
Clearance for Toilet Bowl	14.3"-18.1" / 36.5-46cm
Net Weight	7.2kg
Weight Capacity	100kg
Wheels Diameter	4"
Number of Lockable Wheels	4
Armrests Type	ABS Armrest, Fixed
Footrests	NIL
Frame Material	Aluminium
Foldable Frame	✓
Seat Cover	✓
Toilet Bowl Friendly	✓
Rear Anti-Tipper	X
Front Anti-Tipper	X
Back Lift	X
Safety Belt	Optional

Technical Drawings:



Warranty Coverage

- 1 Year off-site warranty on manufacturing defects
- Wear & tear parts are not covered under warranty



Scan QR code for more information