

bion

Bion, Homecare Bed, H100

Item code: 02525



Specifications:

Max. Weight Capacity	150 kg
Dimensions (L x W)	210 cm x 102 cm
Gross Weight	120 kg
Mattress Weight	15 kg
Mattress Thickness	3 inch
Lowest Height <small>(from ground to top of mattress)</small>	43 cm
Highest Height <small>(from ground to top of mattress)</small>	81 cm
Back Recline Angle	75°
Leg Elevating Angle	30°
Frame Finishing	Powder Coated Steel
Locking System	4 Lockable Castors 4"




Lower Body Elevation



Upper Body Elevation

Technical Drawing:



Adjustable Bed Height

Scan QR code for more information on product and warranty

