

bion

Bion, Pressure Relief Mattress, V450

Item code: 02386



prevention & treatment of mild pressure sores*



static or alternating inflation



adjustable firmness



low vibration



low noise



waterproof mattress cover



CPR rapid deflation





protective zipper cover



replaceable tubes

Product Diagram:



Motor Pump Specifications:

Size (L x W x H):	27 x 13 x 9 cm
Weight:	1.4 kg
Material:	ABS Plastic
Power Supply:	220-240V, 50Hz
Fuse Rating:	T500mAH, 250V
Cycle Time:	12 mins
Pressure Range:	25-140 mmHg

Atmospheric Pressure

80KPa to 106KPa

Temperature

Operation: 5°C to 35°C (41°F to 95°F)
Storage: -15°C to 50°C (5°F to 122°F)
Shipping: -15°C to 70°C (5°F to 158°F)

Humidity

Operation: 10% to 90% non-condensing
Storage: 10% to 90% non-condensing
Shipping: 10% to 90% non-condensing

Pressure Relief Mattress Specifications:

Dimension (L x W x H):	192 x 83 x 11.4 cm
Weight:	4.2 kg
Max. Support Weight	145 kg
Cell Material	TPU
Cover	PU coated Nylon

Scan QR code for more information on product and warranty

