









Bion Code	02318	02437	00646	02426	02311	02435	02353	01224
Series	iLight	iLight	iLight	Comfy	Comfy	Comfy	iLight Detach	iLight Detach
Intended Purpose	General Purpose	General Purpose	General Purpose	Urban Travel	Urban Travel	Urban Travel	Patient-Use	Patient-Use
Unique Feature	Lightweight. Economy.	Lightweight. Economy.	Comfort. Lightweight. Easy-Wash.	Lightweight. Compact.	Lightweight. Compact.	Lightweight. Compact. Anti-tippers. Tension Back.	Flip-up Armrests. Removable leg rest. Anti-tippers.	Flip-up Armrests. Removable leg rest. Anti-tippers.
Image								
Description	Bion, iLight Pushchair, EZ	Bion, iLight Pushchair, 3G	Bion, iLight Pushchair	Bion, Comfy Pushchair, 2G	Bion, Comfy Pushchair, 3G	Bion, Comfy Pushchair, 4G	Bion, iLight Pushchair, Detachable, 16"	Bion, iLight Pushchair, Detachable, 18"
Seat Width	17¾" (45cm)	16½" (42cm)	15¾" (40cm)	15¾" (40cm)	15¾" (40cm)	15¾" (40cm)	15½" (39cm)	18" (46cm)
Seat Depth	16¾" (43cm)	15¾" (39cm)	16" (41cm)	15½" (39.5cm)	15¾" (40cm)	16" (41cm)	17" (43cm)	17" (43cm)
Seat Height	18" (46cm)	18¾" (47.5cm)	19" (48cm)	17½" (44.5cm)	17½" (44.5cm)	17" (43cm)	18.5" (47cm)	18" (46cm)
Backrest Height (from Seat)	15¾" (39cm)	17¾" (44cm)	16" (41cm)	17½" (44.5cm)	16¾" (41.5cm)	16¾" (42.5cm)	15" (38cm)	16¾" (42cm)
Distance between Armrest	16¾" (43cm)	15¾" (40cm)	17¾" (45cm)	17¾" (44cm)	18" (45.5cm)	17½" (44.5cm)	16" (41cm)	18½" (47cm)
Overall Width	25" (65cm)	24¾" (61.5cm)	23½" (60cm)	21¾" (55cm)	21½" (54.5cm)	21¾" (55cm)	21" (53cm)	24" (61cm)
Overall Height	33½" (85cm)	35" (89cm)	35" (89cm)	34¾" (88cm)	34¾" (87.5cm)	34¾" (87cm)	34" (86cm)	33½" (85cm)
Overall Length	40½" (103cm)	39½" (100cm)	39½" (100cm)	39¾" (99.5cm)	39" (99cm)	39¾" (101cm)	40½" (103cm)	40½" (103cm)
Net Weight	10.65kg	10.6kg	11kg	10.5kg	10.4kg	11.2kg	13kg	13kg
Weight Capacity	100kg	100kg	100kg	100kg	100kg	100kg	100kg	100kg
Front Wheel Size	6"	7"	6"	6"	6"	6"	6"	6"
Rear Wheel Size	12" Solid Tyres, Mag Rim	12" Solid Tyres, Mag Rim	12" Solid Tyres, Mag Rim	12" Solid Tyres, Mag Rim	12" Solid Tyres, Mag Rim	12" Solid Tyres, Spokes Rim	12" Solid Tyres, Mag Rim	12" Solid Tyres, Mag Rim
Frame Material	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium
Armrest	Half-Length, Padded, Fixed	Half-Length, PU, Fixed	Half-Length, PU, Fixed	Half-Length, Padded, Fixed	Half-Length, Padded, Fixed	Half-Length, PU, Fixed	Half-Length, Padded, Flip-Up	Half-Length, Padded, Flip-Up
Footrest	Composite Footplate, Fixed	Composite Footplate, Fixed	Composite Footplate, Fixed	Composite Footplate, Fixed	Composite Footplate, Flip-Up	Composite Footplate, Flip-Up	Composite Footplate, Detachable	Composite Footplate, Detachable
Foldable Backrest	✓	✓	✓	✓	✓	✓	✓	✓
Removable Backrest Cover	✗	✗	✓	✗	✗	✓	✓	✓
Removable Seat Cover	✗	✗	✓	✗	✗	✓	✓	✓
User Brakes	✓	✓	✓	✓	✓	✓	✓	✓
Attendant Brakes	✓	✓	✓	✓	✓	✓	✓	✓
Safety Belt	✓	✓	✓	✓	✓	✓	✓	✓
Anti-Tipper	Optional	Optional	No	Optional	Optional	Yes	Yes	Yes
Updated as of	11/10/2018	30/08/2019	6/10/2020	12/07/2019	15/02/2019	30/08/2019	19/08/2020	6/07/2020

Bion Code	02297	02442
Series	Traveller	Traveller
Intended Purpose	Overseas travelling	Overseas travelling
Unique Feature	Lightweight. Compact. Carrying Bag.	Lightweight. Flip-up Armrests Carrying Bag.
Image		
Description	Bion Traveller	Bion Traveller Elite
Seat Width	14" (36cm)	15¾" (39cm)
Seat Depth	14¾" (37.5cm)	14¾" (37.5cm)
Seat Height	18½" (47cm)	18½" (47cm)
Backrest Height (from Seat)	16" (41cm)	15" (38cm)
Distance between Armrest	14¾" (37.5cm)	15¾" (40cm)
Overall Width	17" (43cm)	19" (48cm)
Overall Height	34½" (87.5cm)	34" (86.5cm)
Overall Length	33" (83.75cm)	37" (94cm)
Net Weight	7.3kg	9.6kg
Weight Capacity	90kg	100kg
Front Wheel Size	5"	6"
Rear Wheel Size	8" Solid Tyres, Mag Rim	8" Solid Tyres, Mag Rim
Frame Material	Aluminium	Aluminium
Armrest	Half-Length, Rubber, Flip-Up	Half-Length, PU, Flip-Up
Footrest	Flip-Up,	Flip-Up,
Foldable Backrest	✓	✓
Removable Backrest Cover	✗	✗
Removable Seat Cover	✗	✗
User Brakes	✗	✓
Attendant Brakes	✓	✓
Safety Belt	✓	✓
Anti-Tipper	No	No
Updated as of	18/03/2020	18/03/2020