

OXYGEN CONCENTRATOR

8F-5A



Global sales leader



Medical grade



High oxygen concentration



5L oxygen flow



Failure detection



Ultra-silence



Castors



- Flow Range: 0.5 ~ 5L/min
- Oxygen Concentration: 95.5% ~ 87%
- Maximum Limited Pressure: 70kPa
- Sound Level: 5L/min Sound Pressure ≤ 54dB(A)
- Input Power: 350VA
- Minimum Operating Time: 15minutes

PACKAGING PARAMETERS

Model	QTN (PCS/CTN)	G.W/N.W (KG)	Size (cm)
8F-5A	1	19/16.5	47.5×33×58.5