

Bion, iLight Pushchair

BION Code: 00646



Seat Width	15³/₄" (40cm)
Seat Depth	16" (41cm)
Seat Height	19" (48cm)
Backrest Height (From Seat)	16" (41cm)
Distance between Armrest	17³/₄" (45cm)
Overall Width	23¹/₂" (60cm)
Overall Height	35" (89cm)
Overall Length	39¹/₂" (100cm)
Net Weight	11kg
Weight Capacity	100kg
Front Wheel Size	6"
Rear Wheel Size	12" Solid Tyres, Mag Rim
Frame Material	Aluminium
Armrest	Half-Length, PU, Fixed
Footrest	Composite Footplate, Fixed
Foldable Backrest	✓
Removable Backrest Cover	✓
Removable Seat Cover	✓
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	Optional

1 year warranty on frame only. Wear and tear parts are not covered under warranty.