

## Bion, Blood Pressure Monitor MA350f

Item code: 02406



movement detection



hypertension risk indication



irregular heartbeat detection (IHB)



latex-free patented universal cone cuff


**90M**

90 memories



average of last 3 readings

### Product Diagram:




### Clinically Validated By Multiple Medical Associations

The Bion Blood Pressure Monitor MA350f is a basic blood pressure monitor developed by Bion. It has been clinically validated by multiple medical associations, including BHS (British Hypertension Society) and ESH (European Validation Protocols).

### Specifications:

Device Weight	220 g (without batteries)
Power Source	DC 6V, ≥600mA (Plug size: outer(-) is Ø4.0, inner(+) is Ø1.7) or DC 6V four AAA Batteries
Arm Cuff Circumference	24-40 cm (9.4"~15.7")
Ingress Protection Rating	IP 21
Pressure Measurement Range	30~260mmHg
Pressure Accuracy	±3 mmHg
Pulse Measurement Range	40~199 bpm
Pulse Measurement Accuracy	±5%
Dimensions	129.8mm(L) X 85mm(W) X 68.2mm(H)



Scan QR code for more information on product and warranty

