

Bion Digital Oral Thermometer LB100

Item code: 02410



water resistant



10 sec fast read*



fever alarm



flexible tip



auto shut off



2 years warranty



Fast & Accurate


Bion Digital Oral Thermometer with transparent casing is a water-resistant thermometer that can be used anywhere and anytime.

It has a fast pace temperature reader that can process temperature readings within seconds, allowing users to save a lot of waiting time.

Specifications:

Type	Digital thermometer (not predictive)
Measure Range	32.0°C~42.9°C (90.0°F~109.9°F)
Accuracy	±0.1°C 35.5°C~42.0°C (±0.2°F 95.9°F~107.6°F) ±0.2°C under 35.5°C or over 42.0°C (±0.4°F under 95.9°F or over 107.6°F)
Operating Mode	Direct mode
Display	Liquid crystal display
Memory	For storing the last measured value
Battery	One 1.5 V DC. button battery (size LR41 or SR41, or equivalent)
Battery Life	Approx. 200 hours of continuous operation or 1 year with 3 measurements per day
Dimension	13.9cm x 2.2cm x 1.2cm (L x W x H)
Weight	Approx. 12g including battery
Expected Service Life	Three years
Ambient Operating Range	Temperature: 5°C~40°C (41°F~104°F) Relative Humidity: 15%~95%RH Atmospheric Pressure: 700hPa ~ 1060hPa
Storage & Transportation Condition	Temperature: -20°C~55°C (-4°F~131°F) Relative Humidity: 15%~95%RH Atmospheric Pressure: 700hPa ~ 1060hPa
Ingress Protection Rating	IP 27
Classification	Type BF

Product Diagram:



Scan QR code for more information on product and warranty

