

Bion, Comfy Wheelchair, 4G

BION Code: 02425



Seat Width	15 ³ / ₄ " (40cm)
Seat Depth	15 ³ / ₄ " (40cm)
Seat Height	17 ³ / ₄ " (45cm)
Backrest Height (From Seat)	17" (43cm)
Distance between Armrest	18" (45.5cm)
Overall Width	24" (61cm)
Overall Height	34 ¹ / ₂ " (87.5cm)
Overall Length	40 ¹ / ₄ " (102cm)
Net Weight	11.9kg
Weight Capacity	100kg
Front Wheel Size	6"
Rear Wheel Size	20" x 1 ³ / ₈ " Solid Tyres, Spokes Rim
Frame Material	Aluminium
Armrest	Half-Length, PU, Fixed
Footrest	Composite Footplate, Flip-Up
Foldable Backrest	✓
Removable Backrest Cover	✓
Removable Seat Cover	✓
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	Optional

1 year warranty on frame only. Wear and tear parts are not covered under warranty.