

Bion, Comfy Wheelchair, 2G

BION Code: 02436



Seat Width	15 ³ / ₄ " (40cm)
Seat Depth	15 ¹ / ₂ " (39.5cm)
Seat Height	17 ¹ / ₂ " (44.5cm)
Backrest Height (From Seat)	16 ¹ / ₂ " (42cm)
Distance between Armrest	17 ³ / ₄ " (45cm)
Overall Width	23 ¹ / ₂ " (60cm)
Overall Height	34 ³ / ₄ " (88cm)
Overall Length	38 ¹ / ₂ " (98cm)
Net Weight	10.9kg
Weight Capacity	100kg
Front Wheel Size	6"
Rear Wheel Size	20" x 1 ³ / ₈ " Solid Tyres, Spokes Rim
Frame Material	Aluminium
Armrest	Half-Length, Padded, Fixed
Footrest	Composite Footplate, Flip-Up
Foldable Backrest	✓
Removable Backrest Cover	✗
Removable Seat Cover	✗
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	Optional

1 year warranty on frame only. Wear and tear parts are not covered under warranty.