

bion

Bion, Postur Wheelchair, 18", S300

Item code: 02553



Specifications

Seat Width	18.1" / 46cm
Seat Depth	16.9" / 43cm
Seat Height	19.7" / 50cm
Backrest Height (From Seat)	17.3" / 44cm
Distance between Armrest	20.7" / 52.5cm
Overall Width (Open)	27.2" / 69cm
Overall Height	35" / 89cm
Overall Length (Open)	39.8" / 101cm
Net Weight	13.3kg
Weight Capacity	100kg
Front Wheel Size	6" PU
Rear Wheel Size	22" X 1 1/8" PU Tyres, Spokes Rim
Frame Material	Aluminium
Armrests	Half-Length, PU, Flip-Up
Footrests	Composite Footplate, Detachable
Foldable Backrest	✓
Removable Backrest Cover	✓
Removable Seat Cover	✓
User Brakes	✓
Attendant Brakes	✓
Safety Belt	✓
Anti-Tipper	✓

Postur Seating

Ergonomically designed with your comfort as our top priority with every aspect tailored to support your body's natural alignment.



Warranty Coverage

- 1 Year off-site warranty on manufacturing defects
- Wear & tear parts are not covered under warranty

Dimensions:

