

# bion

## Bion, iLight Pushchair, L110

Item code: 02556



### iLight Just Got Lighter!

The new iLight Pushchair, L110 is the next generation of Bion iLight Pushchair EZ.


A new model that successfully combines functionality and aesthetics while also being lighter and more maneuverable.



### Specifications:

|                             |                                  |
|-----------------------------|----------------------------------|
| Seat Width                  | 17.3" / 44cm                     |
| Seat Depth                  | 16.9" / 43cm                     |
| Seat Height                 | 18.9" / 48cm                     |
| Backrest Height (From Seat) | 15.7" / 40cm                     |
| Distance between Armrest    | 17.5" / 44.5cm                   |
| Overall Width (Open)        | 32.7" / 83cm                     |
| Overall Height              | 35.4" / 90cm                     |
| Overall Length (Open)       | 37.8" / 96cm                     |
| Net Weight                  | 10.3kg                           |
| Weight Capacity             | 100kg                            |
| Front Wheel Size            | 6" PVC                           |
| Rear Wheel Size             | 16" X 1.75" PU Tyres, Spokes Rim |
| Frame Material              | Aluminium                        |
| Armrests                    | Half-Length, PU, Fixed           |
| Footrests                   | Composite Footplate, Fixed       |
| Foldable Backrest           | ✓                                |
| Removable Backrest Cover    | X                                |
| Removable Seat Cover        | X                                |
| User Brakes                 | ✓                                |
| Attendant Brakes            | ✓                                |
| Safety Belt                 | ✓                                |
| Anti-Tipper                 | X                                |

### Dimensions:



### Warranty Coverage

- 1 Year off-site warranty on manufacturing defects
- Wear & tear parts are not covered under warranty

