

How to choose a wheelchair?


Choosing a wheelchair for the first time may seem daunting, but fret not, we got your back. Follow this flowchart to find the perfect Bion wheelchair suitable for your needs!

Our wheelchair and pushchair comparison sheet contains the specifications of all our wheelchair & pushchair models for your convenience.



*elevating leg rest is useful for users with leg injuries, the elevating leg rest allows them to rest their leg extended

** reclining feature helps in preventing bedsores and allows user to change body angle without much hassle







Bion Pushchair Comparison Chart

bion

Item No	2002400	2001700	2003800	2000800	2000900
Series	Comfy	Comfy	iLight	iLight	iLight
Item Name	Bion Comfy Pushchair 4G	Bion Comfy Pushchair 3G	Bion iLight Pushchair LI10	Bion iLight Pushchair	Bion iLight Pushchair EZ
Image					
Seat Width	15¾" / 40 cm	15¾" / 40 cm	17.3" / 44 cm	15¾" / 40 cm	17" / 45 cm
Seat Depth	16" / 41 cm	15¾" / 40 cm	16.9" / 43 cm	16" / 41 cm	16¾" / 43 cm
Seat Height	17" / 43 cm	17½" / 44.5 cm	18.9" / 48 cm	18" / 48 cm	18" / 46 cm
Backrest Height from Seat	16¾" / 42.5 cm	16¾" / 41.5 cm	15.7" / 40 cm	16" / 41 cm	15¾" / 39 cm
Distance between Armrests	17½" / 44.5 cm	18" / 45.5 cm	17.5" / 44.5 cm	17¾" / 45 cm	16¾" / 43 cm
Overall Width	21¾" / 55 cm	21½" / 54.5 cm	32.7" / 83 cm	23½" / 60 cm	25" / 65 cm
Overall Height	34¾" / 87 cm	34¾" / 87.5 cm	35.4" / 90 cm	35" / 89 cm	33½" / 85 cm
Overall Length	39¾" / 101 cm	39" / 99 cm	37.8" / 96 cm	39½" / 100 cm	40½" / 103 cm
Net Weight	11.2 kg	10.4 kg	10.3 kg	11.0 kg	10.65 kg
Weight Capacity	100 kg	100 kg	100 kg	100 kg	100 kg
Front Wheel Size	6"	6"	6" PVC	6"	6"
Rear Wheel Size	12" solid tyres, spokes rim	12" solid tyres, mag rim	16" x 1½" PU tyres, spokes rim	12" solid tyres, mag rim	12" solid tyres, mag rim
Frame Material	Aluminium	Aluminium	Aluminium	Aluminium	Aluminium
Armrest	Half-length, PU, fixed	Half-length, padded, fixed	Half-length, PU, fixed	Half-length, PU, fixed	Half-length, padded, fixed
Footrest	Composite footplate, flip-up	Composite footplate, flip-up	Composite footplate, fixed	Composite footplate, fixed	Composite footplate, fixed
Foldable Backrest	Yes	Yes	Yes	Yes	Yes
Removable Backrest Cover	Yes	No	No	Yes	No
Removable Seat Cover	Yes	No	No	Yes	No
User Brakes	Yes	Yes	Yes	Yes	Yes
Attendant Brakes	Yes	Yes	Yes	Yes	Yes
Safety Belt	Yes	Yes	Yes	Yes	Yes
Anti-Tipper	Yes	Optional	No	No	Optional
Tension Backrest	No	No	No	No	No
Quick Release Wheels	No	No	No	No	No
Price	MYR 879.00	MYR 669.00	MYR 499.00	MYR 719.00	MYR 619.00
Version	Ver20190830_1	Ver20190830_1	Ver20231019_1	Ver20220117_1	Ver20190304_1
Warranty	1 year warranty on frame only. Wear and tear parts are not covered under warranty.	1 year warranty on frame only. Wear and tear parts are not covered under warranty.	1 year warranty on frame only. Wear and tear parts are not covered under warranty.	1 year warranty on frame only. Wear and tear parts are not covered under warranty.	1 year warranty on frame only. Wear and tear parts are not covered under warranty.

Bion Pushchair Comparison Chart

bion

Item No	2001000	2002000	2003600	2000300
Series	iLight	iLight	Postur	Traveller
Item Name	Bion iLight Detachable Pushchair 16"	Bion iLight Detachable Pushchair 18"	Bion Postur Pushchair 18"S310	Bion Traveller
Image				
Seat Width	15½" / 39 cm	18" / 46 cm	18.1" / 46 cm	14" / 36 cm
Seat Depth	17" / 43 cm	17" / 43 cm	16.9" / 43 cm	14¾" / 37.5 cm
Seat Height	18.5" / 47 cm	18" / 46 cm	19.7" / 50 cm	18½" / 47 cm
Backrest Height from Seat	15" / 38 cm	16¼" / 42 cm	17.3" / 44 cm	16" / 41 cm
Distance between Armrests	16" / 41 cm	18½" / 47 cm	21.7" / 55 cm	14¾" / 37.5 cm
Overall Width	21" / 53 cm	24" / 61 cm	26" / 66 cm	17" / 43 cm
Overall Height	34" / 86 cm	33½" / 85 cm	35" / 89 cm	34½" / 87.5 cm
Overall Length	40½" / 103 cm	40½" / 103 cm	39.8" / 101 cm	33" / 83.75 cm
Net Weight	13 kg	13 kg	12.1 kg	7.3 kg
Weight Capacity	100 kg	100 kg	100 kg	80 kg
Front Wheel Size	6"	6"	6" PU	5"
Rear Wheel Size	12" solid tyres, mag rim	12" solid tyres, mag rim	16" x 1¾" PU tyres, spokes rim	8" solid tyres, mag rim
Frame Material	Aluminium	Aluminium	Aluminium	Aluminium
Armrest	Half-length, padded, flip-up	Half-length, padded, flip-up	Half-length, PU, flip-up	Half-length, rubber, flip-up
Footrest	Composite footplate, detachable	Composite footplate, detachable	Composite footplate, detachable	Flip-up
Foldable Backrest	Yes	Yes	Yes	Yes
Removable Backrest Cover	Yes	Yes	Yes	No
Removable Seat Cover	Yes	Yes	Yes	No
User Brakes	Yes	Yes	Yes	No
Attendant Brakes	Yes	Yes	Yes	Yes
Safety Belt	Yes	Yes	Yes	Yes
Anti-Tipper	Yes	Yes	Yes	No
Tension Backrest	No	No	No	No
Quick Release Wheels	No	No	No	No
Price	MYR 1,149.00	MYR 1,149.00	MYR 1,149.00	MYR 539.00
Version	Ver20200819_1	Ver20200706_1	Ver20231020_1	Ver20220111_1
Warranty	1 year warranty on frame only. Wear and tear parts are not covered under warranty.	1 year warranty on frame only. Wear and tear parts are not covered under warranty.	1 year warranty on frame only. Wear and tear parts are not covered under warranty.	1 year warranty on frame only. Wear and tear parts are not covered under warranty.