



# 2023/2024 PRODUCT GUIDE



Protecting Clinicians, Patients, and the Practice.

Palmero is now a proud member of HuFriedyGroup





# Don't Take Chances, Take Control with Palmero.

When it comes to infection prevention and control, you can't afford to take risks. For over 60 years, Palmero Healthcare has kept the dental community safe and secure with a complete portfolio of protective products. We are Palmero Healthcare and we don't take chances, we take responsibility.

Palmero Healthcare has been committed to bringing a wide range of infection prevention products such as hospital-level disinfectants, barrier products, and safety products to the dental community.

We take pride in knowing that our legacy will continue to expand and flourish as a proud member of HuFriedyGroup. With quality and reliable dental products, Palmero's goal is to become the recognized brand and leader in infection prevention and safety products.

## Mission:

To provide infection prevention and safety products that protect clinicians, patients, and the practice.

## National and International Memberships:



Proud Member



DENTAL TRADE ALLIANCE



**Certifications:** Certified ISO 9001:2015, Certified System ISO 9001 QMI-SAI Global

## VISIT US ON:



## Palmero Healthcare LLC

120 Goodwin Place

Stratford, CT 06615

[www.palmerohealth.com](http://www.palmerohealth.com)

1-800-344-6424 or 203-377-6424

# TABLE OF CONTENTS

CLEAN

PG. 4

DISINFECT

PG. 10

PROTECT

PG. 13

ORGANIZE

PG. 48

SUPPORT

PG. 56

WHAT'S NEW.....	2
WIPES.....	4
HAND ASEPSIS.....	5
EVACUATION.....	6
SCRUBBERS/BRUSHES.....	8
INSTRUMENT SCRUBBING SYSTEM.....	9
LENS CLEANING.....	9
SURFACE DISINFECTANTS.....	10
HAND ASEPSIS.....	12
BARRIERS.....	13
EYEWEAR & SHIELDS.....	17
EYEWEAR MAINTENANCE & ACCESSORIES.....	37
HEADLIGHTS.....	39
X-RAY APRONS & HANGERS.....	40
HOLDERS.....	48
BUR BLOCKS.....	50
BUR BLOCKS BOX & ORGANIZER.....	53
DISPENSERS/STORAGE.....	54
BOWLS & SPATULAS.....	56
OSHA COMPLIANCE LABELS.....	59
BIB HOLDERS, LAB APRONS, & PATIENT THROWS.....	60
MIRRORS.....	61
OFFICE MAINTENANCE.....	61

## WHAT'S NEW



### ProVision® Secure™ with Strap

- Distortion-free wraparound anti-fog and anti-static polycarbonate lens
- Vent design at top and bottom of foam insulation allows ventilation for comfort while protecting against spray, splashes and splatter
- Detachable, adjustable strap firmly secures eyewear in place

See page 25 for more details



### HVEsolo™ Disposable Evacuation Tips 50/bag

- Our short (2.5" length), lightweight shaft is uniquely designed for solo practitioner use of the high-volume evacuation (HVE) in non-dominant hand
- Reduce aerosol contamination and effectively eliminate fluid and debris from the oral cavity for improved visibility, patient acceptance and clinician comfort

See page 8 for more details



### Dynamic Disposables™ Safety Eyewear

- Our office packs for grey and orange bonding lenses is part of our newest addition to our Dynamic Disposables safety eyewear line
- A lightweight, disposable, cost-effective alternative eye protection that reduces eye splash contamination incidents

See page 31 for more details



### Tablet Barrier Sleeves 100-Pack 8" x 11.5"

- Designed to protect sensitive touch screen tablet devices by minimizing disinfection requirements
- Easy application disposable sleeves accelerate operator turnover
- Adhesive strip securely protects touch screen tablet device
- Compatible with iPad® or tablet

See page 14 for more details



### ProVision® Anti-Fog Solution 8 oz. & 1.7 oz.

- An alcohol-free optical foaming cleaning solution that resist scratches, static, and stains.
- Foam formula yields 5 times the applications compared to a traditional pump or squeeze bottle.
- Compatible with all eyewear, face shields, loupes, and scopes.
- Alcohol-free, non-toxic, and non-flammable

See page 9 & 37 for more details



## Dynamic Disposables™ Snapeez™ Face Shields Full & Half

- The next generation of disposable face shields is here
- Our new light-weight, Dynamic Disposables Snapeez face shields provide an easy-to-use snap-on disposable lens feature
- Meets OSHA bloodborne pathogen guidelines

See page 32-33 for more details



## Protectall™ Upholstery Cleaner Wipes

- Protectall™ Upholstery Cleaner Wipes are premoistened wipes containing an advanced water-based silicone formula
- Protectall wipes safely protects, cleans, and restores leather, vinyl, plastic, and rubber upholstery surfaces like x-ray aprons and dental unit chairs
- Protectall is recommended for periodic use after disinfecting upholstery surface with DisCide® Ultra surface disinfectant

See page 4 & 45 for more details



## Cling Shield® Smiley Face X-Ray Apron

- Cling Shield® X-Ray product line is designed to protect your staff and patients from harmful radiation
- Our PROTECTALL™ x-ray apron offers a built-in thyroid neck collar to optimally protect reproductive organs and thyroid gland

See page 40-41 for more details



## ProVision® New Safety Eyewear

- Palmero is the leading safety eyewear brand providing an impressive range of trendy eyewear including safety, bonding, disposable and magnifying styles for clinicians and patients
- We're excited to introduce our newest safety eyewear line with 13 new designs to enhance protection and comfort for your practice, staff, and patients
- All of our ProVision® safety eyewear meets ANSI Z87.1: Eye and face protection standard with UVA and UVB protection

See page 17-30 for more details



## ProVision® Crystal Clear Anti-Fog Solo Wipes 100-Pack

- As an addition to our Eyewear Maintenance product line from lens station to optical wipes and cleaning solutions, we're now offering convenient individually wrapped anti-fog wipes
- Great for eyewear and goggles

See page 5, 9, & 37 for more details

# CLEAN

## Clean Operatory, Clean Conscience.

At Palmero Healthcare, we don't just think about cleaning products, we obsess over them – from trusted cleaners like DisCide® XRA hand sanitizing wipes to specialized products like Vacuum Shock and Clean for evacuation systems. Because, at the end of the day, a clean operatory comes with a clean conscience.

## WIPES



### Protectall™ Upholstery Cleaner Wipes

- Premoistened wipes containing an advanced water-based silicone formula
- Safely protects, cleans, and restores leather, vinyl, plastic, and rubber upholstery surfaces like x-ray aprons and dental unit chairs
- Recommended for periodic use after disinfecting upholstery surface with DisCide® Ultra surface disinfectant
- Evenly wipe down the surfaces with Protectall™ wipes and allow to dry to touch

---

Ref #	Size	Quantity
3800	10" x 12"	50 count canister

---



### Versi-Wipe™ Dry Towelettes

- Can be used as a dry or moist wipe
- To customize your wipes, simply add soap water solution or your preferred solution
- Includes OSHA-compliant, fill-in blank labels for personalized labeling
- Container/lid does not contain latex

---

Ref #	Size	Quantity
3516	6" x 6 ¾"	160 count canister

---



NEW & IMPROVED

### Refreshing Patient Facial Towelettes

- Improve your patient experience with lightly scented facial wipes
- New calming lavender scent, now formulated with natural oils
- Improved softer towelette sheets
- Container/lid does not contain latex

---

Ref #	Size	Quantity
3515	6" x 6 ½"	120 count canister

---



### ProVision® Crystal Clear Optical Cleaning Wipes

- Convenient wipes premoistened with anti-fog, anti-static lens cleaning, and polishing solution
- Cleans dirt, grease, and debris from eyewear
- Container/lid does not contain latex

---

Ref #	Size	Quantity
3534	6" x 6 ¾"	160 count canister

---

NEW



### ProVision® Crystal Clear Anti-Fog Solo Wipes

- Cleans lenses while providing long-lasting, anti-fog protection
- Safe and effective for use on uncoated lenses including polycarbonate, acrylic, and glass
- Dispenser box contains 100 premoistened individually wrapped towelettes
- Great for eyewear and goggles

---

Ref #	Size	Quantity
3534S	8" x 5"	100 wipes/box

---




### Citrisol™ Orange Solvent Towelettes

- Two-ply premoistened orange solvent towelettes with lanolin
- One side is abrasive to loosen dirt and grease with a reverse smooth side to wipe clean
- Quickly removes cements, impression materials, adhesives, grease, gum, tar, and more

---

Ref #	Size	Quantity
69	9 ½" x 12"	70 count canister

---



### DisCide® XRA Hand Sanitizing Wipes

- Kills 99.9% of disease-causing microbes
- Contains aloe and vitamin E to soothe and moisturize the skin
- Fresh wintergreen fragrance
- Dries in seconds to ensure staff compliance with hand asepsis guidelines

---

Ref #	Size	Quantity
3500	6" x 8"	160 count canister

---

## HAND ASEPSIS

### DisCide® Effect Professional Hand Asepsis Soap

- A professional hand asepsis soap with aloe vera, vitamin E, and fresh almond scent to promote healthy skin and staff compliance
- Quickly and easily removes all soils and greases without irritation to even the most delicate hands
- Pearlescent viscous formula produces a dense creamy lather and contains no alkali or other irritants
- Container/lid does not contain latex

---

Ref #	Description	Quantity
3543	DisCide Effect Pump	Quart Bottle Pump Dispenser
3542	DisCide Effect Refill	Gallon Refill

---



## OFFICE MAINTENANCE

### Once-A-Day One Second Spray

- For use in air turbine hand pieces, prophy angles, and contra angles
- Also ideal as a lubricant or cleaner for hand pieces
- Contains hexane, non-corrosive

---

Ref #	Description
702	8.8 oz. aerosol can with extension tube

---






## EVACUATION

# From Screen to Clean

## The Two-Step System Too Simple to Ignore

Keeping your evacuation line clean is an important part of keeping your HVE and saliva ejector running smoothly. Vacuum Shock and Vacuum Clean are an easy two-step system that keeps high/low speed lines clean and fresh. When used in conjunction, the time-release formulation ensures your line is protected around the clock. Our Evacuation Product Line features neutral formulations with a pH of 7 that follows best management practices (BMPs) identified by the EPA to help ensure the efficiency of amalgam separators. These products do not contain bleach, chlorine, iodine, peroxide, or any oxidizers. Don't let a compromised evacuation system make problems for your practice, use the one-two punch that keeps evacuation lines happy.



### E-VAC System Kit

- Convenient, all-in-one Vacuum Shock and Vacuum Clean Evacuation System Maintenance Kit
- Kit includes 6-tablets/jar of Vacuum Shock and 45-tablets/jar of Vacuum Clean

---

Ref #	Quantity
3547S	2 jars

---



### Vacuum Shock™

- First stage in Palmero 2-Step Evacuation System Maintenance Kit
- Time-released tablet; use 1 tablet monthly
- Removes accumulations and keeps vacuum lines clean to restore vacuum pressure and extends pump life
- Compatible with amalgam separator systems
- Neutral (7 pH) formulation follows EPA Guidelines for waterline management
- Does not contain bleach, chlorine, iodine, peroxide, or any oxidizers

---

Ref #	Quantity
3546	6 tablets/jar

---



### Vacuum Clean™

- Second stage in Palmero 2-Step Evacuation System Maintenance Kit
- Time-released tablet; use 1 tablet weekly
- Continuously cleans evacuation system lines
- Restores pressure, and extends pump life
- Keeps vacuum lines functioning at peak efficiency
- Can be used as a cost effective detergent to clean high/low speed lines from tip to trap by dissolving 5 tablets in the Pour and Clean pre-mix bottle (3547P) for premixing concentrate
- Compatible with amalgam separator systems
- Neutral (7 pH) formulation follows EPA Guidelines for waterline management
- Does not contain bleach, chlorine, iodine, peroxide, or any oxidizers

---

Ref #	Quantity
3547	45 tablets/jar

---





## Shock & Clean Starter Kit

- The one-two punch for evacuation system maintenance in one convenient kit
- Kit includes 1 jar of Vacuum Shock, 1 jar of Vacuum Clean and 1 Pour & Clean Bottle for premixing
- Clean tablets can be used as a cost effective detergent to clean high/low speed lines from tip to trap by dissolving 5 tablets in the Pour and Clean pre-mix bottle for premixing concentrate

---

Ref #	Quantity
3546K	1 kit

---



## Pour & Clean Bottle

- Measuring container to prepare rinse treatment for the HVE and low-speed lines
- Convenient pre-measured compartments with easy markers for easy mixing and pouring of solution

---

Ref #	Quantity
3547P	16 oz. container

---




## Shock & Clean Convenient Operatory & Trial Kits

- 2-step system keeps vacuum lines functioning at peak efficiency
- Cleans and maintains your evacuation system
- No fuss, easy-to-use operatory kits premeasured for 16-day schedule
- Convenient packages for operatory storage. Kit includes 6 tubes, each containing 16 Vacuum Clean tablets and 1 Vacuum Shock tablet. (also available in a single trial tube)

---

Ref #	Quantity
35460	6 tubes kit
3546TK	1 trial tube

---



## Evac Screens

- Helps avoid contact with potentially hazardous material
- Keeps debris from lodging in evacuation lines
- Disposable screens color-coded by size and fit

---


Ref #	Size/Description	Quantity
1805	Blue, 2 1/8" A-DEC with plastic canisters	144 screens/box
1807	Green, 1 7/8" Belmont, Biotec, Chayes, Dentech, DentalEZ, Pelton/Crane, Proma and others	144 screens/box

---



## HVESOLO™ Disposable Evacuation Tips

Uniquely designed for clinicians performing power instrumentation (air polishing/ultrasonic) procedures without an assistant



- Prevent disease transmission with HVE as it's been shown to reduce contamination arising from the operative site by more than 90%\*
- 3 vent design and smooth round edges gently retract lip/cheek and protect soft tissues to maximize patient comfort
- Proper venting and a wide scoop tip design maximize suction to clear the field of fluid and debris efficiently
- Reduced shaft length (50%) lessens tension/draft from the hose and allowing a more neutral wrist positioning for clinician comfort
- Beveled interface smoothly inserts into standard HVE adaptor
- Does not contain latex

\* Aerosols and splatter in dentistry HARREL, STEPHEN K. et al. JADA Volume 135, Issue 4, 429 -437

Ref #	Size	Quantity
1808S	2.5" Length	50/bag



## Saliva Ejector Screens

- Universal fit disposable saliva ejector screens
- Keeps debris from being aspirated and lodged in the lines

Ref #	Quantity
1602	25 pack
1602A	50 pack
1602C	75 pack
1602B	100 pack
1602D	250 pack



## Ultra Super Sonic

- Technologically advanced ultrasonic and instrument pre-cleaning tablet
- Engineered to achieve maximum instrument cleaning and proper removal of bacteria, debris, and bioburden
- No mixing required, tablet is self-dissolving
- Contains cleaning properties for efficient ultrasonic cleaning or presoak
- Non-toxic and non-corrosive

Ref #	Quantity
3559	30 tablets/jar



## SCRUBBERS/BRUSHES

- Scrubbers help effectively remove debris from instruments before cleaning and sterilization
- Hand scrub brushes assist staff and clinicians with hand washing and surgical hand scrubbing
- Stiff bristles and sturdy ergonomic handles
- Autoclavable, high-impact plastic



Ref #	Size/Description
1868	Hand Scrub Brush, 3 1/2"L x 1 1/2"W, with 1/2" bristle length, white
1868H	Super Scrubber, 4"L x 2"W, with 3/4" bristle length, white

## INSTRUMENT SCRUBBING SYSTEM



- Mechanical scrubbing system made of autoclavable material
- Built-in suction base to secure the tub in place
- Recommend use with an approved solution such as DisCide® ULTRA disinfectant

Ref #	Quantity
1952	Scrubbing Complete Kit
1952R	Replacement brushes (set of 8)

## LENS CLEANING

### ProVision® Crystal Clear Anti-Fog Solo Wipes



- Cleans lenses while providing long-lasting, anti-fog protection
- Safe and effective for use on uncoated lenses including polycarbonate, acrylic, and glass
- Dispenser box contains 100 premoistened individually wrapped towelettes
- Great for eyewear and goggles

Ref #	Size	Quantity
3534S	8" x 5"	100 wipes/box

### ProVision® Crystal Clear Optical Cleaning Wipes



- Convenient wipes premoistened with anti-fog, anti-static lens cleaning, and polishing solution
- Cleans dirt, grease, and debris from eyewear
- Container/lid does not contain latex

Ref #	Size	Quantity
3534	6" x 6 3/4"	160 count canister

### ProVision® Anti-Fog Solution



- Compatible with all eyewear, face shields, mirrors, loupes, and scopes
- Resists scratches, static, and stains
- Foam formula yields 5 times the applications compared to a traditional pump or squeeze bottle
- Alcohol-free, non-toxic, and non-flammable

Ref #	Size	Quantity
3533	8 oz. (237 ml)	1 bottle
3536	1.7 oz. (51 ml)	1 bottle

### ProVision® Lens Cleaning Station



- Convenient eyewear lens cleaning station with anti-fog and anti-static additives
- Designed for countertop placement or wall mounting with hanging hardware
- Station includes 8 oz. pump spray bottle of optical cleaning solution and 600 dry wipes

Ref #	Size	Quantity
3535	7 3/4" x 8" x 5"	1 station

# DISINFECT



## Accelerate Operator Turnover with DisCide® Ultra Disinfectant

With intermediate level DisCide Ultra, our surface disinfectant and decontaminant cleaner, you can accelerate the way you disinfect clinical settings. In just 1 minute, DisCide Ultra is effective against TB, viruses (HBV, HCV, RSV, HIV-1), bacteria (MRSA, VRE, E. coli), and fungi (C. albicans). Specifically formulated to resist evaporation, it's non-corrosive, non-staining, and leaves no residue. Plus, it meets OSHA's blood-borne pathogen standard. Pre-mixed for the ultimate in convenience, DisCide Ultra is available in identical formulations in both spray and towelettes formats to ensure efficacy when liquid and towelette are used in tandem. DisCide to be sure with fast acting DisCide Ultra.

### SURFACE DISINFECTANTS

#### DisCide® ULTRA

- Hospital-level, 1-step, ready-to-use quaternary ammonium, intermediate high-level alcohol-based disinfectant
- Proven to kill deadly pathogens in 1 minute or less for accelerated operator turnover
- Proprietary formulation is non-corrosive, non-staining, and leaves no residue with a pleasant herbal scent
- Non-woven, lint-free, and durable towelettes are premoistened and ready-to-use
- Meets CDC, OSHA, and infection control guidelines. EPA and FDA registered
- Container/lid does not contain latex
- Available in towelettes, spray, or bulk gallon refill



Ref #	Description/Size	Quantity
10DIS	DisCide Towelettes, 10" x 10"	60 count canister
60DIS	DisCide Towelettes, 6" x 6 3/4"	160 count canister
3564	DisCide Spray	4 oz. bottle
3565Q	DisCide Spray	Quart Sprayer
3565G	DisCide Bulk Refill	Gallon Refill

#### Accessories

3567	Empty Quart Bottle with Sprayer and DisCide Ultra label
3567T	Quart Spray Trigger
3519	Unlabeled Empty Quart Bottle with Sprayer

## DISINFECTANT CHART

	MICROORGANISM TESTED	DISCIDE ULTRA TOWELETTES	DISCIDE ULTRA LIQUIDS
	EPA REGISTRATION NUMBER	10492-4	10492-5
GRAM-POSITIVE BACTERIA	MRSA-METHICILLIN-RESISTANT STAPHYLOCOCCUS AUREUS	1 MIN	1 MIN
	STAPH-STAPHYLOCOCCUS AUREUS	1 MIN	1 MIN
	VRE-VANCOMYCIN-RESISTANT ENTEROCOCCUS FAECALIS	1 MIN	
	VRE-VANCOMYCIN-RESISTANT ENTEROCOCCUS FAECIUM		1 MIN
GRAM-NEGATIVE BACTERIA	ANCINETOBACTER BAUMANNII		1 MIN
	BURKHOLDERIA CEPACIA	1 MIN	1 MIN
	CAMPYLOBACTER JEJUNI	1 MIN	1 MIN
	CARBAPENEM-RESISTANT E COLI		1 MIN
	E COLI-ESCHERICHIA COLI	1 MIN	1 MIN
	ESCHERICHIA COLI O157:H7	1 MIN	1 MIN
	KLEBSIELLA PNEUMONIAE		1 MIN
	PROTEUS MIRABILIS	1 MIN	1 MIN
	PSEUDOMONAS AERUGINOSA	1 MIN	1 MIN
	SALMONELLA ENTERICA	1 MIN	1 MIN
ENVELOPED VIRUSES	AVIAN INFLUENZA A VIRUS-H3N2	30 SEC	30 SEC
	H1N1-PANDEMIC 2003 INFLUENZA A VIRUS (FORMERLY SWINE FLU)	30 SEC	30 SEC
	HBV-HEPATITIS B VIRUS	1 MIN	1 MIN
	HCV-HEPATITIS C VIRUS	1 MIN	1 MIN
	HERPES SIMPLEX VIRUS TYPE 2	1 MIN	1 MIN
	HIV-1 (AIDS VIRUS)	1 MIN	1 MIN
	HUMAN CORONAVIRUS	30 SEC	30 SEC
	INFLUENZA A VIRUS (HONG KONG)	30 SEC	30 SEC
	RSV-RESPIRATORY SYNCYTIAL VIRUS	1 MIN	1 MIN
	VACCINIA VIRUS	1 MIN	1 MIN
NON-ENVELOPED VIRUSES	ADENOVIRUS TYPE 2	1 MIN	1 MIN
FUNGI	CANDIDA ALBICANS	1 MIN	1 MIN
MYCOBACTERIA	TB-MYCOBACTERIUM BOVIS (TUBERCULOSIS)	1 MIN	1 MIN
	TOTAL MICROORGANISMS TESTED	23	26



## DisAseptic® XRQ

- Cost-effective quaternary ammonium-based detergent, disinfectant, cleaner, and deodorizer
- One-step, ready-to-use no mixing or measuring
- Laboratory-proven to kill microorganisms on hard, inanimate, non-porous surfaces
- Reduces the risk of infection and cross-contamination
- EPA registered
- Container/lid does not contain latex



Ref #	Description	Quantity
3503	DisAseptic XRQ Spray	Quart Sprayer
3504	DisAseptic XRQ Bulk Refill	Gallon Refill
3532	DisAseptic XRQ Bulk Refill	5 Gallon Refill (two 2.5 gallon with spigot)

### Accessories

3508	Empty Quart Bottle with Sprayer and DisAseptic XRQ label
3567T	Quart Spray Trigger
3519	Unlabeled Empty Quart Bottle with Sprayer

## HAND ASEPSIS

### DisCide® Effect Professional Hand Asepsis Soap

- A professional hand asepsis soap with aloe vera, vitamin E, and fresh almond scent to promote healthy skin and staff compliance
- Quickly and easily removes all soils and greases without irritation to even the most delicate hands
- Pearlescent viscous formula produces a dense creamy lather and contains no alkali or other irritants
- Container/lid does not contain latex



Ref #	Description	Quantity
3543	DisCide Effect Pump	Quart Bottle Pump Dispenser
3542	DisCide Effect Refill	Gallon Refill

### DisCide® XRA Hand Sanitizing Wipes

- Kills 99.9% of disease-causing microbes
- Contains aloe and vitamin E to soothe and moisturize the skin
- Fresh wintergreen fragrance
- Dries in seconds to ensure staff compliance with hand asepsis guidelines



Ref #	Size	Quantity
3500	6" x 8"	160 count canister

PROTECT

## Stand and Deliver


The Palmero Healthcare portfolio of protective products has the confidence to stand up to your convictions. From eyesight protection to radiation mitigation wear, we'll stand and deliver protection with the Palmero Healthcare line of safety products.

## BARRIERS

### Films



- Polyethylene film with a low-tack adhesive that adheres to most operatory surfaces, leaves no residue, does not contain latex
- Accelerates operatory turnover, makes cleanup between patients fast and easy
- Protects equipment from harmful chemicals by minimizing disinfection requirements and staff and patients from cross-contamination
- Available in a variety of sizes and colors, packed in a corrugated film box that serves as a dispenser



Ref #	Size	Description	Quantity
1803B	4" x 6"	Blue	1200 sheets/roll
1803C	4" x 6"	Clear	1200 sheets/roll
1803M	4" x 6"	Mauve	1200 sheets/roll
1803CP	4" x 6"	Blue Bear Print	1200 sheets/roll
1865C	2 1/2" x 6"	Clear	1200 sheets/roll
1866C	6" x 9"	Clear	800 sheets/roll
1867C	12" x 14"	Clear	514 sheets/roll

### Accessories

Plastic Durable Dispensers (film not included), See ORGANIZE, Page 54



### Keyboard Film

- Perforated, plastic sheets with a light adhesive that covers virtually any keyboard
- Packed in a corrugated film box that serves as a dispenser
- Prevents cross-contamination and leaves behind no residue, does not contain latex

Ref #	Size	Description	Quantity
1802	9" x 19"	Keyboard film	800 sheets/roll

## Protector Sleeves

- Convenient disposable plastic sleeves to protect your instruments/equipment, saves time between patients
- Prevents cross-contamination, does not contain latex
- Easy to apply and remove, quick clean-up after every procedure
- Available in a variety of styles and sizes



### Tablet Barrier Sleeves

- Designed to protect sensitive touch screen tablet devices by minimizing disinfection requirements
- Easy application disposable sleeves accelerate operator turnover
- Adhesive strip securely protects touch screen tablet device
- Compatible with iPad® or tablet

Ref #	Size	Quantity
1804	8" x 11 1/2"	100/box



### T-Light Handle Protector Sleeves

- Fits over all "T" style handles
- Reduces the risk of cross-contamination from repeated touching of lighting

Ref #	Size	Quantity
1914	5 3/4" x 4 1/4"	500/box



### Tubing Protectors

- Keeps tubing clean and safe from contaminants, can also be used on armrests
- Roll can be cut to desired length, comes in 2 width sizes

Ref #	Size	Quantity
1912	2 1/4" wide	1200/roll
1912B	3" wide	1200/roll




### Digital X-Ray Sensor Protectors

- Disposable plastic covers for digital x-ray sensors
- Universal fit for a variety of manufacturer's equipment

Ref #	Size	Film Size	Quantity
1915A	2 1/2" x 10"	2	500/box
1970	1 3/8" x 8"	0/1	500/box
1971	1 5/8" x 8"	2	500/box





## Chair Protectors

- Protects dental units from long-term use of cleaners and disinfectants
- Save time and energy with our easy-to-use and universal-fitting chair protectors



Ref #	Description/Size	Quantity
1901	Half Chair, 14" x 13 ¾" x 23"	225/roll
1902	Full Chair, 39" x 84"	175/roll



## Headrest Protectors

- Reduces cross-contamination and provides cleanliness between patients

Ref #	Size	Quantity
1918A	11 ½" x 3" x 8 ¾"	250/box
1918B	14 ½" x 3" x 9"	250/box



## X-Ray Head Protectors

- Prevents x-ray equipment from chemicals and cross-contamination

Ref #	Size	Quantity
1910	15" x 26"	225/roll
1910B	23" x 28"	225/roll
1920	10" x 6" x 24"	500/box, for use with Nomad x-ray



## Bite Block Protectors

- Clear, convenient plastic barriers for use with panoramic x-ray and Siemens machine bite blocks

Ref #	Size	Quantity
1850	Pano 1" x 2"	1000/box
1850S	Siemens 1 ½" x 3"	500/box



## Delivery System Protectors

- Covers and protects your delivery system from chemicals, cross-contamination, and aerosols

Ref #	Size	Quantity
1911A	19" x 19"	500/box



## Air/Water Syringe Protectors

- Protects non-autoclavable parts from cross-contamination

Ref #	Size	Quantity
1915	2 ½" x 10"	500/box



## Syringe Tip Covers

- Disposable, rotatable plastic syringe tip cover
- Designed to easily slide over all air/water syringes

Ref #	Description	Quantity
1950	White	100/bag
1949	White	250/bag

## Eyewear Collections at a Glance

Traditional Safety Eyewear Meets ANSI Z87.1 standard. Does not contain latex						Disposables Does not contain latex	
Wraparound	Goggle	Magnification	Bifocals	Bonding	Laser	Eyewear	Face Shields
<ul style="list-style-type: none"> <li>• Ergonomic lightweight</li> <li>• Full &amp; minis</li> <li>• Ultra anti-fog</li> </ul>	<ul style="list-style-type: none"> <li>• Wear over prescription glasses</li> <li>• Side shields</li> </ul>	<ul style="list-style-type: none"> <li>• Full lens magnification</li> </ul>	<ul style="list-style-type: none"> <li>• Up to +3.0 diopter</li> </ul>	<ul style="list-style-type: none"> <li>• Red and orange lens</li> <li>• Full &amp; minis</li> </ul>	<ul style="list-style-type: none"> <li>• Yellow and green lens</li> </ul>	<ul style="list-style-type: none"> <li>• Snap-on feature</li> <li>• One frame, multicolored lens options</li> </ul>	<ul style="list-style-type: none"> <li>• Full &amp; half sizes</li> <li>• Foam elastic band or snap-on style</li> </ul>

## EYEWEAR & SHIELDS

### The Insight to Safeguard Your Eyesight

Palmero offers an impressive range of trendy eyewear including safety, bonding, disposable, and magnifying styles. Our safety eyewear meets the ANSI Z87.1 eye and face protection standards with UVA, UVB and blue light protection. We also offer accessories such as clip-on lights, magnifiers, and lens cleaning station and maintenance products.


### Implement an Eye Conservation Program for your Practice

Select safety eyewear that serves dental professionals and patients in your practice on multiple levels to: eliminate blue light exposure; protect against spray, splash and splatter for infection control; protect from high impact flying projectiles that can cause permanent damage to the eye and are comfortable enough to insure compliance.

## A Guide To Eyewear Parts And Terms

**NEW**


Designed to meet the ANSI Z87.1 standard



**ProVision<sup>®</sup> Infinity**, See Page 21

### The Four C's of Safety Eyewear

**Certified + Clarity + Comfort = Compliance**



**CERTIFIED** Meets industry standards and guidelines for lens, cornea and retina protection i.e., ANSI Z87.1, OSHA and CDC. Wraparound design for full coverage. Manufactured from high impact resistant materials to protect from flying debris and projectiles. Specialty bonding lens filters out damaging Blue Light 400 nm-500 nm levels to the retina.



**CLARITY** Excellent optics for visual acuity provided by anti-fog coating, scratch resistant technologies, and superior lens designs.



**COMFORT** Broad range of frame, temple and nose bridge styles to provide an individualized fit for all day security and optimal comfort.



**COMPLIANCE** Willingness and discipline to use safety eyewear for patients and clinicians, in every procedure, every day.

## Wraparound Collection



### ProVision® Econo Wraps™

- Sleek, featherweight polycarbonate wraps, universal fit
- Fog-free, scratch-resistant
- Available in 3 colors
- Blocks up to 95% of UVA and UVB rays
- Meets the ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3601G	Grey frame/grey lens	full size
3601C	Clear frame/clear lens	full size
3600	Clear frame/red bonding lens	full size

#### Bulk, 10-Pack Dispenser

3601GB	10-pack of Ref 3601G	full size
3601CB	10-pack of Ref 3601C	full size

NEW



ProVision®  
Mini Econo Wraps  
(25% smaller)



ProVision®  
Econo Wraps™

### ProVision® Mini Econo Wraps

- Same great features and benefits as full size Econo Wraps
- 25% smaller, ideal for kids and adults with smaller/narrow facial structure
- Available in 3 colors

Ref #	Description	Size/Facial Fit
3607G	Grey frame/grey lens	child/youth
3607C	Clear frame/clear lens	child/youth
3607B	Red frame/red bonding lens	child/youth

NEW



### ProVision® Rainbow™ Mini

- Flexible temples in 6 colors make it easy to have your smaller or pediatric patients wear their safety glasses
- Reduce the effect of harsh operatory lights with scratch-resistant, 100% polycarbonate wraparound grey lens
- The Mini is specifically designed to fit smaller faces
- Blocks up to 99.9% of UVA and UVB rays
- Meets the ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3607GB	12-Pack of 6 temple colors (blue, orange, green, purple, yellow & pink) with grey lenses	child/youth



Clear/grey lens option



Ultra lightweight frame



Soft ventilated nose pad

### ProVision® Allure™

- Virtually unbreakable Pink or Lavender ultra lightweight frame (24 grams)
- Soft ventilated nose pad for a secure, comfortable fit
- Curved lens for wraparound splatter protection
- Fits small/narrow and medium face shapes
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref #	Description	Size/Facial Fit
3604PC	Pink frame/clear lens	small/narrow/medium
3604LC	Lavender frame/clear lens	small/narrow/medium
3604LG	Lavender frame/grey lens	small/narrow/medium

### Did You Know?

Palmero Provision Eyewear provides full eye protection as the polycarbonate material filters out 99.99% of the harmful UV rays to the front of the eye, protecting the cornea and the lens.

- Light is electromagnetic radiation within a certain portion of the electromagnetic spectrum.
- Visible light is usually defined as having wavelengths in the range of 400–700 nanometers (nm), between the infrared (with longer wavelengths) and the ultraviolet UV [with shorter wavelengths].
- UV rays are unseen from the human eye but can cause damage due to long term exposure.
- UV light affects the front of the eye

#### Electromagnetic Spectrum of Wavelengths



#### Dangers of Light to the Eye



UV light affects the front of the eye.  
Blue light affects the back of the eye.



### ProVision® Air™

- Ergonomically designed temples fit the contour of your head for a comfortable, secure, non-slip fit
- Flexible, anti-static, lightweight frame (15 grams)
- Extended lens coverage for great protection and style
- Distortion-free, wraparound polycarbonate lens offers superior peripheral vision and protection
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3608C	Clear frame/clear lens	small/narrow/medium
3608G	Grey frame/grey lens	small/narrow/medium



### ProVision® Element™

- Ergonomically designed temples fit the contour of your head for a comfortable, secure, non-slip fit
- Flexible, lightweight frame (20 grams)
- Extended lens coverage for great protection and style
- Distortion-free wraparound polycarbonate lens offers superior peripheral vision and protection
- Wider frame design allows for broader facial fit
- Anti-static and anti-fog lens
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3606C	Clear frame/clear lens	medium/large
3606G	Grey frame/grey lens	medium/large



Detachable & replaceable lens



Soft rubber ribbed nose bridge



Soft rubber tips



Telescopic & ratcheting temples



Top vented silicone frame



### ProVision® InFINITY™

- Upscale black frame and clear wraparound lens with anti-fog and anti-static features
- Top vented silicone brow guard provides added protection and venting to reduce fogging
- Enhanced with pantoscopic tilt feature provided with telescopic and ratcheting temple arms allowing for customizable fit and optimal comfort
- Soft micro ribbed nose bridge grip for secure, non-slip fit
- Fits broader facial structures
- Detachable and replacement lens in clear and bonding
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard


Ref #	Description	Size/Facial Fit
3613C	Black frame/clear lens	universal
3613R	Clear replacement lens	universal
3613BR	Red bonding replacement lens	universal

### Did You Know?


- Lens tilt improves the way a lens works and contributes to how well a pair of glasses fits. Tilt is dependent on the interaction of the heights of the ears and bridge of the nose.
- Pantoscopic tilt refers to the frame alignment in the up and down position of the frame.
- Palmero's ProVision InFINITY design addresses pantoscopic tilt with adjustable telescopic and ratcheting temple arms for customizable fit and optimal comfort.

#### Pantoscopic Tilt


The angle of the lenses as they sit on a clinician's face



Natural Head Position



Ear Lower





Ear Higher



### ProVision® Clarity™

- Lightweight, frameless, flexible, translucent wrap (24 grams)
- Firm nose pad and soft, rubberized temple tips for comfort and grip
- Curved lens for wraparound splatter protection
- Enhanced with **Ultra Anti-Fog Technology** to provide incredibly clear and uninterrupted vision
- Fits small/narrow and medium face shape
- Blocks up to 99.9% UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref #	Description	Size/Facial Fit
3605L	Clear frame/lavender tips/clear lens	small/narrow/medium
3605B	Clear frame/navy tips/clear lens	small/narrow/medium



### ProVision® Econo Tilts™

- Lightweight, polycarbonate frame with tilting arms ensures a great contoured fit
- Firm nose pad and soft, rubberized temple tips guarantee comfort and grip
- Curved lens provides wraparound splatter protection
- Anti-fog polycarbonate lens
- Available in 2 colors (Rose Quartz and Blue)
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3614CB	Clear frame/blue tips/clear lens	universal
3614CQ	Clear frame/rose quartz tips/clear lens	universal



### ProVision® Contour Wraps™

- Lightweight polycarbonate frame with telescoping arms ensure a great contoured fit
- Available in 3 colors
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3550	Ice blue frame/black tips/ice blue lens	universal
3552	Black frame/grey lens	universal
3553	Clear frame/black tips/clear lens	universal



### Did You Know?

Palmero ProVision Eyewear provides full eye protection with wraparound eyewear and goggle designs of polycarbonate that meet ANSI Z87.1 standards for high impact resistance to protect from flying debris and projectiles.

- According to the CDC, approximately 2,000 eye injuries occur every day at workplaces in the U.S.
- OSHA currently has adopted ANSI Z87.1 standards as compliant with regulatory requirements for eye protection standards.
- ANSI/ISEA Z87.1 protects the human face in these environments while also acting as an up-to-date guideline by acknowledging any significant changes in the eye and face protector industry.





Short temples for secure & broader facial fit



Integrated, comfortable nose bridge



Soft temple tips

NEW



### ProVision® Flexiwrap™

- Lightweight, short spring-flex temple design with extreme flexibility to fit any facial structure
- Provide secure fit for **patients in the supine position**
- Soft touch temple and soft nose bridge for a comfortable fit
- Full wraparound protection
- Blocks up to 99.9% UVA and UVB rays
- Meets ANSI Z87.1 standard

NEW



Ref #	Description	Size/Facial Fit
3609GG	Black frame/grey tips/grey lens	universal
3609GC	Clear frame/grey tips/clear lens	universal
3609OG	Black frame/orange tips/grey lens	universal
3609OC	Clear frame/orange tips/clear lens	universal

NEW



NEW



### ProVision® Bad Dogs™

- Designed to provide an uninterrupted peripheral view and ultimate protection
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref #	Description	Size/Facial Fit
3620	Clear frame/clear lens	universal
3621	Clear frame/indoor-outdoor lens	universal
3622	Grey frame/grey lens	universal
3623	Yellow frame/yellow lens	universal





## ProVision® Sphere-X Wraps™

- Dual-lens construction allows for the best possible peripheral vision, adjustable arms to 4 different lengths
- Nose pieces do not contain latex
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3707	Black frame/clear lens	full size
3708	Platinum frame/clear lens	full size



Telescoping arms



Detachable & replaceable lens



Integrated, comfortable nose bridge



## ProVision® See-Breez™

- Ultra-lightweight frames with sculpted dual-lenses and telescoping arms for a comfortable fit
- Great for children or adults with smaller faces
- Available in 4 colors and replacement lens
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3560BL	Black frame/clear lens	universal
3560B	Blue frame/clear lens	universal
3561BL	Black frame/grey lens	universal
3560PL	Platinum frame/clear lens	universal

### Replacement Lenses

3560R	Clear lens	universal
3561R	Grey lens	universal





### ProVision® Tech Specs™ Safety

- Lightweight design that provides excellent protection and comfort
- Contoured frame style hugs the face with rubberized (non-latex) attachments on the nosepiece and ear rests
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref #	Description	Size/Facial Fit
3710C	Grey frame/clear lens	universal
3711C	Blue frame/clear lens	universal
3710G	Grey frame/grey lens	universal
3712	Grey frame/bonding UV protective lens	universal



### ProVision® Chic™

- Designed to fit smaller, narrow facial structure with colorful rubber tipped temples and soft rubber nose bridge
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref #	Description	Size/Facial Fit
3760CP	Clear frame/pink tips/clear lens	universal
3760CW	Clear frame/white tips/clear lens	universal
3760CG	Clear frame/green tips/clear lens	universal



### ProVision® Secure™ with Strap

- Distortion-free wraparound with anti-fog and anti-static polycarbonate lens offers superior peripheral vision and protection
- Vent design at top and bottom of foam insulation allows ventilation for comfort while protecting against spray, splashes and splatter
- Detachable, adjustable strap firmly secures eyewear in place to minimize slippage and provide all day comfort
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Vent design at top and bottom

Foam insulation for comfort

Detachable and adjustable strap



NEW

Ref #	Description	Size/Facial Fit
3630GC	Black frame/green tips/clear lens	universal

# Goggles Collection

Full coverage design to accommodate users with or without prescription eyewear



Adjustable temples to fit broader facial



Fully integrated top, bottom & side shields protection




Built-in comfortable nose bridge

## ProVision® Eyesaver Sleeks™

- Sleek, ultra-light, economical goggles
- Adjustable temples allow for broader facial fit
- Designed to fit over prescription eyewear
- All around integrated side shields protection
- Scratch-resistant
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3SLK	Grey frame/grey lens	universal
16SLK	Red bonding frame/lens	universal
17SLK	Clear frame/clear lens	universal
18SLK	Green frame/green lens	universal



All goggles in this collection fit over prescription glasses



## ProVision® Eyesavers™

- Ultra-light, economical, durable no-frills eye protection
- Hinged temples allow for convenient storage when not in use
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3S-PHC	Grey	universal
17S	Clear	universal
18S	Green	universal
19S	Yellow	universal

**Bulk Dispenser**  
1968 10-pack of Ref 17S, hardware not included



### ProVision® End-Fog

- Fog-free, scratch-resistant
- Extra anti-fog coating sets End-Fog apart from other economy eyewear
- Secure snug fit and optically correct, crystal-clear lenses for those special procedures
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref #	Description	Size/Facial Fit
3556B	Blue frame/clear lens	universal
3556C	Clear frame/clear lens	universal
3556R	Red frame/clear lens	universal
3556P	Purple frame/clear lens	universal

**Did You Know?**

Palmero ProVision eyewear line provides anti-fog coatings and solutions to meet CDC guidelines for maintaining clarity of vision as stated in the "Guidance for the Selection and Use of Personal Protective Equipment (PPE) in Healthcare Settings."

"Goggles provide barrier protection for the eyes; personal prescription lenses do not provide optimal eye protection and should not be used as a substitute for goggles. Goggles should fit snugly over and around the eyes or personal prescription lenses. Goggles with antifog features will help maintain clarity of vision."

Source: US Centers for Disease Control and Prevention. Personal Protective Equipment. <https://www.cdc.gov/infectioncontrol/pdf/strive/PPE103-508.pdf>



### ProVision® Overshield™ Eyewear

- Designed to wear over prescription eyewear with telescoping and ratcheting arms to ensure a great fit
- One-piece polycarbonate lens with integral side shields wraps around for a wide field of vision and complete protection
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref #	Description	Size/Facial Fit
3602	Black frame/clear lens	universal
3602BOND	Black frame/orange lens	universal



### ProVision® Oasis™

- Attractive over-the-glass style fits over most prescription eyewear with built-in side shield
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3570BOND	Black frame/reddish-orange lens	universal

## Bifocals Collection



Integrated comfortable nose bridge



Soft rubber temple tips



Four available diopters



### ProVision® EZ-Focals™

- Lightweight, super comfortable wraparound frame (24 grams)
- Contoured temples with soft tips for added comfort, grip, and a secure fit
- Soft, self-fitting nose bridge, minimizes slippage
- Available in +1.5, +2.0, +2.5, +3.0 diopters
- Clear anti-fog and anti-static lens
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Diopter	Size/Facial Fit
3770B	Black frame/clear lens	+1.5	universal
3770C	Black frame/clear lens	+2.0	universal
3770D	Black frame/clear lens	+2.5	universal
3770E	Black frame/clear lens	+3.0	universal



### Provision® Bifocal

- Fog-free, dual curvature mono-lens with integrated side shields
- Fully adjustable, lens pivots with telescopic arms
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Diopter	Size/Facial Fit
3701A	Black frame/clear lens	+1.0	universal
3701B	Black frame/clear lens	+1.5	universal
3701C	Black frame/clear lens	+2.0	universal
3701D	Black frame/clear lens	+2.5	universal
3701E	Black frame/clear lens	+3.0	universal



## ProVision® Kool Daddy™ Bifocal

- Fog-free, scratch-resistant
- Sleek, lightweight wraparound design
- Includes a slip-on neoprene protective sleeve with heavy-duty belt clip
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref #	Description	Diopter	Size/Facial Fit
3740A	Black frame/clear lens	+1.0	universal
3740B	Black frame/clear lens	+1.5	universal
3740C	Black frame/clear lens	+2.0	universal
3740D	Black frame/clear lens	+2.5	universal
3740E	Black frame/clear lens	+3.0	universal
3745B	Black frame/grey lens	+1.5	universal
3745C	Black frame/grey lens	+2.0	universal
3745D	Black frame/grey lens	+2.5	universal
3746B	Black frame/indoor-outdoor lens	+1.5	universal
3746C	Black frame/indoor-outdoor lens	+2.0	universal
3746D	Black frame/indoor-outdoor lens	+2.5	universal



## ProVision® Cool Wraps™ Bifocal

- Sleek, lightweight wraparound design
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Diopter	Size/Facial Fit
3730A	Clear frame/clear lens	+1.0	universal
3730B	Clear frame/clear lens	+1.5	universal
3730C	Clear frame/clear lens	+2.0	universal
3730D	Clear frame/clear lens	+2.5	universal
3730E	Clear frame/clear lens	+3.0	universal



## ProVision® Tech Specs™ Bifocal

- Lightweight design that provides excellent protection and comfort
- Contoured frame style hugs the face with rubberized (non-latex) attachments on the nosepiece and ear rests
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Diopter	Size/Facial Fit
3720A	Grey frame/clear lens	+1.0	universal
3720B	Grey frame/clear lens	+1.5	universal
3720C	Grey frame/clear lens	+2.0	universal
3720D	Grey frame/clear lens	+2.5	universal
3720E	Grey frame/clear lens	+3.0	universal

## Magnification Collection



### ProVision® Econo Loupes

- Full 50mm diameter magnifying lens for a fully-enhanced clinical field
- Economic, lightweight alternative to traditional loupes
- Reduces eye strain and increases visual precision
- Wraparound frame and soft flexible nose bridge to minimize slippage and provide all day comfort
- Soft tip flexible temples ensure a secure fit
- Available in +1.5, +2.0 and +2.5 diopters
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref #	Description	Diopter	Size/Facial Fit
3780B	Black frame/clear lens	+1.5	universal
3780C	Black frame/clear lens	+2.0	universal
3780D	Black frame/clear lens	+2.5	universal

## A Guide To ProVision® Econo Loupes Full 50mm Magnifying Lens

**NEW**

Designed to meet the ANSI Z87.1 standard



### Magnification

**Full 50mm  
Magnifying Lens**



**V.S.**

**Bifocal Lens**



### Did You Know?

Palmero's ProVision Econo Loupes was designed with loupes users in mind, featuring a lightweight frame and full magnifying lens that can be used as an alternative to the traditional loupes.

- The use of magnification can help reduce chronic back pain while it improves both the vision and posture of the clinician.
- Investing in traditional loupes can be significant while our Econo Loupes provide an economical magnifier, for a fraction of the cost, providing the needed clarity, magnification, and depth of field to improve clinical outcomes.




## Disposables Collection



### Dynamic Disposables™ Safety Eyewear

- Convenient, disposable, cost-effective alternative eye protection
- Disposable eyewear can be worn over prescription eyewear
- Reusable frames reduces costs
- OSHA blood-borne pathogen standard compliance
- Fog-free, anti-glare, anti-static, and optically clear
- Does not contain latex



Ref #	Description	Quantity
<b>Office/Value Packs</b>		
3900	Office Pack	10 black frames, 20 clear lens & 1 cord
3905C	Office Pack	10 black frames, 20 grey lens & 1 cord
3906C	Office Pack	10 black frames, 20 orange bonding lens & 1 cord
3908	Value Pack	25 black frames, 100 clear lens & 5 cords
<b>Frames Replacement</b>		
3901	Black frames	10 pack
<b>Lens Replacement</b>		
3902	Clear lens	100 pack
3903	Clear lens	25 pack
3905	Grey lens	25 pack
3905B	Grey lens	100 pack
3906	Orange bonding lens	10 pack
3906B	Orange bonding lens	100 pack

### Dynamic Disposables™ Safety Eyewear Pink



Ref #	Description	Quantity
<b>Office/Value Packs</b>		
3915	Office Pack	10 black frames, 20 clear lens with pink stripe & 1 cord
<b>Frames Replacement</b>		
3901	Black frames	10 pack
<b>Lens Replacement</b>		
3916	Clear lens/pink stripe	25 pack
3917	Clear lens/pink stripe	100 pack

### Did You Know?

When you purchase our new Dynamic Disposables PINK Safety Eyewear, 5% of the proceeds will be donated to the National Breast Cancer Foundation, Inc.

**The National Breast Cancer Foundation's (NBCF)** mission is to provide assistance and inspire hope to those affected by breast cancer through early detection, education, and support services. NBCF provides free mammograms and breast cancer education to women in all 50 states through a network of hospitals. NBCF was founded in 1991 and has a Charity Navigator's highest 4-Star rating. To learn more, visit [nationalbreastcancer.org](http://nationalbreastcancer.org).



**CLINICAL  
PROBLEM  
SOLVER**



## Dynamic Disposables™ Snapeez™ Face Shields

- As an addition to our traditional face shields line, the Snapeez is the next generation of face shields with an easy-to-use snap-on disposable lens feature
- Ergonomic frame wraps around the head for a more secure fit, designed to be worn with or without prescription eyewear, bridge-free nose means a more comfortable fit
- Self-sealing face shields prevent exposure to bloodborne pathogens
- Economical Reusable Frame reduces the cost per use
- Easy-to-Use and Convenient
- Disposable and cost effective face shields fit over eye glasses and face masks
- Ultra-Clear™ Technology provides anti-fog, anti-static, and optically clear
- Meets OSHA bloodborne pathogen guidelines
- Does not contain latex



Ref #	Description	Size	Quantity
3918	Office Pack	7.25"	10 frames, 20 full shields & 1 cord/box
3919	Full Shields	7.25"	25 full face shields/box
3920	Office Pack	5.5"	10 frames, 20 half shields & 1 cord/box
3921	Half Shields	5.5"	25 half face shields/box



## Dynamic Disposables™ Traditional Face Shields

- Disposable, lightweight full-face shields offer complete facial coverage below the chin
- Protects the eyes, nose, and mouth from aerosols, sprays, and splatters
- Soft elastic and foam forehead band provides a comforting snug fit
- Meets OSHA bloodborne pathogen guidelines
- Does not contain latex



Ref #	Description	Size	Quantity
1964	Clear full face shield	8" length	24/box

	Traditional Face Shield	NEW Snapeez™ Face Shield
Products		
Unique Features	<p><b>Foam Design</b> Soft, stretch elastic strap for a secure and comfortable fit</p>	<p><b>Wraparound Frame</b> Ergonomic frame wraps around the head for a more secure fit, designed to be worn with or without prescription eyewear, bridge-free nose means more comfortable fit</p> <p><b>Easy Snap-On Attachment</b> Self-sealing face shields prevent exposure to bloodborne pathogens</p> <p><b>Economical Reusable Frame</b> Reduce the cost per use</p>
Common Features	<p><b>Easy-to-Use and Convenient</b></p> <p><b>Disposable &amp; Cost Effective</b> Cost effective face shields fit over eye glasses and face masks</p> <p><b>Ultra-Clear Technology</b> Anti-fog, anti-static, and optically clear</p> <p><b>Reliability &amp; Safety</b> Meets OSHA bloodborne pathogen guidelines</p> <p><b>Does not contain latex</b></p>	
Length	8" Full Face Shield	7.25" Full Face Shield 5.5" Half Face Shield

### Did You Know?

ProVision Safety Eyewear and Goggles wraparound lens designs and Dynamic Disposable Face Shields meet OSHA and CDC guidelines and protect against splashes, sprays and splatter from bloodborne pathogens.

“When splashes, sprays, splatters, or droplets of blood or OPIM pose a hazard to the eyes, nose or mouth, then masks in conjunction with eye protection (such as goggles or glasses with solid side shields) or chin-length face shields must be worn.”

Source: US Department of Labor - Occupational Safety and Health Administration.  
Personal Protective Equipment (PPE) Reduces Exposure to Bloodborne Pathogens  
[https://www.osha.gov/OshDoc/data\\_BloodborneFacts/bbfact03.pdf](https://www.osha.gov/OshDoc/data_BloodborneFacts/bbfact03.pdf)

## Bonding Collection

Palmero takes eye conservation seriously by offering a comprehensive line of disposable and reusable eyewear/shields for dental curing light procedures. Dental curing light intensity has increased and studies show many offices do not practice proper eye protection protocols from exposure to the dental curing light. Offer eye protection for everyone in the operator who could be exposed to the bright light with the appropriate orange (blue light blocking) safety glasses to prevent ocular damage.

Here are four factors to consider when selecting the proper eyewear:

- All bonding eyewear is tested and proven to protect blue light rays between 385-500 nm which exceeds the blue light wavelengths of 450-495 nm
- Select from orange or red polycarbonate lenses for maximum protection
- Filters up to 95% of blue, UVA, and UVB rays
- Comfortable, lightweight, fog-free, and scratch-resistant



### ProVision® Bonding Wraps™

- Durable lightweight eyewear with sleek styling and a close fit
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3600	Red bonding wraps	full size



### ProVision® Mini Bonding Wraps™

- Mini version of our popular full size sleek, bonding eyewear with the same great features and benefits
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3607B	Red bonding wraps	child/youth



### Dynamic Disposables™ Bonding Frames and Lenses

- Disposable lens with reusable ergonomically framed eyewear that can be worn over prescription glasses
- Meets OSHA bloodborne pathogen standard



Ref #	Description	Quantity
3901	Black frames	10 pack
3906	Orange bonding lens	10 pack
3906B	Orange bonding lens	100 pack
3906C	10 black frames, 20 orange bonding lens & 1 cord*	Office Pack

*\*see page 31 for more information*



### ProVision® Infinity™

- Replaceable bonding lens for ProVision Infinity Safety Eyewear (see page 21)
- Meets ANSI Z87.1 standard

Ref #	Description
3613BR	Red bonding replacement lens



### ProVision® Bond Oasis™

- Attractive over-the-glasses style with built-in side shields
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3570BOND	Black frame/reddish-orange lens	universal



### ProVision® Econo Bonding

- One piece goggle design with detachable side shields
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
1100	Red Econo Bonding	universal



### ProVision® Eyesaver Sleeks™ Bonding

- Sleek, ultra-light, economical goggles with adjustable temples to allow for broader facial fit
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
16SLK	Red bonding frame/lens	universal



### ProVision® Tech-Specs™ Bonding

- Stylish, impact-and-stress resistant material
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3712	Grey frame/red lens	universal



### ProVision® Bond OverBond™

- Designed to wear over prescription eyewear with telescoping and ratcheting arms to ensure a great fit
- Meets ANSI Z87.1 standard

Ref #	Description	Size/Facial Fit
3602BOND	Black frame/orange lens	universal



### ProVision® Bonding Lens Shield


- Easily adapts to most pistol or wand type VLCs (visible light cure) to help protect eyes from the harmful blue light used in curing and bleaching procedures
- Constructed of an unbreakable polycarbonate with a protection of 530nm

Ref #	Description
125	3" diameter

## Did You Know?

Palmero ProVision Eyewear portfolio includes specialized bonding eyewear materials that filter out the harmful Blue Light to the back of the eye, protecting the optic nerve and the retina.

- Blue light is being used for multiple dental products and procedures including Bonding, Curing, and Whitening agents.
- Blue Light Radiation has been proven to cause moderate to severe retinal damage to the Dental Professionals and even patients who experience long term exposure over a period of time without wearing any type of optical protection.
- Using eyewear filters to block UV and the blue wavelength of light from 380-495nm is an important eye conservation measure for dental professionals and patients.
- Blue light affects the back of the eye, the retina.



Color	Wavelength
Violet	380-450 nm
Blue	450-495 nm
Green	495-570 nm
Yellow	570-590 nm
Orange	590-620 nm
Red	620-750 nm

## Laser Collection

Palmero offers 3 premium laser safety eyewear products.

- Meets ANSI Z87.1 and Z136.1 standard for laser protection
- Optical Density OD >5+
- Filter Specifications 800-830 nm
- Includes upscale storage kit with zippered hard case, draw string bag and lanyard



### ProVision® Diode Wraparound Laser Eyewear

- Wraparound black frame with yellow/green polycarbonate filter lens
- Rubber nose pad and foam pad insert along brow provides comfortable fit while the wraparound design ensures your eyes are protected
- Visible Light Transmission (VLT) 68%

Ref # 3589YG



### ProVision® Diode Alexandrite Design Laser Eyewear

- Silver/black frame with green polycarbonate filter lens
- Rubber nose pad and rubberized temple bars provides comfortable grip while the wraparound design ensures your eyes are protected
- Visible Light Transmission (VLT) 50%

Ref # 3588DA



### ProVision® Diode Uni-Fit Wrap Laser Eyewear

- Lightweight, universal-fit, black wrap frame
- Yellow/green polycarbonate filters lens
- Rubber nose pad and wraparound design ensures your eyes are protected
- Visible Light Transmission (VLT) 68%

Ref # 3590YG

Ref #	Wavelength	Optical Density (OD)	Lens Color	Visible Light Transmission (VLT)
3589YG	800-830 nm	>5+	Yellow/Green	68%
3588DA	800-830 nm	>5+	Green	50%
3590YG	800-830 nm	>5+	Yellow/Green	68%



## Did You Know?

When selecting laser eyewear there are 4 important factors to consider when choosing the proper laser eyewear.

- Consult your laser manual to determine the protection level needed based on laser's output and the recommended Optical Density (OD)
- Select a filter whose specifications match the manual and provides the highest Visible Light Transmission (VLT)
- Find a frame that provides a comfortable fit
- Provide laser eye protection for clinicians, staff and patients

## Eyewear Maintenance



### ProVision® Crystal Clear Anti-Fog Solo Wipes

- Cleans lenses while providing long-lasting, anti-fog protection
- Safe and effective for use on uncoated lenses including polycarbonate, acrylic, and glass
- Dispenser box contains 100 premoistened individually wrapped towelettes
- Great for eyewear and goggles

Ref #	Size	Quantity
3534S	8" x 5"	100 wipes/box



### ProVision® Crystal Clear Optical Cleaning Wipes

- Convenient wipes premoistened with anti-fog, anti-static lens cleaning and polishing solution
- Cleans dirt, grease, and debris from eyewear
- Container/lid does not contain latex

Ref #	Size	Quantity
3534	6" x 6 3/4"	160 wipes/canister



### ProVision® Anti-Fog Solution

- Compatible with all eyewear, face shields, mirrors, loupes, and scopes
- Resists scratches, static, and stains
- Foam formula yields 5 times the applications compared to a traditional pump or squeeze bottle
- Alcohol-free, non-toxic, and non-flammable

Ref #	Size	Quantity
3533	8 oz. (237 ml)	1 bottle
3536	1.7 oz. (51 ml)	1 bottle



### ProVision® Lens Cleaning Station

- Convenient eyewear lens cleaning station with anti-fog and anti-static additives
- Designed for countertop placement or wall mounting with hanging hardware
- Station includes 8 oz. pump spray bottle of optical cleaning solution and 600 dry wipes

Ref #	Size	Quantity
3535	7 3/4" x 8" x 5"	1 station

## Did You Know?

Palmero has a full line of lens care products and disinfectants to meet CDC guidelines for eyewear cleaning and disinfection.

- Reusable protective eyewear should be cleaned with soap and water when visibly soiled as well as disinfected between patients.

Source: US Centers for Disease Control and Prevention. Personal Protective Equipment.  
<https://www.cdc.gov/oralhealth/infectioncontrol/faqs/personal-protective-equipment.html>

## Eyewear Accessories



### Dynamic Disposables™ Side Shields

- Impact-resistant, slip-on plastic clear side shields
- Quick and easy to apply on virtually any frame arms
- Turn your glasses into safety eyewear
- Available in rigid or flexible plastic material



Ref #	Description
3611	Rigid shields
3612	Flexible shields

Eyewear not included. Shields will not project against severe impact hazards or chemical splash.



### Dynamic Disposables™ Eyewear Cords

- Flexible, multicolor classic string style cord, 24" length cords
- Eliminate lost time searching for misplaced eyewear
- Allows eyewear to hang conveniently around neck when not in use

Ref #	Description	Quantity
3904	Multicolor classic string style cords	5 pack

Eyewear not included



### Bryte-Syte™ Clip-On Headlight 4000

- Universal clip on mount, unthethered
- ¼-watt LED delivers bright white light, 17 lumens
- Comfortable, lightweight light & clip with on/off switch control, weighing only 1 oz.
- Secure clip design for all major TTL or Flip Up Loupes
- Runs on two 3-volt CR2032 lithium coin battery (included)
- Provides 100-hour battery life

Ref #	Description
4000	Low-Power Clip-On Headlight 4000





### Bryte-Syte™ Clip-On Headlight 4006

- Lightweight, universal clip on mount weighing only 1 oz.
- Corded with belt-clip battery pack with variable intensity control
- 1-watt LED delivers bright white light, 500 lumens
- Secure clip design for all major TTL & Flip Up Loupes
- Runs on 3 AAA Batteries (one set included)

---

Ref #	Description
4006	High-Power Clip-On Headlight 4006

---



### Bryte-Syte™ High Power Headband Headlight

- Extremely mobile option to deliver light into the field
- No lengthy wires, battery pack integrated in headband, weighs 8.2 oz.
- Adjustable crown and lateral headband for comfortable fit with variable intensity control
- Provides 20-hour battery life, 100,000 hour LED life
- 1-watt LED delivers bright white light, 500 lumens
- Runs on 3 AAA Batteries (one set included)

---

Ref #	Description
4008	High-Power Headband Light

---

#### Did You Know?

Palmero's ProVision eyewear accessories include light weight headlights that deliver brilliant lines of sight light into the field.

## X-RAY APRONS & HANGERS

### Cling Shield®

- FDA registered, medical-grade textured vinyl backing that clings in place ensuring coverage by preventing slippage
- Made in the USA and custom-made to order
- Available in 2 material types of Vinyl or Nylon and Lead-Free (.3MM) or Lead Lined (.3MM or .5MM)
- Available in a variety of colors and selected silk screen prints

#### Did You Know:

- Do not fold aprons as it can cause cracking within the lining. See Palermo's assortment of hangers to protect your investment safely.
- All lead aprons have a useful life as their shielding properties wear down gradually. Proper maintenance includes either checking via fluoroscopy or replace if they show any signs of damage.
- The EPA states that recycling the lead apron is the preferred method since it keeps lead out of landfills. Contact Palermo for our environmentally friendly X-Ray Apron Recycle/Trade-in program.

## Cling Shield® Child Protectall™ Apron with Integrated Thyroid Collar

NEW

**PROTECTALL™** design has a built in thyroid collar with velcro clasp

Available in vinyl or nylon in 10+ colors including new Smiley Face apron for children

**Cling Shield®** medical-grade textured vinyl backing prevents sliding and provides stability



Integrated nylon loops allow for easy hanging and proper storage

When not in use, store your apron properly on a rack or apron hanger specifically designed for aprons

High density lead and lead free materials effectively shield radiation

Easy to clean materials that can be disinfected using **DisCide® Ultra** surface disinfectant and cleaned with **Protectall™ Upholstery Cleaner Wipes** to safely protect your investment \*\*



MADE IN  
USA

\*\*With proper use. Data on file.

### Did You Know?



RECYCLE/TRADE-IN  
PROGRAM

Buy 1, Recycle 1,  
Receive 1 apron hanger

### X-Ray Apron Recycling Program\*

1. **Purchase** a new Palermo x-ray apron.
2. **Collect** old or unwanted x-ray apron of any brand.
3. **Submit** invoice as proof of purchase to obtain and receive an RA# via Fax 203.377.8988 or email customerservice@palmerohealth.com. Include ship to address, phone, and email.
4. **Contact** Customer Service for a Return of Authorization Number (RA#) to register for recycling.
5. **Ship** your old apron with RA# for responsible disposal to: Palermo Healthcare LLC, Attn: X-Ray Recycle/Trade-In Program, 120 Goodwin Place, Stratford, CT 06615.
6. **Acquire** a new x-ray apron hanger. Once your recycled apron is received, we'll send you an \$21.85 retail value hanger for your participation!

\*Continental U.S. only



Smiley

### Child Protectall™ Apron with Neck Collar

- Neck collar permanently attached. Covers from thyroid to below the gonadal area. Sized for petite adults or children. Velcro clasp on neck collar. Loops attached for easy hanging and storage.
- 20" x 20"

Ref #	Description	Thickness	Color	Weight
27USMILEY	Lead-Free	.3MM	Slate Blue ONLY	3 lbs
27SMILEY	Lead-Lined	.3MM	Slate Blue ONLY	3 lbs
27UHAPPYFACE	Lead-Free	.3MM	Blue and White Stripes ONLY	4 lbs
27HAPPYFACE	Lead-Lined	.3MM	Blue and White Stripes ONLY	4 lbs
27U	Lead-Free	.3MM		7 lbs
27	Lead-Lined	.3MM		7 lbs
27M	Lead-Lined	.5MM		7 lbs



Happy Face

### Adult Protectall™ Apron, with Neck Collar

- Neck collar permanently attached. Covers from thyroid to below the gonadal area. Velcro clasp on neck collar. Loops attached for easy hanging and storage.
- 22 ¼" x 25 ½"

Ref #	Description	Thickness	Weight
24U	Lead-Free	.3MM	5 lbs
24	Lead-Lined	.3MM	6 lbs
24M	Lead-Lined	.5MM	12 lbs



Smiley

### Child Apron, No Collar

- Covers from over shoulders to below the gonadal area. Loops attached for easy hanging and storage.
- 20" x 19 ½"

Ref #	Description	Thickness	Color	Weight
22USMILEY	Lead-Free	.3MM	Slate Blue ONLY	3 lbs
22SMILEY	Lead-Lined	.3MM	Slate Blue ONLY	4 lbs
22UHAPPYFACE	Lead-Free	.3MM	Blue and White Stripes ONLY	3 lbs
22HAPPYFACE	Lead-Lined	.3MM	Blue and White Stripes ONLY	4 lbs
22U	Lead-Free	.3MM		3 lbs
22	Lead-Lined	.3MM		4 lbs
22M	Lead-Lined	.5MM		7 lbs





### Adult Apron, No Collar

- Covers from over the shoulders to below the gonadal area. Loops attached for easy hanging and storage.
- 22 ¼" x 24 ¼"


Ref #	Description	Thickness	Weight
20U	Lead-Free	.3MM	4 lbs
20	Lead-Lined	.3MM	5 lbs
20M	Lead-Lined	.5MM	8 lbs



### Child Pano Dual Apron, No Collar

- Covers in front from shoulders to below gonadal area. Back is the same length as front. Velcro clasp at shoulder.
- 19 7/8" x 19 1/2"

Ref #	Description	Thickness	Color	Weight
26CUSMILEY	Lead-Free	.3MM	Slate Blue ONLY	7.50 lbs
26CSMILEY	Lead-Lined	.3MM	Slate Blue ONLY	9 lbs
26CUHAPPYFACE	Lead-Free	.3MM	Blue and White Stripes ONLY	7.50 lbs
26CHAPPYFACE	Lead-Lined	.3MM	Blue and White Stripes ONLY	9 lbs
26CU	Lead-Free	.3MM		7.50 lbs
26C	Lead-Lined	.3MM		9 lbs
26CM	Lead-Lined	.5MM		11 lbs



### Adult Pano Dual Apron, No Collar

- Covers in front from shoulders to below gonadal area. Back is the same length as front. Velcro clasp at shoulder.
- 22 ½" x 26 ½"

Ref #	Description	Thickness	Weight
26U	Lead-Free	.3MM	8 lbs
26	Lead-Lined	.3MM	10 lbs
26M	Lead-Lined	.5MM	15 lbs



### Adult Lap Apron

- Covers pelvic and gonadal area. Extra long straps can be tied in the front or back.
- 24" x 18 1/8"

Ref #	Description	Thickness	Weight
35LU	Lead-Free	.3MM	2.50 lbs
35L	Lead-Lined	.3MM	4 lbs
35LM	Lead-Lined	.5MM	4.50 lbs



### Adult Deluxe 3/4 Pano Dual Apron, No Collar

- Covers from shoulders to below gonadal area. Back covers from shoulders to mid-lumbar area. Velcro closing for easy size adjustment.
- 22 1/2" x 30"

Ref #	Description	Thickness	Weight
23U	Lead-Free	.3MM	7.50 lbs
23	Lead-Lined	.3MM	9 lbs
23M	Lead-Lined	.5MM	11 lbs



### Adult Pano Cape Apron, No Collar

- Lays over the shoulders, partially down the front and back with Velcro closure on either front or back. Loops attached for easy hanging and storage.
- 23 1/2" x 7 1/2"

Ref #	Description	Thickness	Weight
29U	Lead-Free	.3MM	2.50 lbs
29	Lead-Lined	.3MM	4 lbs
29M	Lead-Lined	.5MM	4.50 lbs



### Adult Pano Stole Apron, No Collar

- Lays over the shoulders to mid-thorax area in front. Back attaches with crisscross Velcro fasteners for easy size adjustment.
- 30" x 12 1/2"

Ref #	Description	Thickness	Weight
29SU	Lead-Free	.3MM	2.50 lbs
29S	Lead-Lined	.3MM	4 lbs
29SM	Lead-Lined	.5MM	4.50 lbs



## Adult Tech Apron

- Covers and protects technician from shoulders to below the knees. Ties behind neck and waist straps wrap around the back and tie in front of technician waist.
- 23 1/8" x 35"

Ref #	Description	Thickness	Weight
21U	Lead-Free	.3MM	6 lbs
21	Lead-Lined	.3MM	7 lbs
21M	Lead-Lined	.5MM	12 lbs



## Adult Coat Apron

- Full-length frontal protection. Wraps around back over shoulder blades. Ties at upper and lower back. 23" x 37 1/2"
- Back is 10" x 17 3/4"

Ref #	Description	Thickness	Weight
34U	Lead-Free	.3MM	7 lbs
34	Lead-Lined	.3MM	8 lbs
34M	Lead-Lined	.5MM	12 lbs



## Child Neck/Thyroid Collar

- For the protection of the thyroid area. Velcro closing. One size fits all. Use in conjunction with Child Patient Apron (page 41).
- 17" x 4"
- Not available in Silk Screen Facings

Ref #	Description	Thickness	Weight
25CU	Lead-Free	.3MM	9 oz
25C	Lead-Lined	.3MM	1 lb
25CM	Lead-Lined	.5MM	1.5 lbs



## Adult Neck/Thyroid Collar

- For the protection of the thyroid area. Velcro closing. One size fits all. Use in conjunction with Adult Patient Apron (page 42).
- 22 5/8" x 4 1/4"
- Not available in Silk Screen Facings

Ref #	Description	Thickness	Weight
25U	Lead-Free	.3MM	9 oz
25	Lead-Lined	.3MM	1 lb
25M	Lead-Lined	.5MM	1.5 lbs

### Did You Know?

Aprons should be stored flat or on an apron hanger. To avoid creasing, it should never be folded. Aprons should also be visually examined at frequent intervals and evaluated for potential cracks in the protective material no less than annually.



**Buy 1, Recycle 1,  
Receive 1 apron hanger**

Contact Palmero for our environmentally friendly X-Ray Apron Recycle/Trade-in program\*. Visit [palmerohealth.com](http://palmerohealth.com), call 800.344.6426 or email [customerservice@palmerohealth.com](mailto:customerservice@palmerohealth.com).

*\*Continental U.S. only*

# Hold-It™ X-Ray Apron Hangers and Accessories



## Standard Apron Hanger

- White-coated steel hanger
- Compatible with adult and child Patient and Protectall™ style aprons
- Hanging hardware included

---

Ref #	Description
28	15" x 2"

---



## Economy Apron Hanger

- White-coated steel hanger
- Not compatible with Dual, Coal, or Tech aprons
- Hanging hardware included

---

Ref #	Description
31	15" x 2"

---



## CC Style Apron Hanger

- Sturdy chrome-coated hanger
- Compatible with adult and child Patient and Protectall™ style aprons
- Hanging hardware included

---

Ref #	Description
32	15"

---



## Economy X-Ray Apron Hanger

- White-coated steel hanger
- Compatible with adult and child Patient and Protectall™ style aprons, Coat and Dual style aprons

---

Ref #	Description
33	20 7/8" x 10 1/4"

---



## Deluxe Coat Panoramic Apron Hanger

- Sturdy chrome-coated hanger
- Compatible with Panoramic aprons, Coat and Dual style aprons

---

Ref #	Description
30	21" x 9 1/4"

---



## Protectall™ Upholstery Cleaner Wipes

- Premoistened wipes containing an advanced water-based silicone formula
- Safely protects, cleans, and restores leather, vinyl, plastic, and rubber upholstery surfaces like x-ray aprons and dental unit chairs
- Recommended for periodic use after disinfecting upholstery surface with DisCide® Ultra surface disinfectant
- Evenly wipe down the surfaces with Protectall™ wipes and allow to dry to touch

---

Ref #	Size	Quantity
3800	10" x 12"	50 count canister

---

## CLING SHIELD® APRONS

		LINING	
		LEAD-FREE	
		.3MM	
APRONS	SIZE	WT.(LBS)	REF #
CHILD PROTECTALL™ APRON, WITH NECK COLLAR	20" X 20"	3.00	27U +(COLOR)
ADULT PROTECTALL™ APRON, WITH NECK COLLAR	22¼" X 25¼"	5.00	24U +(COLOR)
CHILD APRON, NO COLLAR	20" X 19½"	3.00	22U +(COLOR)
ADULT APRON, NO COLLAR	22¼" X 24¼"	4.00	20U +(COLOR)
CHILD PANO DUAL APRON, NO COLLAR	19⅞" X 19½"	7.50	26CU +(COLOR)
ADULT PANO DUAL APRON, NO COLLAR	22½" X 26½"	8.00	26U +(COLOR)
ADULT LAP APRON	24" X 18⅛"	2.50	35LU +(COLOR)
ADULT DELUXE PANO DUAL APRON, NO COLLAR	22½" X 30"	7.50	23U +(COLOR)
ADULT PANO CAPE APRON, NO COLLAR	23½" X 7½"	2.50	29U +(COLOR)
ADULT PANO STOLE APRON, NO COLLAR	30" X 12½"	2.50	29SU +(COLOR)
ADULT TECH APRON	23⅓" X 35"	6.00	21U +(COLOR)
ADULT COAT APRON	23" X 37½" (FRONT), 10" X 17¾" (BACK)	7.00	34U +(COLOR)

## CLING SHIELD® COLLARS

CHILD NECK/THYROID COLLAR	17" X 4"	9 OZ.	25CU +(COLOR)
ADULT NECK/THYROID COLLAR	22⅝" X 4¼"	9 OZ.	25U +(COLOR)

### Vinyl Colors

With Black Binding



Black

Steel  
Grey

White

Taupe

Mauve

Slate  
Blue

Royal  
Blue

Teal

Forest  
Green

### Tough Nylon Colors

With Black Binding



Navy

Aqua

Gold

Burgundy

Red

With Grey Binding



Light Grey



LINING				MATERIAL & COLORS				
LEAD-LINED				NYLON*	VINYL*			
.3MM		.5MM		COLORS	COLORS	SILK SCREENING DESIGN **		
WT.(LBS)	REF #	WT.(LBS)	REF #	ALL	ALL	BEAR	TOOTH	HAPPY FACE
4.00	27 +(COLOR)	7.00	27M +(COLOR)	✓	✓	✓	✓	✓
6.00	24 +(COLOR)	12.00	24M +(COLOR)	✓	✓	✓	✓	✓
4.00	22 +(COLOR)	7.00	22M +(COLOR)	✓	✓	✓	✓	✓
5.00	20 +(COLOR)	8.00	20M +(COLOR)	✓	✓	✓	✓	✓
9.00	26C +(COLOR)	11.00	26CM +(COLOR)	✓	✓	✓	✓	✓
10.00	26 +(COLOR)	15.00	26M +(COLOR)	✓	✓	✓	✓	✓
4.00	35L +(COLOR)	4.50	35LM +(COLOR)	✓	✓	N/A	N/A	N/A
9.00	23 +(COLOR)	11.00	23M +(COLOR)	✓	✓	N/A	N/A	N/A
4.00	29 +(COLOR)	4.50	29M +(COLOR)	✓	✓	N/A	N/A	N/A
4.00	29S +(COLOR)	4.50	29SM +(COLOR)	✓	✓	N/A	N/A	N/A
7.00	21 +(COLOR)	12.00	21M +(COLOR)	✓	✓	N/A	N/A	N/A
8.00	34 +(COLOR)	12.00	34M +(COLOR)	✓	✓	N/A	N/A	N/A

1 LB.	25C +(COLOR)	1.5 LBS.	25CM +(COLOR)	✓	✓	N/A	N/A	N/A
1 LB.	25 +(COLOR)	1.5 LBS.	25M +(COLOR)	✓	✓	N/A	N/A	N/A

**Vinyl Colors**

With Grey Binding



Grey

With Beige Binding



Beige

With Blue Binding



Blue



**BEAR** prints black on white vinyl with black binding



**TOOTH** prints blue on white vinyl with blue binding

With Blue Binding



Happy Face Blue and White Stripes

With Black Binding



Smiley Slate Blue



**HAPPY FACE** prints yellow on blue and white vinyl with blue binding



**SMILEY** prints gradient yellow on slate blue vinyl with black binding (Ref# 27, 22, and 26C only)

\* = For Vinyl or Nylon material aprons, specify the color description (i.e. 27U Black) with the part number.

\*\* = For silk screening designs, add "BEAR," "TOOTH," or "HAPPYFACE" (i.e. 27U after the Lead-Free and Lead-Lined part number).

N/A = Not available in Silk Screening Designs

# ORGANIZE



## Let's Get It Together

A well-organized operatory is the backbone of every well-run practice. At Palmero Healthcare, we offer the systems and dispensers to put the products you need in the places you want them. From wall-mounted holders, and assorted dispensers to bur blocks and industry-specific containers, your search for organization ends here. Let's get together and get organized.

## HOLDERS

### Hold-It™ Holder

- Durable holder can be mounted easily in accessible areas to prevent clutter
- Save valuable counter space, directly mounts to a cabinet or wall
- Fits all Palmero towelette containers and many other manufacturer's brands
- Mounting hardware included



### Hold-It™ Canister Holder

- Compatible with DisCide® ULTRA towelettes canister (see Page 10)

Ref #	Description
109	6" x 5 5/8" holder



### Hold-It™ Locking Canister Holder

Ref #	Description
109L	7 1/4" x 5 5/8" holder
109LK	Padlock and locking key





## Hold-It™ Round Quart Bottle Holder

Ref #	Description
3514	5 ¾" x 3 ¾"



## Hold-It™ Square Quart Bottle Holder

- Compatible with DisCide® ULTRA quart bottle (see Page 10)

Ref #	Description
3517	4 ¾" x 4" x 5 ½"




## Hold-It™ Face Mask Holder

Ref #	Description
1856	5" x 6" x 5"



## Hold-It™ Glove Box Holder

Ref #	Description
1858	Single/double box holder, 11" x 4" x 4"
1858T	Triple box holder, 12 ¼" x 4 ½" x 15"



## Hold-It™ Tissue Box Holder

Ref #	Description
1859R	Rectangular tissue box holder, 9 ¾" x 2 ½" x 4 ¼"
1859S	Square tissue box holder, 5" x 5" x 4 ¼"



## Hold-It™ Box Holder

Ref #	Description
1854	9" x 4 ½" x 6 ½"

## BUR BLOCKS

### Clear Plexiglass Bur Blocks

- An attractive, clear plexiglass bur block with 16-, 27-, 36-, 50-, 72-, or 100-hole configurations to hold and organize burs
- Available in friction grip, latch, or combo (half friction grip/half latch)
- Non autoclavable



### 16-Hole Bur Blocks

- 2 13/16" x 1 ¼" x ½", 7/16" holes spacing

Ref #	Description
1516FG	Friction grip
1516LATCH	Latch style



### 27-Hole Bur Blocks

- 2 13/16" x 1 7/8" x 5/16", 7/16" holes spacing

Ref #	Description
1516AFG	Friction grip
1516ALATCH	Latch style



### 36-Hole Bur Blocks

- 5" x 2" x 1/2", 10/16" holes spacing

Ref #	Description
1516BFG	Friction grip
1516BHALF	Combo
1516BLATCH	Latch style



### 50-Hole Bur Blocks

- 5 5/16" x 2" x 1/2", 5/16" holes spacing

Ref #	Description
1516CFG	Friction grip
1516CHALF	Combo
1516CLATCH	Latch style



### 72-Hole Bur Blocks

- 5" x 4" x 1/2", 5/16" holes spacing

Ref #	Description
1516DFG	Friction grip
1516DHALF	Combo
1516DLATCH	Latch style



### 100-Hole Bur Blocks

- 5 1/4" x 3 3/4" x 1/2", 5/16" holes spacing

Ref #	Description
1516EFG	Friction grip
1516EHALF	Combo
1516ELATCH	Latch style

## Anodized Aluminum Bur Blocks

- Autoclavable, dry-heat sterilized, chemiclaved, or placed in an ultrasonic cleaner
- Preattached cover secures the burs in place, 1/4" spacing between holes
- Available in friction grip, latch, or combo on selected 10-, 15-, 21-, or 30-hole bur blocks



### 10-Hole Bur Blocks

- 1 5/8" x 1 1/8"

---

Ref #	Description
T10CHALF	Combo Bur Block 10-hole, 5 friction grip & 5 latch
T10CLATCH	Latch Bur Block 10-hole

---



### 15-Hole Bur Blocks

- 1 5/8" x 1 1/8"

---

Ref #	Description
T15CFG	Friction Grip Bur Block 15-hole
T15CLATCH	Latch Bur Block 15-hole

---



### 21-Hole Bur Blocks

- 2 13/16" x 1 1/4"

---

Ref #	Description
T20CFG	Friction Grip Bur Block 21-hole
T20CLATCH	Latch Bur Block 21-hole
T20CHALF	Combo Bur Block 21-hole, 10 friction grip & 10 latch

---



### 30-Hole Bur Blocks

- 2 7/8" x 1 1/4"

---

Ref #	Description
T30CFG	Friction Grip Bur Block 30-hole
T30CLATCH	Latch Bur Block 30-hole

---

## BUR BLOCKS BOX & ORGANIZER



### Bur Blocks/Clear Hinged Box

- Clear boxes, to store or organize multiple items in your dental office
- Fits Palermo's 16-Hole and 27-Hole Bur Blocks
- Non-autoclavable



Ref #	Size
76	2 7/8" x 1 3/16" x 1"
77	2 7/8" x 2" x 2"



### Rotary Bur Blocks Organizer

- Clear plastic with a dome-shaped cover
- Turns 360° for easy access to burs
- 11/16" spacing between holes, 3.5"W x 2.25"H
- Non-autoclavable
- Available in friction grip, latch, or combo mix

Ref #	Description
40	60-hole, friction grip
40A	62-hole, combo (18 friction grip/44 latch)



### Mesh Bur Holder

- Stainless steel mesh bur holder closes tightly, allowing you to clean and sterilize burs, earrings, crowns, piercing studs, or any small objects
- Safely secures small objects during ultrasonic cleaning or sterilization, safe for steam or dry heat

Ref #	Description
1861	6" overall length, 1 1/2" mesh cup diameter (closed position)

## DISPENSERS/STORAGE



NEW

### Film Dispensers

- Durable plastic film dispenser to organize your barrier films
- Available in clear acrylic or white polystyrene

See PROTECT section on Page 13 for Barrier Films



Ref #	Description	Size
1855A	Clear Acrylic Dispenser	7 1/2"L x 5 1/2"W x 7 1/4"H fits 4" x 6" and 2 1/2" x 6" film
1855	White Polystyrene Dispenser	7 1/4"L x 5 1/4"W x 7"H fits 4" x 6" film

Film not included



### Cotton Roll Dispenser

- Hinged lid for easy loading
- Pull-out drawer dispenses up to three rolls while keeping bulk rolls free from contamination
- Weighted non-skid base, glutaraldehyde safe
- Non-autoclavable, available in 3 colors



Available Colors

Ref #	Description	Color
1600WHITE	4" x 2" x 4"	White
1600BLUE	4" x 2" x 4"	Blue
1600BEIGE	4" x 2" x 4"	Beige



### Denture/Retainer Boxes

- A high-impact, one-piece plastic container to hold mouthguards, retainers, bleaching laminates, bridgework, or partials
- Identification label included
- Designed with contoured corners for easy cleaning
- 12-pack assorted color with ID labels



Ref #	Description	Colors	Size
1962	Denture box	blue, yellow, beige, & mauve	3 3/4" x 2 1/2" x 2"
1963	Retainer box	neon pink, neon green, light blue, & yellow	3" x 2 3/4" x 1"





Available Colors

### Cotton Roll Holder

- Base contains over 50 stainless steel pins to securely hold an entire package of cotton rolls
- Hinged lid closes to keep remaining rolls clean when not in use
- Phenol and glutaraldehyde safe
- Steam autoclavable, available in 4 colors
- Specify white, blue, beige, or mauve

Ref #	Description	Color
1426WHITE	2 ¼" x 3 ⅛" diameter	White
1426BLUE	2 ¼" x 3 ⅛" diameter	Blue
1426BEIGE	2 ¼" x 3 ⅛" diameter	Beige
1426MAUVE	2 ¼" x 3 ⅛" diameter	Mauve



Available Colors

### Sponge Dispenser

- Spring action releases pads for one-at-a-time removal, great for gauze or cotton pads, can be used for dry or moistened pads
- Leakproof – holds all sterilizing and germicidal solutions including phenols, weighted base prevents tipping
- Steam autoclavable, available in 4 colors
- Specify white, blue, beige, or mauve

Ref #	Description	Color
1403WHITE	2" x 2" x 2 ½"	White
1403BLUE	2" x 2" x 2 ½"	Blue
1403BEIGE	2" x 2" x 2 ½"	Beige
1403MAUVE	2" x 2" x 2 ½"	Mauve

### Pump Dispenser

- Dispense the right amount of liquid by pressing lightly on the cap
- Swing-away cap springs back to prevent evaporation and spillage
- Unbreakable plastic, 8 oz. capacity



46

Ref #	Description
46	Clear
47	Brown
48	Clear with Alcohol print



47



48

# SUPPORT



## Operatory Essentials

Sometimes it's the little things that can make the biggest impact. Maybe it's the mirrors that reflect on how patients see your practice. Maybe it's the OSHA labels that confirm your commitment to compliance. Or maybe it's the mixing supplies, waste management provisions, patient throws or any of the countless operatory essentials, that you count on every day. All you need to know is that if you need it in your operatory, we've got it in inventory. We are Palmero Healthcare and we'll make sure the little things never become a big deal.



### Candeez™ Flexible Mixing Bowls

- Fun, scented or unscented, flexible mixing bowls
- Unique design and material allows for inside-out folding, does not contain latex
- Made from an extremely pliable PVC allowing for easier and more efficient mixing of impression materials
- Large stabilizing base, ultra-smooth surface for bubble-free mixing
- Use with most automatic alginate mixers
- Available in various sizes and translucent colors



**Please Note:** Bowls are NOT autoclavable. Spatulas are autoclavable to 270°.



#### Scented Kit

Ref #	Color
1532BPS	Bubblegum/Pink
1532TOS	Tangerine/Orange
1532LGS	Lime/Green
1532MAS	Mint/Aquamarine

Kit Contains:  
1 each of Ref 1520 (spatula), 1530 (medium bowl) and 1531 (large bowl)

#### Scented Bowls

Medium 4" x 2 1/2"		Large 4 1/2" x 3"	
Ref #	Color	Ref #	Color
1530BP	Bubblegum/Pink	1531BP	Bubblegum/Pink
1530TO	Tangerine/Orange	1531TO	Tangerine/Orange
1530LG	Lime/Green	1531LG	Lime/Green
1530MA	Mint/Aquamarine	1531MA	Mint/Aquamarine



#### Unscented Bowls

Medium 4" x 2 1/2"		Large 4 1/2" x 3"	
Ref #	Color	Ref #	Color
1530R	Red	1531R	Red
1530P	Purple	1531P	Purple
1530B	Blue	1531B	Blue
1530G	Green	1531G	Green
1530Y	Yellow	1531Y	Yellow
1530C	Clear	1531C	Clear



Clean Candeez bowl with a damp cloth or use Citrisol™ Orange Solvent Towlettes (page 5) for heavier residue



## Candeez™ Ergonomic Spatulas

- Sturdy, high-grade polypropylene handle with flexible tip optimizes mixing ability
- Designed specifically for easy mixing and application
- Autoclavable to 270°, does not contain latex
- Available in scented or unscented, colors, and sizes for every fit



### Scented - 7 ½"L

Ref #	Description	Color
1520BP	Bubblegum scent	Pink
1520LG	Lime scent	Green
1520TO	Tangerine scent	Orange
1520MA	Mint scent	Aquamarine



### Unscented - 7 ½"L

Ref #	Description
1520W	White
1520GY	Grey
1520B	Blue
1520P	Purple
1520M	Mauve
1520Y	Yellow

### Special Sizes

Ref #	Description	Size	Color
1520LW	Left-Handed design	7 ½"L	White
1522W	Jumbo	11"L	White



## Mixing Sticks and Spatulas

- Does not contain latex, disposable plastic mixing stick/spatula
- Double-ended cement spatula for easy mixing and applying (Ref# 1519A)
- Autoclavable to 270°

### Square Ends

Ref #	Description	Size	Quantity
Ref 1517	Mixing Stick	4 1/2"L	100/bag



### Round Ends

Ref #	Description	Size	Quantity
1519	Alginate Spatula	7 1/2"L	1/bag
1519A	Cement Spatula	6 3/8"L	1/bag
1519L	Alginate Spatula, Extra Long Handle	11 1/4"	1/bag



## Double-Sided Measure Scoop

- Does not contain latex, double-sided alginate tool is essential for proper measuring, especially when dry and wet ingredients need to be kept separate prior to mixing
- Dimensions: 6 1/2"L x 3/4"H, 1 3/4" dia. large scoop & 1" dia. small scoop

Ref #	Description	Color
1525	Translucent scoop	Green



## Sure-Sta™ Mixing Pads

- White, non-skid foam base pads
- Laminated, high-quality bond paper, tear-resistant, and impervious to solutions
- Available in various sizes, approximately 70 sheets/pad

Ref #	Description	Quantity
1514A	3" x 6"	4 pads/bag
1514AS	3" x 6"	1 pad/bag
1514B	6" x 6"	2 pads/bag
1514BS	6" x 6"	1 pad/bag
1514C	3" x 3"	8 pads/bag
1514CS	3" x 3"	1 pad/bag
1514D	6" x 8"	2 pads/bag
1514DS	6" x 8"	1 pad/bag



## Mouth Props

- Autoclavable, silicone mouth props reduce gagging and mouth fatigue
- Ideal for adults or children, finger channel for easy insertion
- Deep corrugation for secure stabilization
- Does not contain latex


Ref #	Description	Size
101	Small child (#01)	1 9/16" x 1 3/8" x 9/16"
101A	Child (#02)	1 3/4" x 1 1/8" x 5/8"
101B	Adult (#10)	1 7/8" x 1 1/4" x 3/4"



## CaRing Tension Relief Band

- Flexible and convenient way to soothe sore, achy fingers
- Stimulate circulation, releasing tension built up after a long day's work
- Contain tiny triangular shaped points that apply light pressure to your fingers
- Promotes healthy hand care and maintenance
- Improve overall stress management

Ref #	Description	Quantity
3810	Tension Relief Band	5 Rings/pack



## OSHA Compliance Labels

- OSHA Compliance Label systems for all healthcare facility needs
- Easily identify all hazardous substances, waste, and radiation areas

Ref #	Description	Size	Quantity
1951	Disinfected sticker	1 5/8" x 1"	250/roll
1957	Caution Radiation	3" x 3"	5/pack
1958	Hazardous Waste	2 1/4" x 3 1/4"	100/roll
1959	Hazardous material	2 1/4" x 4"	100/roll
1960	Biohazard warning	3" x 3"	100/roll
3549	OSHA value pack kit, contains labels 1957, 1958, 1959 & 1960		



## Dis-Pose Infectious Waste Bags

- Designed for safe collection, handling, and disposal of infectious waste
- Meets OSHA and CDC recommended guidelines
- Available in 3 sizes, 100 bags/box, ties included

Ref #	Description	Size
1961A	10-gallon waste bag	24" x 24"
1961B	16-gallon waste bag	24" x 31"
1961C	33-gallon waste bag	41" x 41"



Available Colors

## Chain Napkin Bib Holders

- Autoclavable chrome-plated alligator clips with durable spring and sharp tooth jaws
- Beautifully finished with ball chain
- Available in 6 colors (blue, white, yellow, green, red, and black)

Ref #	Description	Size	Color
310BLUE	Chain (autoclavable)	15" long	Blue
310WHITE	Chain (autoclavable)	15" long	White
310YELLOW	Chain (autoclavable)	15" long	Yellow
310GREEN	Chain (autoclavable)	15" long	Green
310RED	Chain (autoclavable)	15" long	Red
310BLACK	Chain (autoclavable)	15" long	Black



Available Colors

## Chill-Free Napkin Bib Holders

- Chrome-plated alligator clips with durable spring and sharp tooth jaws
- Finished with clear plastic "chill-free" sleeves
- Available in 6 colors (blue, white, yellow, green, red, and black)

*Do not use dry heat sterilization or steam autoclave*

Ref #	Description	Size	Color
311BLUE	Chill-free sleeve	15" long	Blue
311WHITE	Chill-free sleeve	15" long	White
311YELLOW	Chill-free sleeve	15" long	Yellow
311GREEN	Chill-free sleeve	15" long	Green
311RED	Chill-free sleeve	15" long	Red
311BLACK	Chill-free sleeve	15" long	Black



## Lab Apron

- Soft, supple, durable lab aprons are waterproof
- Stain-resistant and washable
- Clear color with tie closures

Ref #	Description	Size
610	Ankle length, clear color	45 1/2" x 34 1/2"



## Patient Throw

- Soft, supple, durable throws are waterproof, stain-resistant, and washable
- Velcro closure adjusts to fit snugly and comfortably
- Reinforced towel clips attached to collar
- Available in nylon or plastic material, selected lengths, and colors

### Nylon

Ref #	Description	Size	Color
612	Knee length	27" x 44"	blue/yellow
613	Waist length	27" x 30"	blue/yellow/beige

### Plastic

Ref #	Description	Size	Color
850	Ankle length	36" x 44"	blue/green/clear
851	Knee length	27" x 44"	blue/clear
852	Waist length	27" x 30"	blue/clear





### Hand Mirror

- Durable plastic, double-sided patient hand mirror
- Standard mirror on one side, 2.5X magnifying mirror on the flip side for close-up viewing

Ref #	Description
45	5" diameter



### Disposable Mouth Mirrors Bulk

- Highly reflective light-weight, single-use mouth mirrors for examinations and soft tissue retraction procedures.
- **Improved** wider handle design provides slip-free and comfortable grip
- Clear, distortion-free, 7/8" diameter mirrors
- Ultra-thin reflecting surface for enhanced imaging
- Ideal for high-volume situations where interim sterilization is not practical

Ref #	Description	Quantity
93	Bulk Pack	72/box



### Disposable Mouth Mirrors Individually Wrapped

- Same great features and benefits as our bulk pack
- Ideal for distribution to patients for in home oral care procedures

Ref #	Description	Quantity
94	Individually Wrapped	24/box

## Disposable Mouth Mirrors

**Wider handle design**  
provides slip-free and comfortable grip

**INDIVIDUALLY WRAPPED**  
Ideal for high-volume situations where interim sterilization is not practical and for patients distribution for in home oral care procedures

6"

**NEW**



### Once-A-Day One Second Spray

- For use in air turbine hand pieces, prophy angles, and contra angles
- Also ideal as a lubricant or cleaner for hand pieces
- Contains hexane, non-corrosive

Ref #	Description
702	8.8 oz. aerosol can with extension tube

## TERMS AND CONDITIONS



**ORDERS AND QUOTATIONS:** Contact Palmero Healthcare for a quotation at 1-800-344-6424. All orders may be placed via fax 203-377-8988. Visit our website at [www.palmerohealth.com](http://www.palmerohealth.com)



**MINIMUM INVOICE:** \$50.00 minimum invoice order amount. There will be a \$10.00 handling charge for all orders not meeting the minimum \$50.00 order requirement.



**SHIPPING TERMS:** All shipments are F.O.B. Stratford, CT. Drop shipments will be made prepaid, unless otherwise specified, with applicable transportation charges added to the invoice. Common carrier shipments will be sent freight collect, and U.P.S. or Federal Express shipments will be prepaid and added to the invoice. If no shipping instructions are specified with the order, our method of transportation will be deemed acceptable. Palmero reserves the right to make partial shipments, unless specifically requested otherwise. If a special request is made, an additional handling charge will be included on the invoice to cover the resulting processing and carrier charges.



**PAYMENT TERMS:** All invoices are payable net 30 days from date of invoice. Palmero reserves the right to refuse orders, suspend credit, and cancel unfulfilled orders when credit conditions of a customer warrant such action.

**ACCEPTANCE:** All orders are subject to Palmero's approval. Palmero Healthcare reserves the right to alter the design or form of any product or discontinue any product prior to the date of delivery without incurring any obligation.

### SHIPMENTS, DAMAGES, AND SHORTAGES

**IN TRANSIT:** When signing for receipt of goods, all visible damages and/or carton shortages must be indicated on the carrier delivery receipt. The customer shall inspect all items upon arrival and submit in writing any claim of a shortage or non-conformance within 7 days after receipt of shipment. Responsibility for an error in shipment passes from the shipper to the consignee if the merchandise, or any portion thereof, is used prior to notification.



**MERCHANDISE RETURN POLICY:** A Return Authorization Number (RAN) must be requested prior to returning any merchandise to Palmero Healthcare. Issue of a RAN should not be interpreted as a final credit to the customer account. Palmero reserves the right to evaluate incoming returns prior

to issuing a credit to the customer. Merchandise must be returned in resalable condition (unless originally defective); in original complete packages (no partial packages accepted for credit); forwarding shipping costs are prepaid; a handling/restocking fee of 20% to 30% will be assessed to merchandise which is not defective. Percentage will be based on the condition of the merchandise returned. Non-returnable products received by Palmero will be returned to the customer and shipping costs will be billed to the customer. Discontinued Items by PHC will not be accepted for return. Specialty products such as X-Ray Aprons are date sensitive and any returns received beyond 9 months of the stamped date will not be issued credit. The X-Ray Aprons will be disposed of in the proper recycling method as specified by Palmero Healthcare's policy.

**PRODUCT WARRANTY:** Palmero's products are guaranteed to be free of material defects. Palmero is committed to a high standard of quality and craftsmanship. This standard will be employed in the manufacture of each product. If our products do not meet these standards, our customer may return the product to be replaced, at our discretion, and at no charge.

**WARRANTY STIPULATIONS:** Palmero Healthcare warrants that it has good title to the products sold to the buyer; and that its goods will conform to buyer's specifications, if any. Except as otherwise provided herein, all other warranties, whether expressed or implied, and including, but not limited to, the warranty for merchantability and fitness for a particular purpose are hereby excluded. In all events, buyer's sole and exclusive remedy shall be limited to the amount of the purchase price of any products on which a claim is made as its exclusive remedy for breach of the warranty expressed herein and for negligence and otherwise. Under no circumstances shall Palmero be liable to buyer or any other person for consequential or incidental damages.

**CHOICE OF LAW:** All disputes arising out of all sales by Palmero to buyer shall be governed by the State of Connecticut and litigated in Fairfield County, Connecticut. Buyer hereby consents to the jurisdiction of the courts residing in Fairfield County, Connecticut.

### PALMERO HEALTHCARE LLC

120 GOODWIN PLACE

STRATFORD, CT 06615

[WWW.PALMEROHEALTH.COM](http://WWW.PALMEROHEALTH.COM)

1-800-344-6424 OR 203-377-6424







# 2023/2024 PRODUCT GUIDE

palmerohealth.com / 1-800-344-6424



Palmero is now a proud member of HuFriedyGroup

