



VACUUM SHOCK™ & CLEAN™



From Screen to Clean. The Two-Step Evacuation System Maintenance Too Simple to Ignore.

Keeping your evacuation line clean is an important part of keeping your HVE and saliva ejector running smoothly. Vacuum Shock and Vacuum Clean are an easy two-step system that keeps high/low speed lines clean and fresh.

Vacuum Shock breaks up any accumulated debris while Vacuum Clean keeps your lines fresh and functioning at peak efficiency. When used in conjunction, the time-release formulation ensures your line is protected around the clock. Don't let a compromised evacuation system create problems for your practice, use the one-two punch that keeps evacuation lines happy.



Our Evacuation Product Line follows Best Management Practices for waterline management by EPA Guidelines. These products do not contain bleach, chlorine, iodine, peroxide, or any oxidizers.

	Vacuum Shock	Vacuum Clean
Products:		
Two-Step of Evacuation System Maintenance Kit:	Stage 1 1 tablet monthly	Stage 2 1 tablet weekly
Features & Benefits:	<ul style="list-style-type: none"> Keeps vacuum lines clean and extends pump life 	<ul style="list-style-type: none"> Restores pressure, keeps vacuum lines clean, extends pump life Self-activating tablet releases powerful cleaning enzymes through the week Allows consistent vacuum suction and asepsis
Alternative Use:	n/a	Suction Lines

VACUUM SHOCK™ & CLEAN™

An easy 2-Step Evacuation System Maintenance that keeps vacuum lines functioning at peak efficiency

Don't let a compromised evacuation system make problems for your practice, use the one-two punch that keeps the evacuation lines happy!



Find the Fit for Your Practice



Vacuum Shock™

Ref# 3546, 6 tablets/jar

- First stage in Palmero 2-Step Evacuation System Maintenance Kit
- 1 tablet monthly



Vacuum Clean™

Ref# 3547, 45 tablets/jar

- Second stage in Palmero 2-Step Evacuation System Maintenance Kit
- 1 tablet weekly



E-VAC System Set™

Ref# 3547S, 1 Set

- Convenient all-in-one Vacuum Shock and Vacuum Clean Evacuation System Maintenance Kit



Shock & Clean Operatory and Trial Kits

Ref 3546O, 1 Kit

Ref 3546TK, 1 Trial Tube

- No fuss easy to use operatory kits premeasured for 16 day schedule
- Convenient packages for operatory storage, kit includes 6 tubes, each containing 16 Vacuum Clean tablets and 1 Vacuum Shock tablet, also available in a single trial tube

For Optimal Results, use these evacuation maintenance products...



Shock & Clean Starter Kit

Ref# 3546K, 1 Set

- The one-two punch for evacuation system maintenance in one convenient kit
- Kit includes 1 jar of Vacuum Shock, 1 jar of Vacuum Clean and 1 Pour and Clean Bottle for concentrate premixing



Pour and Clean Bottle

Ref 3547P, 1-16 oz. Container

- Measuring container to prepare rinse treatment for the HVE and low-speed lines
- Convenient pre-measured compartments clearly marked for easy mixing and pouring of solution



Saliva Ejector Screens

Ref 1602, 25-Pack

Ref 1602A, 50-Pack

Ref 1602C, 75-Pack

Ref 1602B, 100-Pack

Ref 1602D, 250-Pack

- Disposable and aseptic
- Keeps debris from being separated and lodged in the lines
- Universal fit



Evac Screens

Ref 1805 (blue), 144 screens/bag

Ref 1807 (green), 144 screens/bag

- Helps avoid contact with potentially hazardous material
- Keeps debris from lodging in evacuation lines
- Disposable screens color-coded by size and fit
- Available in two sizes
 - Blue 2 1/8" A-DEC with plastic canisters
 - Green 1 7/8" Belmont, Biotec, Chayes, Dentech, DentalEZ, Pelton/Crane, Proma and others



For more information, visit palmerohealth.com, call 800.344.6424 or email customerservice@palmerohealth.com.