

PROVISION® NEW SAFETY EYEWEAR

The Insight to Safeguard Your Eyesight

Palmero's Safety Eyewear is a sight to behold. From broad protective products consisting of full face shields and goggles, to more procedure specific systems like laser safety and bonding glasses, our protective eyewear offers the ultimate in comfort, convenience, and affordability. Our safety eyewear meets the ANSI Z87.1 eye and face protection standards with UVA and UVB protection. Give your eyesight the protection it deserves with Palmero Healthcare Safety Eyewear.

Implement an Eye Conservation Program for your Practice

Select safety eyewear that serves dental professionals and patients in your practice on multiple levels to: eliminate blue light exposure; protect against spray, splash and splatter for infection control; protect from high impact flying projectiles that can cause permanent damage to the eye and are comfortable enough to ensure compliance.

The Four C's of Safety Eyewear

Certified + Clarity + Comfort = Compliance



CERTIFIED Meets industry standards and guidelines for lens, cornea, and retina protection i.e., ANSI Z87.1, OSHA and CDC. Wraparound design for full coverage. Manufactured from high impact resistant materials to protect from flying debris and projectiles. Specialty bonding lens filters out damaging Blue Light 400nm-500nm levels to the retina.



CLARITY High Excellent optics for visual acuity provided by anti-fog coating, scratch resistant technologies, and superior lens designs.



COMFORT Broad range of frame, temple, and nose bridge styles to provide an individualized fit for all day security and optimal comfort.



COMPLIANCE Willingness and discipline to use safety eyewear for everyone, in every procedure, every day.



Eyewear Collections at a Glance

Traditional Safety Eyewear Meets ANSI Z87 Standards, Does not contain latex						Disposables Does not contain latex	
Wraparound	Goggles	Bifocals	Bonding	Magnification	Laser	Face Shields	Eyewear
<ul style="list-style-type: none"> Ergonomic lightweight Full & minis Ultra anti-fog 	<ul style="list-style-type: none"> Wear over prescription glasses Side shields 	<ul style="list-style-type: none"> Up to +3.0 diopter 	<ul style="list-style-type: none"> Red and orange lens 	<ul style="list-style-type: none"> Full lens magnification 	<ul style="list-style-type: none"> Yellow and green lens 	<ul style="list-style-type: none"> Full & half sizes Foam elastic band or snap-on style 	<ul style="list-style-type: none"> Snap-on feature One frame, multicolored lens options



Protecting Clinicians, Patients, and the Practice

PROVISION[®] SAFETY EYEWEAR



Introducing Palmero's newest safety eyewear to enhance protection and comfort for your practice, staff, and patients!

Wraparound Collection

ProVision[®] Allure[™]

- Virtually unbreakable Pink or Lavender ultra lightweight frame (24 grams)
- Soft ventilated nose pad for a secure, comfortable fit
- Curved lens for wraparound splatter protection
- Fits small/narrow and medium face shapes
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref#	Frame/Lens	Eyewear Type	Facial Fit
3604PC	Pink/Clear	Clinician/Patient	Small/Narrow/Medium
3604LC	Lavender/Clear	Clinician/Patient	Small/Narrow/Medium
3604LG	Lavender/Grey	Patient	Small/Narrow/Medium



Clear/grey lens option

Ultra lightweight frame

Soft ventilated nose pad

ProVision[®] Air[™]

- Ergonomically designed temples fit the contour of your head for a comfortable, secure, non-slip fit
- Flexible, anti-static, lightweight frame (15 grams)
- Extended lens coverage for great protection and style
- Distortion-free, wraparound polycarbonate lens offers superior peripheral vision and protection
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref#	Frame/Lens	Eyewear Type	Facial Fit
3608C	Clear/Clear	Clinician/Patient	Small/Narrow/Medium
3608G	Grey/Grey	Patient	Small/Narrow/Medium



Clear/grey lens option

Ultra lightweight Frame

Integrated nose bridge

ProVision[®] Element[™]

- Ergonomically designed temples fit the contour of your head for a comfortable, secure, non-slip fit
- Flexible, lightweight frame (20 grams)
- Extended lens coverage for great protection and style
- Distortion-free wraparound polycarbonate lens offers superior peripheral vision and protection
- Wider frame design allows for broader facial fit
- Anti-static and anti-fog lens
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref#	Frame/Lens	Eyewear Type	Facial Fit
3606C	Clear/Clear	Clinician/Patient	Medium/Large
3606G	Grey/Grey	Patient	Medium/Large



Anti-fog lens in clear and grey

Ultra lightweight & wider frame for broader facial fit

Integrated nose bridge



Protecting Clinicians, Patients, and the Practice

PROVISION® SAFETY EYEWEAR



Introducing Palmero's newest safety eyewear to enhance protection and comfort for your practice, staff, and patients!

Wraparound Collection

ProVision® Clarity™

- Lightweight, frameless, flexible, translucent wrap (24 grams)
- Firm nose pad and soft, rubberized temple tips for comfort and grip
- Curved lens for wraparound splatter protection
- Enhanced with **Ultra Anti-Fog Technology** to provide incredibly clear and uninterrupted vision
- Fits small/narrow and medium face shape
- Blocks up to 99.9% UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref#	Frame/Lens/Tips	Type	Facial Fit
3605L	Clear/Clear/Lavender	Clinician/Patient	Small/Narrow/Medium
3605B	Clear/Clear/Navy	Clinician/Patient	Small/Narrow/Medium



Soft rubber color tips



Ultra anti-fog coating technology



Integrated, comfortable nose bridge

ProVision® Infinity™

- Upscale black frame and clear wraparound lens with anti-fog and anti-static features
- Top vented silicone brow guard provides added protection and venting to reduce fogging
- Enhanced with pantoscopic tilt feature provided with telescopic and ratcheting temple arms allowing for customizable fit and optimal comfort
- Soft micro ribbed nose bridge grip for secure, non-slip fit
- Fits broader facial structures
- Detachable and replacement lens in clear and bonding
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Detachable & replaceable lens



Soft rubber ribbed nose bridge



Soft rubber tips




Telescopic & ratcheting temples




Top vented silicone frame

Pantoscopic Tilt


The angle of the lenses as they sit on the clinician's face



Natural Head Position



Ear Lower



Ear Higher

Ref#	Description	Frame/Lens	Type	Facial Fit
3613C	ProVision Infinity	Black/Clear	Clinician	Universal
3613R	Replacement Lens	Clear Lens	Clinician	Universal
3613BR	Replacement Lens	Red Lens	Clinician	Universal



Protecting Clinicians, Patients, and the Practice

PROVISION[®] SAFETY EYEWEAR



Introducing Palmero's newest safety eyewear to enhance protection and comfort for your practice, staff, and patients!

Wraparound Collection

ProVision[®] Econo Tilts[™]


- Lightweight, polycarbonate frame with tilting arms ensures a great contoured fit
- Firm nose pad and soft, rubberized temple tips guarantee comfort and grip
- Curved lens provides wraparound splatter protection
- Anti-fog polycarbonate lens
- Available in 2 colors (Rose Quartz and Blue)
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref#	Frame/Lens/Tips	Type	Facial Fit
3614CB	Clear/Clear/Blue	Clinician/Patient	Universal
3614CQ	Clear/Clear/Rose Quartz	Clinician/Patient	Universal



Integrated comfortable nose bridge



Soft rubber temple tips




Ratcheting temples

ProVision[®] Mini Bonding Wraps[™]

- Mini version of our popular full size sleek, bonding eyewear with the same great features and benefits
- Tested and proven to protect blue light rays between 385-500 nm which exceeds the blue light wavelengths of 450-495 nm
- Red polycarbonate lenses for maximum protection
- Extended lens coverage for great protection and style
- Filters up to 95% of blue, UVA, and UVB rays
- Comfortable, lightweight, fog-free, and scratch-resistant
- Meets ANSI Z87.1 standard



Ref#	Frame/Lens	Eyewear Type	Facial Fit
3607B	Red/Red	Patient	Small/Narrow/Medium



Bonding red lens



Ultra lightweight frame



Integrated nose bridge

For more information, visit palmerohealth.com, call 800.344.6424 or email customerservice@palmerohealth.com.



Protecting Clinicians, Patients, and the Practice

PROVISION® SAFETY EYEWEAR



Introducing Palmero's newest safety eyewear to enhance protection and comfort for your practice, staff, and patients!

Wraparound Collection

ProVision® Rainbow™ Mini

- Put your younger patients at ease with our brightly colored and economical patient eyewear
- Flexible temples in 6 colors make it easy to have your small patients wear their safety glasses
- Reduce the effect of harsh operatory lights with scratch-resistant, 100% polycarbonate wraparound grey lens
- The Mini is specifically designed to fit smaller faces including pediatric patients
- Blocks up to 99.9% of UVA and UVB rays
- Scratch resistant lens designed to meet ANSI Z87.1 Safety Standards of impact protection
- Sold in a box of 12; 2 of each color (blue, orange, green, purple, yellow & pink)



Grey lens

Ultra lightweight frame

Integrated nose bridge

Ref#	Frame/Lens	Eyewear Type	Facial Fit
3607GB	6 colors/grey	Patient	Small/Narrow/Medium

ProVision® Flexiwrap™

- Lightweight, innovative spring-flex temple design with extreme flexibility to fit any facial structure
- Short temples for a secure comfortable fit and **great for patients in the supine position**
- Soft touch temple and soft nose bridge
- Full wraparound protection
- Blocks up to 99.9% UVA and UVB rays
- Meets ANSI Z87.1 standard



Short temples for secure & broader facial fit



Integrated comfortable nose bridge



Soft temple tips

Ref#	Frame/Tips/Lens	Type	Facial Fit
3609GG	Black/Grey/Grey	Patient	Universal
3609GC	Black/Grey/Clear	Clinician/Patient	Universal
3609OG	Black/Orange/Grey	Patient	Universal
3609OC	Black/Orange/Clear	Clinician/Patient	Universal

For more information, visit palmerohealth.com, call 800.344.6424 or email customerservice@palmerohealth.com.



Protecting Clinicians, Patients, and the Practice

PROVISION® SAFETY EYEWEAR



Introducing Palmero's newest safety eyewear to enhance protection and comfort for your practice, staff, and patients!

Wraparound Collection

ProVision® Secure™ with Strap

- Distortion-free wraparound anti-fog and anti-static polycarbonate lens offers superior peripheral vision and protection
- Vent design at top and bottom of foam insulation allows ventilation for comfort while protecting against spray, splashes and splatter
- Detachable, adjustable strap firmly secures eyewear in place to minimize slippage and provide all day comfort
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref#	Frame/Tips/Lens	Type	Facial Fit
3630GC	Black/Green/Clear	Clinician	Universal



Vent design and anti-fog and anti-static polycarbonate lens



Foam insulation for comfort and protection



Detachable, adjustable strap for security

Goggles Collection

ProVision® Eyesaver Sleeks™

- Sleek, ultra-light, economical goggles
- Adjustable temples allow for broader facial fit
- Designed to fit over prescription eyewear
- All around integrated side shields protection
- Scratch-resistant
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref#	Frame/Lens	Type	Facial Fit
3SLK	Grey/Grey	Patient	Universal
16SLK	Red bonding frame/lens	Clinician/Patient	Universal
17SLK	Clear/Clear	Clinician/Patient	Universal
18SLK	Green/Green	Patient	Universal



Adjustable temples to fit broader facial



Fully integrated top, bottom & side shields protection



Built-in comfortable nose bridge

For more information, visit palmerohealth.com, call 800.344.6424 or email customerservice@palmerohealth.com.



Protecting Clinicians, Patients, and the Practice

PROVISION® SAFETY EYEWEAR



Introducing Palmero's newest safety eyewear to enhance protection and comfort for your practice, staff, and patients!

Bifocals Collection

ProVision® EZ-Focals™

- Lightweight, super comfortable wraparound frame (24 grams)
- Contoured temples with soft tips for added comfort, grip, and a secure fit
- Soft, self-fitting nose bridge, minimizes slippage
- Available in +1.5, +2.0, +2.5, +3.0 diopters
- Clear anti-fog and anti-static lens
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref#	Diopter	Frame/Lens	Type	Facial Fit
3770B	+1.5	Black/Clear	Clinician	Universal
3770C	+2.0	Black/Clear	Clinician	Universal
3770D	+2.5	Black/Clear	Clinician	Universal
3770E	+3.0	Black/Clear	Clinician	Universal



Integrated comfortable nose bridge



Soft rubber temple tips



Four available diopters

Magnification Collection

ProVision® Econo Loupes

- Full 50mm diameter magnifying lens for a fully-enhanced clinical field
- Economic, lightweight alternative to traditional loupes
- Reduces eye strain and increases visual precision
- Wraparound frame and soft flexible nose bridge to minimize slippage and provide all day comfort
- Soft tip flexible temples ensure a secure fit
- Available in +1.5, +2.0 and +2.5 diopters
- Blocks up to 99.9% of UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref#	Diopter	Frame/Lens	Type	Facial Fit
3780B	+1.5	Black/Clear	Clinician	Universal
3780C	+2.0	Black/Clear	Clinician	Universal
3780D	+2.5	Black/Clear	Clinician	Universal



Full 50mm diameter magnifying lens



Soft rubber temple tips



Integrated comfortable nose bridge

For more information, visit palmerohealth.com, call 800.344.6424 or email customerservice@palmerohealth.com.



Protecting Clinicians, Patients, and the Practice

PROVISION[®] EYEWEAR MAINTENANCE

Eyewear Maintenance

ProVision[®] Crystal Clear Anti-Fog Solo Wipes

- Cleans lenses while providing long-lasting, anti-fog protection
- Safe and effective for use on lenses including polycarbonate, acrylic, and glass
- Dispenser box contains 100 premoistened individually wrapped towelettes
- Anti-fog, streak-free, fresh scent
- Great for eyewear and goggles

Ref#	Description	Quantity
3534S	ProVision Crystal Clear Anti-Fog Solo Wipes 8" x 5"	100-pack box



ProVision[®] Crystal Clear Optical Cleaning Wipes

- Convenient wipes premoistened with anti-fog, anti-static lens cleaning, and polishing solution
- Cleans dirt, grease, and debris from eyewear
- Container/lid does not contain latex

Ref#	Description	Quantity
3534	ProVision Crystal Clear Optical Cleaning Wipes 6" x 6 3/4"	160 count canister



ProVision[®] Anti-Fog Solution

- Compatible with all eyewear, face shields, mirrors, loupes, and scopes
- Resists scratches, static, and stains
- Foam formula yields 5 times the applications compared to a traditional pump or squeeze bottle
- Alcohol-free, non-toxic, and non-flammable

Ref#	Description	Product Qty	Case Qty
3533	ProVision Anti-Fog Solution 8 oz.	1 bottle	16 bottles
3536	ProVision Anti-Fog Solution 1.7 oz.	1 bottle	24 bottles



ProVision[®] Lens Cleaning Station

- Convenient eyewear and lens cleaning station
- Designed for countertop placement or wall mounting with hanging hardware
- Station includes 8 oz. pump spray bottle of optical cleaning solution and 600 dry wipes

Ref#	Description	Product Qty
3535	ProVision Lens Cleaning Station 7 3/4" x 8" x 5"	1 station



For more information, visit palmerohealth.com, call 800.344.6424 or email customerservice@palmerohealth.com.



Protecting Clinicians, Patients, and the Practice