

PROVISION[®] MINI SAFETY EYEWEAR

The Insight to Safeguard Your Eyesight

Palmero's Safety Eyewear is a sight to behold. From broad protective products consisting of full face shields and goggles, to more procedure specific systems like laser safety and bonding glasses, our protective eyewear offers the ultimate in comfort, convenience, and affordability. Our safety eyewear meets the ANSI Z87.1 eye and face protection standards with UVA and UVB protection. Give your eyesight the protection it deserves with Palmero Healthcare Safety Eyewear. Mini specifically designed to fit smaller faces including pediatric patients.

Implement an Eye Conservation Program for your Practice

Select safety eyewear that serves dental professionals and patients in your practice on multiple levels to: eliminate blue light exposure; protect against spray, splash and splatter for infection control; protect from high impact flying projectiles that can cause permanent damage to the eye and are comfortable enough to insure compliance.

The Four C's of Safety Eyewear

Certified + Clarity + Comfort = Compliance

CERTIFIED Meets industry standards and guidelines for lens, cornea, and retina protection i.e., ANSI.Z87.1, OSHA and CDC. Wraparound design for full coverage. Manufactured from high impact resistant materials to protect from flying debris and projectiles. Specialty bonding lens filters out damaging Blue Light 400nm-500nm levels to the retina.

CLARITY High Excellent optics for visual acuity provided by anti-fog coating, scratch resistant technologies, and superior lens designs.

COMFORT Broad range of frame, temple, and nose bridge styles to provide an individualized fit for all day security and optimal comfort.

COMPLIANCE Willingness and discipline to use safety eyewear for everyone, in every procedure, every day.



Eyewear Collections at a Glance

Traditional Safety Eyewear Meets ANSI.Z87 Standards, Does not contain latex						Disposables Does not contain latex	
Wraparound	Goggles	Bifocals	Bonding	Magnification	Laser	Face Shields	Eyewear
<ul style="list-style-type: none"> Ergonomic lightweight Full & minis Ultra anti-fog 	<ul style="list-style-type: none"> Wear over prescription glasses Side shields 	<ul style="list-style-type: none"> Up to +3.0 diopter 	<ul style="list-style-type: none"> Red and orange lens 	<ul style="list-style-type: none"> Full lens magnification 	<ul style="list-style-type: none"> Yellow and green lens 	<ul style="list-style-type: none"> Full & half sizes Foam elastic band or snap-on style 	<ul style="list-style-type: none"> Snap-on feature One frame, multicolored lens options

PROVISION® MINI SAFETY EYEWEAR

Introducing Palmero's newest safety eyewear to enhance protection and comfort for your pediatric patients!

Wraparound Collection

Available
Jan-2019



ProVision® Rainbow™ Mini

- Put your younger patients at ease with our brightly colored and economical patient eyewear
- Flexible temples in 6 colors make it easy to have your small patients wear their safety glasses
- Reduce the effect of harsh operatory lights with scratch-resistant, 100% polycarbonate wraparound grey lens
- The Mini is specifically designed to fit smaller faces including pediatric patients
- Blocks up to 99.9% of UVA and UVB rays
- Scratch Resistant lens designed to meet ANZI Z787.1 Safety Standards of impact protection
- Sold in a box of 12; 2 of each color (blue, orange, green, purple, yellow & pink)



Grey lens



Ultra lightweight
frame



Integrated
nose bridge



Ref#	Frame/Lens	Eyewear Type	Facial Fit
3607GB	6 colors/grey	Patient	Small/Narrow/Medium

ProVision® Mini Bonding Wraps™

- Mini version of our popular full size sleek, bonding eyewear with the same great features and benefits
- Tested and proven to protect blue light rays between 385-500 nm which exceeds the blue light wavelengths of 450-495 nm
- Red polycarbonate lenses for maximum protection
- Extended lens coverage for great protection and style
- Filters up to 95% of blue, UVA, and UVB rays
- Comfortable, lightweight, fog-free, and scratch-resistant
- Meets ANSI Z87.1 standard

Available
Jan-2019



Bonding
red lens



Ultra lightweight
frame



Integrated
nose bridge

Ref#	Frame/Lens	Eyewear Type	Facial Fit
3607B	Red/Red	Patient	Small/Narrow/Medium

PROVISION® MINI SAFETY EYEWEAR

Palmero offers a full line of safety eyewear specifically designed to protect your pediatric patients!

Wraparound Collection

ProVision® Flexiwrap™

- Lightweight, innovative spring-flex temple design with extreme flexibility to fit any facial structure
- Short temples for a secure comfortable fit
- Soft touch temple and soft nose bridge
- Full wraparound protection
- Blocks up to 99.9% UVA and UVB rays
- Meets ANSI Z87.1 standard



Ref#	Frame/Tips/Lens	Type	Facial Fit
3609GG	Black/Grey/Grey	Patient	Universal
3609GC	Black/Grey/Clear	Clinician/Patient	Universal
3609OG	Black/Orange/Grey	Patient	Universal
3609OC	Black/Orange/Clear	Clinician/Patient	Universal



Short temples for secure & broader facial fit



Integrated comfortable nose bridge



Soft temple tips

ProVision® Mini Econo Wraps™

- Same great features and benefits as full size Econo Wraps
- 25% smaller, ideal for kids and adult with smaller/narrow facial structure
- Sleek, featherweight polycarbonate wraps, universal fit
- Fog-free, scratch-resistant
- Available in 2 colors
- Blocks up to 95% of UVA and UVB rays
- Meets the ANSI Z87.1 standard



Ref#	Frame/Lens	Eyewear Type	Facial Fit
3607G	Grey/Grey	Patient	Small/Narrow/Medium
3607C	Clear/Clear	Clinician/Patient	Small/Narrow/Medium



Clear/grey lens option



Ultra lightweight frame



Integrated nose bridge


ProVision® Miniwrap™

- Scaled-down version of a popular style
- Adjustable temples with integrated comfortable nose bridge
- One-piece polycarbonate lens with integral side shields wraps around for wide field of vision and complete protection
- Blocks up to 95% of UVA and UVB rays
- Meets ANSI Z87.1 standard

Ref#	Frame/Lens	Eyewear Type	Facial Fit
3705P	Pink/Clear	Clinician/Patient	Child/Youth



Integrated side shields



Integrated nose bridge



Telescopic temple

PROVISION® EYEWEAR MAINTENANCE

Eyewear Maintenance

ProVision® Crystal Clear Anti-Fog Solo Wipes

- Cleans lenses while providing long-lasting, anti-fog protection
- Safe and effective for use on uncoated lenses including polycarbonate, acrylic, and glass
- Dispenser box contains 100 premoistened individually wrapped towelettes
- Anti-fog, streak-free, fresh scent
- Great for eyewear, goggles, monitors, and tablets

Ref#	Description	Quantity
3534S	ProVision™ Crystal Clear Anti-Fog Solo Wipes 8" x 5"	100-pack box

ProVision® Crystal Clear Optical Cleaning Wipes

- Convenient wipes premoistened with anti-fog, anti-static lens cleaning, and polishing solution
- Cleans dirt, grease, and debris from eyewear
- Container/lid does not contain latex

Ref#	Description	Quantity
3534	ProVision Crystal Clear Optical Cleaning Wipes 6" x 6 ¾"	160 count canister

ProVision® Anti-Fog Solution

- Compatible with all eyewear, face shields, mirrors, loupes, and scopes
- Resists scratches, static, and stains
- Foam formula yields 5 times the applications compared to a traditional pump or squeeze bottle
- Alcohol-free, non-toxic, and non-flammable

Ref#	Description	Product Qty	Case Qty
3533	ProVision Anti-Fog Solution 8 oz.	1 bottle	16 bottles
3536	ProVision Anti-Fog Solution 1.7 oz.	1 bottle	24 bottles

ProVision® Lens Cleaning Station

- Convenient eyewear and lens cleaning station
- Designed for countertop placement or wall mounting with hanging hardware
- Station includes 8 oz. pump spray bottle of optical cleaning solution and 600 dry wipes

Ref#	Description	Product Qty
3535	ProVision Lens Cleaning Station 7 ¾" x 8" x 5"	1 station



Coming Soon!

