

DYNAMIC DISPOSABLES™

The New Face of Safety and Compliance Disposable Face Shields with a Reputation for Protection



The CDC suggests the use of a face shield when additional protection is needed to safeguard the entire face from infectious agents.

Dynamic Disposables Face Shield is a lightweight, disposable, cost effective product that provides extremely comfortable coverage. As an addition to our traditional line, we are proud to introduce the next generation of face shields with an easy-to-use snap-on disposable lens feature.

Unlike some systems that use an 'eye glass style' frame to support the shield, we deliver a more-comfortable, more-secure fit with our light weight frame which does not require a nose bridge. The full and half shields both use our Ultra-Clear™ Technology to prevent fogging, glare, and static electric buildup. The crystal clear plastic allows for excellent front and peripheral vision with protection you can trust.

Choose your protection from our original snug-fitting elastic band design or Palmero's new Snapeez Face Shields!



	Traditional Face Shields	NEW Snapeez Face Shields
PRODUCTS:		
UNIQUE FEATURES:	Foam Design Soft, stretch elastic strap for a secure and comfortable fit	Wraparound Frame Ergonomic frame wraps around the head for a more secure fit, designed to be worn with or without prescription eyewear, bridge-free nose means more comfortable fit Easy Snap On Attachment Self-sealing face shields prevent exposure to bloodborne pathogens Economical Reusable Frame Reduce the cost per use
COMMON FEATURES:	Easy-to-Use and Convenient Disposable & Cost Effective Cost effective face shields fit over eye glasses and face masks Ultra-Clear Technology Anti-fog, anti-static, and optically clear Reliability & Safety Meets OSHA bloodborne pathogen guidelines Does not contain latex	
LENGTH:		
Full Face Shield	8"	7.25"
Half Face Shield	6"	5.5"

DYNAMIC DISPOSABLES™

Take Comfort in Compliance


Our disposable full and half face shields offer the ultimate in PPE protection. The comfortable elastic fastener and ergonomic frame designs keep the lightweight shield secure without putting undo pressure on the nose or face. Surround frames with snap on shields are now also available. Both designs are extremely easy to use and remove, and fit comfortably over glasses or goggles. Because they are disposable, your shield is guaranteed to be scratch free for optimal visibility. For the ultimate in compliance, comfort, and convenience, Palmero Healthcare face shields clearly deliver superior protection for a stress free procedure.

Find the Fit for Your Practice

Traditional Full and Half Face Shields

- Easy-to-use. Cost Effective. Disposable.
- Integrated soft, stretch elastic stripe for a secure and comfortable fit.

REF#	DESCRIPTION	QUANTITY
1964	Disposable Full Face Shield	24 full face shields
1966	Disposable Half Face Shield	24 half face shields



Snapeez Face Shields Office Pack NEW

- Save money with bulk ordering.
- Includes latex free frame, face shields, and cord.

REF#	DESCRIPTION	QUANTITY
3918	Office Pack, Snapeez Full Face Shield	10 black frames, 20 full face shields, 1 cord
3920	Office Pack, Snapeez Half Face Shield	10 black frames, 20 half face shields, 1 cord



Snapeez Face Shield Replacement Shields Clear NEW

- Ultra-Clear™ Technology eliminates fog, glare, and static electricity.
- Works with reusable frame to save money.

REF#	DESCRIPTION	QUANTITY
3919	Clear full face shields	25 pack
3921	Clear half face shields	25 pack



Replacement Frames

- Ergonomic design eliminates pressure points.
- Wrap around design won't slide off head.

REF#	DESCRIPTION	QUANTITY
3901	Black Frames	10 pack



Accessories

- Can be worn over prescription eyewear.
- Allows eyewear to hang conveniently around neck when not in use.

REF#	DESCRIPTION	QUANTITY
3904	22" Classic braided cords	5 pack, assorted
3913	24" Plastic cords	10 pack, black



For more information, visit palmerohealth.com, call 800.344.6424 or email customerservice@palmerohealth.com.