

# CLING SHIELD® X-RAY APRONS

## Introducing Palmero's New Smiley Face Design!



Palmero's Cling Shield® X-Ray product line is designed to protect your staff and patients from harmful radiation.

Our PROTECTALL™ x-ray apron offers a built-in thyroid neck collar to optimally protect reproductive organs and thyroid gland. Each of our x-ray aprons and neck collars are made in the U.S.A.

We're excited to introduce our new smiley face design for our x-ray aprons! Young and old will recognize this fun design that has become so popular with social media. The new design has the same great quality and features as our classic smiley face. So whether you're new school or old school, both designs are available for your selection.

### Features at a Glance

- FDA registered, medical-grade textured vinyl Cling Shield® backing that clings in place ensuring coverage by preventing slippage and providing stability
- High density lead and Palmero's Lead Free materials effectively shield radiation
- Made from easy to clean vinyl material that can be sanitized using **DisCide® Ultra** surface disinfectant and clean and preserve your investment with **Protectall™ Upholstery Cleaner Wipes**
- Integrated nylon hanging loops allow proper storage to avoid damaging the linings
- Custom-made to order in the U.S.A.



| Cling Shield® X-Ray Collection |                 |                      |           |                      |          |            |            |            |                 |
|--------------------------------|-----------------|----------------------|-----------|----------------------|----------|------------|------------|------------|-----------------|
| Protectall With Collar         | Apron No Collar | Dual Apron No Collar | Lap Apron | Deluxe 3/4 Pano Dual | Pan Cape | Pano Stole | Tech Apron | Coat Apron | Thyroid Collar  |
|                                |                 |                      |           |                      |          |            |            |            |                 |
| Child and Adult                | Child and Adult | Child and Adult      | Adult     | Adult                | Adult    | Adult      | Adult      | Adult      | Child and Adult |

For more information, visit [palmerohealth.com](http://palmerohealth.com), call 800.344.6424 or email [customerservice@palmerohealth.com](mailto:customerservice@palmerohealth.com).



Protecting Clinicians, Patients, and the Practice

# CLING SHIELD® X-RAY APRONS

Find the right fit for your practice

## Protectall™ Apron with Neck Collar

- Neck collar permanently attached
- Covers from thyroid to below the gonadal area
- Sized for petite adults or children
- Velcro clasp on neck collar
- Loops attached for easy hanging and storage



| Ref#                   | Description | Thickness | Size                | Color                             | Material*   |
|------------------------|-------------|-----------|---------------------|-----------------------------------|-------------|
| 27USMILEY – <b>NEW</b> | Lead-Free   | .3MM      | 20" x 20" Child     | Slate Blue/Black Binding          | Vinyl       |
| 27UHAPPYFACE           | Lead-Free   | .3MM      | 20" x 20" Child     | Blue & White Stripes/Blue Binding | Vinyl       |
| 27SMILEY – <b>NEW</b>  | Lead-Lined  | .3MM      | 20" x 20" Child     | Slate Blue/Black Binding          | Vinyl       |
| 27HAPPYFACE            | Lead-Lined  | .3MM      | 20" x 20" Child     | Blue & White Stripes/Blue Binding | Vinyl       |
| 24U                    | Lead-Free   | .3MM      | 22 ¼" x 25 ½" Adult | Multiple Color Options            | Vinyl/Nylon |
| 24                     | Lead-Lined  | .3MM      | 22 ¼" x 25 ½" Adult | Multiple Color Options            | Vinyl/Nylon |

## Neck/Thyroid Collar

- For the protection of the thyroid area.
- Velcro closing, one size fits all
- Use in conjunction with Child or Adult Patient Apron

| Ref# | Description | Size               | Material*   |
|------|-------------|--------------------|-------------|
| 25CU | Lead-Free   | 17" x 4" Child     | Vinyl/Nylon |
| 25C  | Lead-Lined  | 17" x 4" Child     | Vinyl/Nylon |
| 25U  | Lead-Free   | 22 ¾" x 4 ¼" Adult | Vinyl/Nylon |
| 25   | Lead-Lined  | 22 ¾" x 4 ¼" Adult | Vinyl/Nylon |



## Hold-It™ Standard and Economy Apron Hangers

- Do not fold aprons as it can cause cracking within the lining. Use Palmero's assortment of hangers to protect your investment in safety
- White-coated steel hanger
- Compatible with adult and child Patient and Protectall™ style aprons
- Hanging hardware included



| Ref# | Description           | Size     | Color |
|------|-----------------------|----------|-------|
| 28   | Standard Apron Hanger | 15" x 2" | White |
| 31   | Economy Apron Hanger  | 15" x 2" | White |

## Go Green with our eco-friendly X-Ray Apron Recycling/Trade-In Program

### Buy 1, Recycle 1, Receive 1\*

1. **Purchase** a new Palmero x-ray apron.
2. **Collect** old or unwanted x-ray apron of any brand.
3. **Submit** invoice as proof of purchase to obtain and receive an RA# via Fax 203.377.8988 or email [customerservice@palmerohealth.com](mailto:customerservice@palmerohealth.com). Include ship to address, phone, and email.
4. **Contact** Customer Service for a Return of Authorization Number (RA#) to register for recycling.
5. **Ship** your old apron with RA# for responsible disposal to: Palmero Healthcare LLC Attn: X-Ray Recycle/Trade-In Program, 120 Goodwin Place, Stratford, CT 06615.
6. **Acquire** a new x-ray apron hanger. Once your recycled apron is received, we'll send you an \$20.00 retail value hanger for your participation!

\*Offer only available for continental U.S.



Buy 1, Recycle 1, Receive 1 apron hanger\*

