

# CLING SHIELD® X-RAY APRONS

## Palmero's Radiation Mitigation Wear

Palmero's Cling Shield™ X-Ray product line is designed to protect your staff and patients from harmful radiation.

Our PROTECTALL x-ray apron offers a built-in thyroid neck collar to optimally protect reproductive organs and thyroid gland. Each of our x-ray aprons and neck collars are made in the U.S.A.

We offer 10 unique styles and sizes from x-ray aprons to thyroid neck collars in a variety of colors and selected silk screen prints as well as two material types of Vinyl or Nylon and Lead-Free [.3MM] or Lead-Lined [.3MM].

### Features at a Glance

- Medical-grade textured vinyl Cling Shield™ backing that clings in place ensuring coverage by preventing slippage and providing stability
- High density lead and lead-free materials effectively shield radiation
- Made from easy to clean vinyl material that can be sanitized using DisCide® Ultra surface disinfectant
- Integrated nylon hanging loops allow proper storage to avoid damage of the lining
- Custom-made to order in the U.S.A.



### X-Ray Collections at a Glance

Cling Shield™ X-Ray Aprons and Collars									
Protectall With Collar	Apron No Collar	Dual Apron No Collar	Lap Apron	Deluxe 3/4 Pano Dual	Pan Cape	Pano Stole	Tech Apron	Coat Apron	Thyroid Collar
Child and Adult	Child and Adult	Child and Adult	Adult	Adult	Adult	Adult	Adult	Adult	Child and Adult

To view Palmero full x-ray product line, visit [www.palmerohealth.com](http://www.palmerohealth.com) and go to CATALOG on the navigation menu for online access or download version.

# CLING SHIELD® X-RAY APRONS

Find the right fit for your practice

## Protectall Apron with Neck Collar

- Neck collar permanently attached
- Covers from thyroid to below the gonadal area
- Velcro clasp on neck collar
- Loops attached for easy hanging and storage

Ref#	Description	Size	Material*
27U	Lead-Free	20" x 20" Child	Vinyl/Nylon
27	Lead-Lined	20" x 20" Child	Vinyl/Nylon
24U	Lead-Free	22 ¼" x 25 ½" Adult	Vinyl/Nylon
24	Lead-Lined	22 ¼" x 25 ½" Adult	Vinyl/Nylon



## Apron, No Collar

- Covers from over shoulders to below the gonadal area
- Loops attached for easy hanging and storage

Ref#	Description	Size	Material*
22U	Lead-Free	20" x 19 ½" Child	Vinyl/Nylon
22	Lead-Lined	20" x 19 ½" Child	Vinyl/Nylon
20U	Lead-Free	22 ¼" x 24 ¼" Adult	Vinyl/Nylon
20	Lead-Lined	22 ¼" x 24 ¼" Adult	Vinyl/Nylon



## Child Pano Dual Apron, No Collar

- Covers in front from over shoulders to below the gonadal area
- Back is the same length as front
- Velcro clasp at shoulder

Ref#	Description	Size	Material*
26CU	Lead-Free	19 7/8" x 19 ½" Child	Vinyl/Nylon
26C	Lead-Lined	19 7/8" x 19 ½" Child	Vinyl/Nylon
26U	Lead-Free	22 ½" x 26 ½" Adult	Vinyl/Nylon
26	Lead-Lined	22 ½" x 26 ½" Adult	Vinyl/Nylon



## Neck/Thyroid Collar

- For the protection of the thyroid area.
- Velcro closing, one size fits all
- Use in conjunction with Child or Adult Patient Apron

Ref#	Description	Size	Material*
25CU	Lead-Free	17" x 4" Child	Vinyl/Nylon
25C	Lead-Lined	17" x 4" Child	Vinyl/Nylon
25U	Lead-Free	22 ¾" x 4 ¼" Adult	Vinyl/Nylon
25	Lead-Lined	22 ¾" x 4 ¼" Adult	Vinyl/Nylon



\* Visit [www.palmerohealth.com](http://www.palmerohealth.com) for the complete x-ray ordering table and colors on Pages 34-35 of the Product Guide.

- For Vinyl or Nylon material aprons, specify the color description (i.e. 27UBLACK) with the part number.

- For silk screening designs, add "BEAR," "TOOTH," or "HAPPYFACE" (i.e. 27U after the Lead-Free and Lead-Lined part number).

Please note not all styles are available in Silk Screening Designs.

# CLING SHIELD® X-RAY APRONS

Find the right fit for your practice

Do not fold aprons as it can cause cracking within the lining. Use Palmero's assortment of hangers to protect your investment in safety.

## Standard Apron Hanger

- White-coated steel hanger
- Compatible with Patient and Protectall style aprons
- Hanging hardware included

Ref#	Description	Size	Color
28	Standard Apron Hanger	15" x 2"	White



## Economy Apron Hanger

- White-coated steel hanger
- Not compatible with Dual, Coat, or Tech aprons
- Hanging hardware included

Ref#	Description	Size	Color
31	Economy Apron Hanger	15" x 2"	White



## CC Style Apron Hanger

- Sturdy chrome-coated hanger
- Compatible with Patient and Protectall style aprons
- Hanging hardware included

Ref#	Description	Size	Color
32	CC Style Apron Hanger	15"	Chrome



## Economy X-Ray Apron Hanger

- White-coated steel hanger
- Compatible with Patient, Protectall, Coat, and Dual style aprons

Ref#	Description	Size	Color
33	Economy X-Ray Apron Hanger	20 7/8" x 10 1/4"	White



## Deluxe Coat Panoramic Apron Hanger

- Sturdy chrome-coated hanger
- Compatible with Panoramic, Coat, and Dual style aprons

Ref#	Description	Size	Color
30	Deluxe Coat Panoramic Apron Hanger	21" x 9 1/4"	Chrome



# CLING SHIELD® X-RAY APRONS

Go green with our eco-friendly and environmentally friendly X-Ray Apron Recycling/Trade Program



RECYCLE/TRADE-IN  
PROGRAM

Buy 1, Recycle 1,  
Receive 1 apron hanger

1. **Purchase** a new Palmero x-ray apron.

2. **Collect** old or unwanted X-ray apron of any brand.

3. **Contact** Customer Service for a Return of Authorization Number (RA#) to register for recycling.

4. **Submit** invoice with RA# via Fax 203.377.8988 or email [customerservice@palmerohealth.com](mailto:customerservice@palmerohealth.com). Include ship to address, phone, and email.

5. **Ship** your old apron with RA# for responsible disposal to:  
Palmero Healthcare LLC Attn: X-Ray Recycle/Trade-In Program 120  
Goodwin Place Stratford, CT 06615

6. **Acquire** a new x-ray apron hanger, an \$18.00 retail value for your participation!



For more information, visit [palmerohealth.com](http://palmerohealth.com),  
call 800.344.6424 or email [customerservice@palmerohealth.com](mailto:customerservice@palmerohealth.com).



*Protecting Clinicians, Patients, and the Practice*