

2~8°C UNDER-COUNTER PHARMACY REFRIGERATOR



Model: HYC-118A



AUDIBLE & VISUAL ALARMS

Audible buzzer and visible flashing light



3 ALARM FUNCTIONS

High/low temperature, door ajar, sensor error



AIR COOLING / AUTO-DEFROST

Better cabinet temperature uniformity, auto defrost

PRODUCT SPECIFICATIONS

TECHNICAL INFORMATION

Type	Pharmacy Refrigerator
Cabinet Type	Undercounter
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise Level (Db)	30

PERFORMANCE

Temperature Range (°C)	2 ~ 8
------------------------	-------

CONTROL

Controller	Microprocessor
Display	LED

ELECTRICAL DATA

Power Supply (V/Hz)	110/60
Power (W)	160
Electrical Current (A)	1.05



DIMENSIONS

Capacity (L/Cu.ft)	118/4.2
Net/Gross Weight (Approx)	46/51 kg 101.4/112.2 lbs
Interior Dimensions (W/D/H)	512 / 415 / 638 mm
Exterior Dimensions (W/D/H)	597 / 645 / 835 mm
Packing Dimensions (W/D/H)	675 / 680 / 930 mm

FUNCTIONS

High/Low Temperature	Y
Remote Alarm	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

ACCESSORIES

Caster	NA
Foot	Y
Test Hole	Y
Shelves/Basket Quantity	3 / 1
USB Interface	NA