

# 2~8°C PHARMACY REFRIGERATOR



## Model: HYC-290



### ADJUSTABLE SHELVES

Adjustable to any height for your requirements



### AIR COOLING / AUTO-DEFROST

Better cabinet temperature uniformity, auto defrost



### USB INTERFACE

Data of temperature history can be downloaded via the built-in USB port, optional print feature

## PRODUCT SPECIFICATIONS



### TECHNICAL INFORMATION

Type	Pharmacy Refrigerator
Cabinet Type	Upright
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise Level (Db)	43

### PERFORMANCE

Temperature Range (°C)	2 ~ 8
------------------------	-------

### CONTROL

Controller	Microprocessor
Display	LED

### ELECTRICAL DATA

Power Supply (V/Hz)	110/60
Power (W)	340
Electrical Current (A)	2.2

### DIMENSIONS

Capacity (L/Cu.ft)	290/10.2
Net/Gross Weight (Approx)	94/116 kg   207.2/255.7 lbs
Interior Dimensions (W/D/H)	530 / 555 / 1080 mm
Exterior Dimensions (W/D/H)	665 / 710 / 1665 mm
Packing Dimensions (W/D/H)	735 / 800 / 1800 mm

### FUNCTIONS

High/Low Temperature	Y
Remote Alarm	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

### ACCESSORIES

Caster	Y
Foot	Y
Test Hole	Y
Shelves/Inner Doors	4 / -
USB Interface	Y