SAFETY DATA SHEET



Accel® INTERVention Wipes One Step Surface Cleaner and Disinfectant

Revision: 2021-12-07

Version:

1. IDENTIFICATION

Product name:

SDS #: Recommended use: Uses advised against:

Manufacturer, importer, supplier: US Headquarters Diversey, Inc. 1300 Altura Rd., Suite 125 Fort Mill, SC 29708 Phone: 1-888-352-2249 SDS Internet Address: https://sds.diversey.com

Emergency telephone number:

Accel® INTERVention Wipes One Step Surface Cleaner and Disinfectant MS0301550 • This product is intended to be used neat. Uses other than those identified are not recommended

> Canadian Headquarters Diversey Canada, Inc. 6150 Kennedy Road Unit 3 Mississauga, Ontario L5T 2J4 Phone: 1-800-668-7171

1-800-851-7145; 1-651-917-6133 (Int'l)

2. HAZARDS IDENTIFICATION

Classification for the undiluted product

This product is not classified as hazardous according to OSHA 29CFR 1910.1200 (HazCom 2012-GHS) and Canadian Hazardous Products Regulations (HPR) (WHMIS 2015-GHS).

Hazard Statements None required. Precautionary Statements None required.

<u>Health hazards not otherwise classified (HHNOC)</u> - Not applicable <u>Physical hazards not otherwise classified (PHNOC)</u> - Not applicable

RTU

3. COMPOSITION/INFORMATION ON INGREDIENTS

Classified Ingredients

Ingredient(s)	CAS #	Weight %
Benzyl alcohol	100-51-6	1 - 5%
Hydrogen peroxide	7722-84-1	> 0.1 - < 1%

4. FIRST AID MEASURES

Undiluted Product:

Eyes: Rinse with plenty of water. If irritation occurs and persists, get medical attention. **Skin:** No specific first aid measures are required. **Inhalation:** No specific first aid measures are required. **Ingestion:** IF SWALLOWED: Call a Poison Center (1-800-851-7145) or doctor/physician if you feel unwell.

<u>Most Important Symptoms/Effects:</u> No information available. Immediate medical attention and special treatment needed Not applicable.

5. FIRE-FIGHTING MEASURES

Specific methods: Suitable extinguishing media: Specific hazards: No special methods required The product is not flammable. Extinguish fire using agent suitable for surrounding fire. None known.

Special protective equipment for firefighters: As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear. **Extinguishing media which must not be used for safety reasons:** No information available.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions: Environmental precautions and clean-up methods:

Put on appropriate personal protective equipment (see Section 8.). Clean-up methods - large spillage. Absorb spill with inert material (e.g. dry sand or earth), then place in a chemical waste container. Use a water rinse for final clean-up.

7. HANDLING AND STORAGE

Handling: Avoid contact with eyes. FOR COMMERCIAL AND INDUSTRIAL USE ONLY. **Storage:** Keep tightly closed in a dry, cool and well-ventilated place.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Exposure Guidelines:

Ingredient(s)	CAS #	ACGIH	OSHA
Hydrogen peroxide	7722-84-1	1 ppm (TWA)	1 ppm (TWA)
			1.4 mg/m³ (TWA)

Undiluted Product:

Engineering measures to reduce exposure:

Good general ventilation should be sufficient to control airborne levels.

Personal Protective Equipment

It is the responsibility of the employer to determine the potential risk of exposure to hazardous chemicals for employees in the workplace in order to determine the necessity, selection, and use of personal protective equipment.

Eye protection:	No personal protective equipment required under normal use conditions.
Hand protection:	No personal protective equipment required under normal use conditions.
Skin and body protection:	No personal protective equipment required under normal use conditions.
Respiratory protection:	No personal protective equipment required under normal use conditions. If aerosols, mists, or vapors are not adequately controlled by ventilation, use appropriate respiratory protection to avoid
Hygiene measures:	over-exposure. Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Physical StateLiquidEvaporation Rate:No information availableOdor threshold:No information available.Melting point/freezing point (°C):Not determinedAutoignition temperature:No information availableSolubility in other solvents:No information availableDensity:1.01Kg/L

Color: Clear Odor: Product specific Surfactant Boiling point/range: Not determined Decomposition temperature: Not determined Solubility: Completely Soluble Relative Density (relative to water): 1.01 Vapor density: No information available Bulk density: No information available Flash point (°F): > 200 °F > 93 °C Viscosity: 0 VOC: 0 % * Flammability (Solid or Gas): Not applicable Sustained combustion: Not applicable Explosion limits: - upper: Not determined - lower: Not determined

Vapor pressure: No information available. Partition coefficient (n-octanol/water): No information available Elemental Phosphorus: 0.12 % by wt. pH: ≈ 3 Corrosion to metals: Not corrosive to metals

* - Title 17, California Code of Regulations, Division 3, Chapter 1, Subchapter 8.5, Article 2, Consumer Products, Sections 94508

10. STABILITY AND REACTIVITY

Reactivity: Stability: Hazardous decomposition products: Materials to avoid: Conditions to avoid: Not Applicable The product is stable None reasonably foreseeable. Do not mix with any other product or chemical unless specified in the use directions. None known.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure: Skin contact, Inhalation, Eye contact

Delayed, immediate, or chronic effects and symptoms from short and long-term exposure

BW).

Skin contact: Unlikely to be irritant in normal use.
Eye contact: May be mildly irritating to eyes.
Ingestion: No information available.
Inhalation: No information available.
Sensitization: No known effects.
Target Organs (SE): None known
Target Organs (RE): None known

Numerical measures of toxicity
ATE - Oral (mg/kg):
ATE - Dermal (mg/kg):
ATE - Inhalatory, gases (mg/l):
ATE - Inhalatory, mists (mg/l):
ATE - Inhalatory, vapors (mg/l):

>20

The estimated acute oral LD50 for the undiluted product is estimated to be >5000 (mg/kg

12. ECOLOGICAL INFORMATION

Ecotoxicity: No information available.

Persistence and Degradability: No information available.

Bioaccumulation: No information available.

Mobility: No information available.

Other adverse effects No information available.

13. DISPOSAL CONSIDERATIONS

Do not contaminate water, food, or feed by storage or disposal.

Waste from residues / unused products (undiluted product):

This product, as sold, if discarded or disposed, is not a hazardous waste according to Federal regulations (40 CFR 261). Under RCRA, it is the

responsibility of the user of the product to determine, at the time of disposal, whether the waste solution meets RCRA criteria for hazardous waste. Dispose in compliance with all Federal, state, provincial, and local laws and regulations.

Pesticide Storage: Refer to product label.

Pesticide Disposal: Refer to product label.

Container Disposal: Refer to product label.

RCRA Hazard Class (undiluted product): Not Regulated. Contaminated Packaging: Do not re-use empty containers.

14. TRANSPORT INFORMATION

DOT/TDG/IMDG: The information provided below is the full transportation classification for this product. This description does not account for the package size(s) of this product, that may fall under a quantity exception, according to the applicable transportation regulations. When shipping dangerous goods, please consult with your internal, certified hazardous materials specialist to determine if any exceptions can be applied to your shipment.

DOT (Ground) Bill of Lading Description: NOT REGULATED

IMDG (Ocean) Bill of Lading Description: NOT REGULATED

15. REGULATORY INFORMATION

International Inventories at CAS# Level

DSL

All components are listed or otherwise exempt

CERCLA/ SARA

Ingredient(s)	CAS #	Weight %	CERCLA/SARA RQ (lbs)	Section 302 TPQ (lbs)	Section 313
Hydrogen peroxide	7722-84-1	> 0.1 - < 1 %		1000	

Canadian Regulations

This product is a Canadian DIN registered product and is subject to certain labeling requirements for disinfectants under national regulations. These requirements differ from the classification criteria and hazard information required for safety data sheets (SDS), and workplace labels of non-disinfectant products by the HPR (WHMIS 2015). This SDS was prepared in accordance with the HPR (WHMIS 2015) criteria, which may result with differences between the label and SDS language.

DIN No.: 02448637

16. OTHER INFORMATION

NFPA (National Fire Protection Association) Rating Scale: (Low Hazard) 0 - 4 (Extreme Hazard)

Health 0 Flammability 0 Instability 0 Special Hazards -

Revision: 2021-12-07 Version:

Reason for revision: Prepared by: Not applicable North American Regulatory Affairs

Additional advice:

Does not contain an added fragrance

Notice to Reader: This document has been prepared using data from sources considered technically reliable. It does not constitute a warranty, express or implied, as to the accuracy of the information contained within. Actual conditions of use and handling are beyond seller's control. User is responsible to evaluate all available information when using product for any particular use and to comply with all Federal, State, Provincial and Local laws and regulations.