



### ADJUSTABLE SHELVES

Particle resistant epoxy coated wire shelves, adjustable height.

### FORCED AIR COOLING WITH AUTO DEFROST

Better cabinet temperature uniformity and defrost

### REMOTE ALARM FUNCTIONS

High/Low temperature and sensor error.

### STANDARD FEATURES

Self closing door, probe test hole and Interior LED lighting

The MDF-20U110 is the most versatile Pharmacold freezer model yet. It is equipped with a reversible door which can be mounted on the left or right side. It has also been specifically designed to fit under a counter or in tight spaces.

Unlike other medical freezers its size the air intake and discharge vents are located in the front of the freezer. Giving it the ability to be placed directly against a wall, no need for added space at the back for air flow. The Pharmacold MDF-20U110 is an asset to any medical institution or laboratory.

### TECHNICAL INFORMATION

Type	Pharmacy Freezer
Cabinet Type	Under Counter
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	R600a

### PERFORMANCE

Temperature Range Celsius	-10 to -25
Accuracy Celsius	0.1
Uniformity Celsius	+ /- 2

### CONTROL

Controller	Microprocessor Thermostat
Display	LED

### ELECTRICAL DATA

Power Supply (V/HZ)	110/60
Power (W)	138
Electrical Current (A)	1.2

### DIMENSIONS

Capacity (L/Cu.ft)	113/4
Net/Gross Weight (Approx)	51 lbs
Exterior Dimensions (W/D/H)	23.5 / 22/ 31.75 inches
Interior Dimensions (W/D/H)	19.75 / 16.75 / 22.5 inches

### FUNCTIONS

High / Low Temperature	Y
Sensor Error	Y

### ACCESSORIES

Caster	Optional
Foot	4
Probe Hole	Standard
Shelves/Basket Quantity	3/-
Door Lock	Standard
Remote Alarm Contact	Standard