



### ADJUSTABLE SHELVES

Particle resistant epoxy coated wire shelves, adjustable height.

### FORCED AIR COOLING WITH AUTO DEFROST

Better cabinet temperature uniformity and defrost

### REMOTE ALARM FUNCTIONS

High/Low temperature, power failure, door ajar / door open alert, sensor error, low battery alarms

### STANDARD FEATURES

Backup battery for alarms, self closing door, probe test hole, Interior LED lighting

The NAB-60 is the most versatile Pharmacold model yet. It is equipped with a reversible door which can be mounted on the left or right side. The Pharmacold NAB-60 is an asset to any medical institution or laboratory.

## TECHNICAL INFORMATION

Type	Pharmacy Refrigerator
Cabinet Type	Under Counter
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free

## PERFORMANCE

Temperature Range Celsius	2 - 8
Accuracy Celsius	0.1
Uniformity Celsius	+ /- 2

## CONTROL

Controller	Microprocessor Thermostat
Display	LED

## ELECTRICAL DATA

Power Supply (V/HZ)	120/60
Power (W)	96
Electrical Current (A)	0.8

## DIMENSIONS

Capacity (L/Cu.ft)	60L / 2.1 cu ft
Net/Gross Weight (Approx)	83.8/88.2 (lbs)
Interior Dimensions (W/D/H)	16.25" / 15" / 16.5" inches
Exterior Dimensions (W/D/H)	20.13" / 22.13" / 26" inches

## FUNCTIONS

High / Low Temperature	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

## ACCESSORIES

Caster	N
Foot	4
Probe Hole	Standard
Shelves/Basket Quantity	1.5/-
Door Lock	Standard
Remote Alarm Contact	Standard