

# AAMI 2 OVERHEAD ISOLATION GOWN

THUMBLOOP WITH ELASTIC CUFF | FULL BACK | YELLOW

REF PM4-20662

## ABOUT THE GOWN

Made from strong, fluid resistant and breathable SMS material, this gown is ideal for situations with potential for low or moderate levels of fluid exposure. The overhead neck closure, full back design and thumbloops with elastic cuffs provide dependable coverage against contaminants.

## PROTECTION STANDARDS



Impact Penetration:  $\leq 1.0$  g  
Hydrostatic Pressure:  $\geq 20$  cm

Variation in lot results exist.  
All results meet or exceed  
AAMI PB70 requirements.

## FEATURES & BENEFITS

### PRIMAGARD® SMS TECHNOLOGY

PRIMED SMS gowns are engineered for the ultimate balance of breathability and fluid resistance to keep clinicians protected from the outside and comfortable on the inside.

### OVERHEAD NECK CLOSURE


Allows quick and easy donning and doffing, minimizing contact with harmful fluids.

### FULL BACK DESIGN

Provides enhanced coverage clinicians can depend on to prevent contact with contaminants.

### THUMBLOOP WITH ELASTIC CUFF

Keeps gown sleeves in place under gloves to protect vulnerable skin from harmful substances and bacteria.



## QUANTITIES

GOWNS/BAG	BAGS/CASE	GOWNS/CASE
10	10	100

WEIGHT  
Medium

BREATHABILITY  
●●●○○

FLUID RESISTANCE  
●●●○○



PRIMED Medical Products Inc.  
200, 2003-91 Street SW, Edmonton, AB, Canada, T6X 0W8  
TF +1.877.877.4633 | F +1.780.497.7670 | W PRIMED.ca | E info@PRIMED.ca

NA-REVA1