



BETTER FIT BETTER PROTECTION




CA-N95 FLAT- FOLD RESPIRATOR (N95 EQUIVALENT, NOT NIOSH N95)



 **4 LAYERS**
PROTECTION

PFE & BFE
≥ 95%



-  Made by Canadian for Canadian
-  Innovative 3D shaped design that provides an enhanced seal around the face area, and helps prevent airborne bacteria contamination.
-  The expanded area surface increases airflow for easier breathability.

-  4 available sizes to help you find your perfect fit.
-  A double layer of highly breathable meltblown effectively filters harmful particles and viruses by ≥ 95%.
-  Extra soft fabric provides long-lasting comfort.

SUPPORT CANADIAN MADE PPE

TECHNICAL DATA SHEET

CA-N95 FLAT- FOLD RESPIRATOR

(N95 EQUIVALENT, NOT NIOSH N95)

ZIPLOCK BAG PACKAGING

COMPACT
REUSABLE
WATERPROOF



RETAIL ZIPLOCK BAG PACKAGING SAMPLE



DESCRIPTION & CLASSIFICATION

PRODUCT NAME	CA-N95 RESPIRATOR
STYLE	FLAT FOLD
STRAP TYPE	FLAT EARLOOP
SIZE	XS, S, M, L
COLOURS	BLACK, WHITE, BLUE
MANUFACTURED IN	ONTARIO, CANADA


PACKAGING SPECIFICATION


RESPIRATOR WEIGHT	0.02 lbs / 9 g
RESPIRATORS PER BAG	10 PC
WEIGHT PER BAG	0.18 lbs / 82g
PACKAGE DIMENSIONS	3.17" X 1.5" X 9.5"
RESPIRATORS PER CASE	1250
BAGS PER CASE	125
WEIGHT PER CASE	22.5 lbs / 10.25 kg
CASE DIMENSIONS	20.25" X 16" X 15"
RESPIRATORS PER SKID	25000 PC
CASE PER SKID	20 CASES
SKID DIMENSIONS	48" X 40" X 60"

PERFORMANCE REQUIREMENTS		
TEST	HEALTH CANADA REQUIREMENTS	CA-N95
PARTICLE FILTER EFFICIENCY	≥95 %	PASS
INHALATION RESISTANCE	≤35 MMH ₂ O (343 PA)	PASS
EXHALATION RESISTANCE	≤25 MMH ₂ O (245 PA)	PASS
MECHANICAL HEADSTRAP STRENGTH	≥20 (NEWTONS)	PASS

	XS (SUITABLE FOR CHILDREN AGES 4-11)			
	WIDTH	179 MM / 7.05"		
	COLOURS	WHITE	BLACK	BLUE
	ID#	FF-WHT-XS	FF-BLK-XS	FF-BLU-XS

	S (SUITABLE FOR OLDER CHILDREN)			
	WIDTH	189 MM / 7.44"		
	COLOURS	WHITE	BLACK	BLUE
	ID#	FF-WHT-S	FF-BLK-S	FF-BLU-S

	M (SUITABLE FOR ADULTS WITH SMALLER FACES)			
	WIDTH	199 MM / 7.83"		
	COLOURS	WHITE	BLACK	BLUE
	ID#	FF-WHT-M	FF-BLK-M	FF-BLU-M

	L (SUITABLE FOR MOST ADULTS)			
	WIDTH	215 MM / 8.46"		
	COLOURS	WHITE	BLACK	BLUE
	ID#	FF-WHT-L	FF-BLK-L	FF-BLU-L