


# MOBIDERM STOCKINGS

## POSITION OF THE ANATOMICAL MARKERS

- G - is located on the gluteal fold (for thigh-high).....
- F - equidistant point between G and E .....
- E - corresponds to the middle of the patella .....
- D - is located on the head of the fibula (3-4 cm under the patella).....
- C - corresponds to the widest point of the calf.....
- B1 - corresponds to the start of the calf muscles.....
- B - to be measured at the narrowest point of the ankle .....
- Y - goes through the instep and the tuberosity of the calcaneus.....
- A - goes through the head of the metatarsal bones (at the widest point).....



All markers are positioned on the external face of the leg



Circumference of the head of the metatarsal bones



Circumference at the instep



Ankle circumference (at the narrowest point)



Circumference of the thigh (at the gluteal fold)



Height measurement (position of measuring tape on the external face of the leg)



- ORDER (by default)
- QUOTATION
- RENEWAL



RETAILER IDENTIFICATION

Customer code:

Patient's surname: .....

Patient's first name: .....

Gender :  M  F  Child

Patient's height: .....

Case No. for renewal

1st treatment

Date : ..... Quantity: .....

I authorize my health care professional to collect my data and to communicate them to Thuasne company as part of the processing of my made to measure medical device in accordance with Law No 78-17 of 6 January 1978, and European Regulation No 2016/679/ EU of 27 April 2016, I have rights including in particular the rights of access, rectification, portability and deletion of my data. I can exercise these rights by contacting the health care professional to who I ordered my medical device.

Patient signature

**If possible, please enclose photos of the limb to be fitted.  
Please draw in the contours of the garment on the diagram and cross unnecessary measures.**

**Models**

- Below-knee
- Thigh-high

- |                          |                          |
|--------------------------|--------------------------|
| LEFT                     | RIGHT                    |
| <input type="checkbox"/> | <input type="checkbox"/> |
| <input type="checkbox"/> | <input type="checkbox"/> |

**Distal options**

- Open toe
- Closed toe

**Foam options**

- All small blocks
- All big blocks
- Big blocks on the leg and small blocks on the feet

**Proximal options**


- Anti-slip with silicone dots 3 cm
- Velcro opening

**Comments**

.....

.....

.....



**Please contact your regular Thuasne distributor**