

# PuraGraft Collagen Wound Dressings



## PuraMEND adapt

Resorbable collagen membrane packaged privately for PuraGraft

ZERO memory, memorable handling.

Type I bovine collagen

Its zero memory characteristics allow it to conform to any surgical defect when it becomes hydrated.

### Features

- Adapts and adheres effortlessly to surgical site
- Guides healing of bone and soft tissue
- Predictable resorption time (26-38 weeks)

<b>PGBA15</b>	<b>15x20</b>
<b>PGBA20</b>	<b>20x30</b>
<b>PGBA30</b>	<b>30x40</b>



## PuraMEND maintain

Resorbable collagen membrane packaged Privately for PuraGraft

- Type I bovine collagen.
- Resorption time, 26-38 weeks, for sustained function in bone repair.
- Higher mechanical strength for membrane stabilization.
- Each piece can be easily trimmed for a final fit.
- May be placed with either side against the tissue.

<b>PGBM15</b>	<b>15x20</b>
<b>PGBM20</b>	<b>20x30</b>
<b>PGBM30</b>	<b>30x40</b>



*Call (877) 540-3258*