

LAB FILTRATION FOR AIR MONITORING

Looking to filter out unnecessary expenses??

Let us do the heavy lifting!

VISIT US:

www.delta-sci.com

Connect with us now! 905-629-4545

infoedelta-sci.com

Deal For You Until
Sept 30th
FREE
SAMPLES

50% OFF

of equal or lower cost

Delta Scientific Laboratory Products 346 Watline Avenue Mississauga, ON Canada

Share your filter usage, and we'll find affordable alternatives, freeing up your budget for other critical investments



Opticlear Mixed Cellulose Esters (MCE) Membrane Filters

- Airborne Asbestos Monitoring
- Sampling and Analysis of Si & Pb
- Sampling of bulk-air for TSP (Total Suspended Particulate Matter)



Glass & Quartz **Thimbles**

• High temperature extraction or dust monitoring.



Binder-free Glass Microfiber Filters (Grade-GB100IR)

• Fiilters with very low trace metal content of As, Pb, Cd used in high and low volume air samplers for aerosols, airborne dust and metal contaminants.





Quartz Fiber Filters

• Filter media for air pollution analysis to sample acid gases in stacks and flues at high temperatures >500°C



PTFE Composite Filters (PG60)

• Measurement of dust concentration in the air



37mm 2PC/3PC Cassette

• Fungi Collection, Microbiological Analysis & Asbestos Monitoring.

